

## Certificate of Analysis for NR-52902

# Vector pMCSG53 Containing the SARS-Related Coronavirus 2, Wuhan-Hu-1 Non-Structural Protein 8 Gene

### Catalog No. NR-52902

#### **Product Description:**

The non-structural protein 8 (nsp8) gene from severe acute respiratory syndrome-related coronavirus 2 (SARS-CoV-2), Wuhan-Hu-1 (GenBank: MN908947) was codon optimized and cloned into the pMCSG53 plasmid. pMCSG53 is an Escherichia coli (E. coli) expression vector that contains an N-terminal hexa-histidine tag, followed by a tobacco etch virus (TEV) protease recognition site prior to the insert coding sequence, resulting in the expression of a cleavable histidine-tagged protein. It also contains tRNA genes covering rare codons for arginine (AGG/AGA) and isoleucine (AUA) to improve expression in E. coli. The beta-lactamase gene, bla, provides transformant selection through ampicillin resistance in E. coli. The deposited plasmid was transformed into One Shot™ TOP10 E. coli (Invitrogen™ C404003), grown in Luria-Bertani broth with ampicillin (50 µg per mL) for 1 day at 37°C in an aerobic atmosphere, extracted using a Plasmid Plus Maxi Kit (QIAGEN® 12965) and vialed in TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0).

Lot: 70036144 Manufacturing Date: 21MAY2020

TEST	SPECIFICATIONS	RESULTS
Next-Generation DNA Sequencing	~ 5390 base pairs	5382 base pairs <sup>1</sup>
Genotypic Analysis		
Sequencing of nsp8 gene (~ 594 base pairs)	≥ 99% sequence identity to depositor's sequence	100% sequence identity to depositor's sequence <sup>2</sup>
Sequencing of pMCSG53 vector (~ 4790 base pairs)	N-terminal hexa-histidine tag confirmed	N-terminal hexa-histidine tag confirmed
	N-terminal TEV protease cleavage site confirmed	N-terminal TEV protease cleavage site confirmed
Antibiotic Resistance		
Ampicillin (encoded by beta-lactamase gene bla) <sup>3</sup>	bla sequence present	bla sequence present
Agarose Gel Electrophoresis		
Digestion with Sapl and BamHI (pre-vial)	~ 3 kb and ~ 2.5 kb	~ 3 kb and ~ 2.5 kb (Figure 1)
Concentration by Qubit™ Measurement	≥ 2 µg/mL	0.3 μg in 20 μL per vial (15 μg/mL)
Amount per Vial	Report results	0.3 μg per vial
OD <sub>260</sub> /OD <sub>280</sub> Ratio	1.7 to 2.1	1.9
Effective Bacterial Transformation		
Invitrogen™ One Shot™ TOP10 <i>E. coli</i>	≥ 50 colonies per ng	266 colonies per ng

<sup>&</sup>lt;sup>1</sup>The sequence was assembled pre-vial using the depositor's predicted sequence as the reference sequence. The complete plasmid sequence and map are provided on the BEI Resources webpage.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org Tel: 800-359-7370

Tel: 800-359-7370 Fax: 703-365-2898


<sup>&</sup>lt;sup>2</sup>The NR-52902 insert was codon optimized but is 100% identical with the SARS-CoV-2, Wuhan-Hu-1 NSP8 protein within the ORF1ab polyprotein (GenPept: QHD43415).

<sup>&</sup>lt;sup>3</sup>The antibiotic ampicillin degrades quickly during growth. Bacterial stationary phase should be minimized during plasmid expansion to avoid plasmid loss and increased antibiotic concentrations may be necessary.



## Certificate of Analysis for NR-52902

Figure 1: Agarose Gel of Undigested and Restriction Enzyme Digested NR-52902



Lane 1: Invitrogen™ TrackIt™ 1 Kb Plus DNA Ladder

Lane 2: NR-52902 digested Lane 3: NR-52902 undigested

/Heather Couch/ Heather Couch

21 AUG 2020

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by ATCC® and the contributor to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

 $\mbox{ATCC}^{\mbox{\tiny{\$}}}$  is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: <a href="mailto:contact@beiresources.org">contact@beiresources.org</a>
Tel: 800-359-7370

Fax: 703-365-2898