

Genomic RNA from SARS-Related Coronavirus 2, Isolate USA-CA4/2020

Catalog No. NR-52508

This reagent is the property of the U.S. Government.

Product Description:

Genomic RNA was extracted from a preparation of cell lysate and supernatant from *Cercopithecus aethiops* kidney epithelial cells (Vero E6; ATCC® CRL-1586™) infected with severe acute respiratory syndrome-related coronavirus 2 (SARS-CoV-2), isolate USA-CA4/2020 (BEI Resources NR-52386 lot 70034885), using QIAamp® Viral RNA Mini Kit (Qiagen® 52904). The viral genomic RNA is in a background of cellular nucleic acid and carrier RNA.

Lot: 70035259


Manufacturing Date: 28MAY2020

TEST	SPECIFICATIONS	RESULTS
Sequencing of Species-Specific Region (860 nucleotides)	≥ 98% identity with SARS-CoV-2, isolate USA-CA4/2020 (GenBank: MT027063.1)	100% identity with SARS-CoV-2, isolate USA-CA4/2020 (GenBank: MT027063.1)
Next-Generation Sequencing (NGS) of Complete Genome Using Illumina® iSeq™ 100 Platform (Refer to Appendix I for NGS information)	≥ 98% identity with SARS-CoV-2, isolate USA-CA4/2020 (GenBank: MT027063.1)	99.99% identity with SARS-CoV-2, isolate USA-CA4/2020 (GenBank: MT027063.1)
Functional Activity by RT-PCR Amplification¹	~ 950 base pair amplicon	~ 950 base pair amplicon (Figure 1)
Genome Copy Number Using BioRad QX200 Droplet Digital PCR (ddPCR™) System (Post vial; 6 replicates)	Report results	2.67 × 10 ⁸ genome equivalents/mL
Virus Inactivation (14 Days, 2 Passages) 10% of total yield inoculated on Vero E6 cells and evaluated for cytopathic effect ² Passage 1 Passage 2	No viable virus detected No viable virus detected	No viable virus detected No viable virus detected
Virus Inactivation (14 Days, 2 Passages) 100% of total yield inoculated on Vero E6 cells and evaluated for cytopathic effect ² Passage 1 Passage 2	No viable virus detected No viable virus detected	No viable virus detected No viable virus detected

¹Amplified using iTaq™ Universal SYBR Green One-step Kit (Bio-Rad® 172-5151) with 5 µL of NR-52508 in a 50 µL reaction

²Supernatant was tested by qPCR to confirm absence of CPE and no evidence of replicative RNA was detected. Samples from both passages were tested by qPCR at the end of day 14 of the passages.

Figure 1: Functional Activity of NR-52508 by RT-PCR Amplification



Lane 1: Invitrogen™ TrackIt™ 1 Kb Plus DNA Ladder
Lane 2: PCR product from 1 µL of NR-52508

/Heather Couch/
Heather Couch

08 OCT 2020

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.



ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

APPENDIX I: NGS Information for NR-52508 lot 70035259

Sequence analysis resulted in the discovery of three SNPs when compared to GenBank MT027063.1 (see Table below). Quality scores over 60 indicate it is improbable that the variant call is incorrect.

Position in NR-52508_70035259 Sequence	Position in MT027063.1	Reported MT027063.1 Sequence	Identified Alternative Base	Quality	Variant Type	Length of Variant	Frequency of Variant
22195	22206	A	G	53	SNP	1	0.1889764
23596	23607	G	T	221	SNP	1	0.8170213
26187	26198	C	T	221	SNP	1	0.7892562