

Vector pHDM Containing the SARS-Related Coronavirus 2, Wuhan-Hu-1 Spike Glycoprotein

Catalog No. NR-52514

Product Description:

The vector for the spike (S) glycoprotein gene from severe acute respiratory syndrome-related coronavirus 2 (SARS-CoV-2), Wuhan-Hu-1 (GenBank: [NC_045512](#)) was designed by codon optimizing the S sequence for mammalian expression and subcloned into the pHDM vector under the CMV promoter. The Kozak sequence in pHDM has been corrected. NR-52514 contains the beta-lactamase gene, *bla*, to provide transformant selection through ampicillin resistance in *Escherichia coli* (*E. coli*). The deposited plasmid was transformed into One Shot™ TOP10 *E. coli* (Invitrogen™ C404003), grown in Luria-Bertani broth with ampicillin (50 µg per mL) for 1 day at 37°C in an aerobic atmosphere, extracted using a Plasmid Plus Maxi Kit (QIAGEN® 12963) and vialled in TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0).

Lot: 70035472

Manufacturing Date: 29APR2020


TEST	SPECIFICATIONS	RESULTS
Next-Generation DNA Sequencing	~ 8380 base pairs	8386 base pairs ¹
Genotypic Analysis Sequencing of S glycoprotein insert (~ 3820 base pairs)	≥ 99% sequence identity to depositor's sequence	100% sequence identity to depositor's sequence ²
Antibiotic Resistance Ampicillin (encoded by beta-lactamase gene <i>bla</i>) ³	<i>bla</i> sequence present	<i>bla</i> sequence present
Agarose Gel Electrophoresis (pre-vial) Digestion with <i>SapI</i>	~ 5 kb and ~ 4 kb	~ 5 kb and ~ 4 kb (Figure 1)
Concentration by PicoGreen® Measurement	≥ 2 µg/mL	0.7 µg in 100 µL per vial (7 µg/mL)
Amount per Vial	Report results	0.7 µg per vial
OD₂₆₀/OD₂₈₀ Ratio	1.7 to 2.1	2.0
Effective Bacterial Transformation Invitrogen™ One Shot™ TOP10 <i>E. coli</i>	≥ 50 colonies per ng	129 colonies per ng

¹The sequence was assembled pre-vial using the depositor's predicted sequence as the reference sequence. The complete plasmid sequence and map are provided on the BEI Resources webpage.

²The NR-52514 insert was codon optimized for mammalian expression, but otherwise is 100% identical to the SARS-CoV-2, Wuhan-Hu-1 S protein (GenPept: YP_009724390).

³The antibiotic ampicillin degrades quickly during growth. Bacterial stationary phase should be minimized during plasmid replication to avoid plasmid loss and increased antibiotic concentrations may be necessary.

Figure 1: Agarose Gel of Undigested and Restriction Enzyme Digested NR-52514



Lane 1: Invitrogen™ TrackIt™ 1 Kb Plus DNA Ladder

Lane 2: NR-52514 digested

Lane 3: NR-52514 undigested

/Heather Couch/

Heather Couch

04 NOV 2020

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by ATCC® and the contributor to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

