

NIH AIDS Reagent Program

20301 Century Boulevard Building 6, Suite 200 Germantown, MD 20874 USA

Phone: 240 686 4740 Fax: 301 515 4015 aidsreagent.org

DATA SHEET

Reagent: * HHV7 JI Infected Sup-T1 Cells

Catalog Number: 12184

Lot Number: 102710

Release Category:

1x10⁶cells/vial Provided:

Sup-T1 cells infected with HHV7 JI strain. Cell Type:

Growth

NOTE: It is recommended that you have a culture of uninfected cells (Sup T-1, Cat# Characteristics: 100) prepared when thawing these cells. To expand this culture, mix the infected cells with uninfected cells at a 1:5 ratio of infected to uninfected. Infected cells tend to cause

clumps of cells with the infected cells in the center. Be sure to gently pipette the cells a

few times to break up clumps to allow a more robust expansion.

Description: Sup-T1 cells infected with HHV7 JI strain.

Special

Characteristics:

The virus is highly cell-associated and as such little virus released to the supernatant.

Recommended

Storage:

Liquid Nitrogen

Contributor: Dr. Dharam Ablashi, HHV-6 Foundation

References: Ablashi DV, Berneman ZN, Kramarsky B, Whitman J Jr, Asano Y, Pearson GR. Human

herpesvirus-7 (HHV-7): current status. Clin Diagn Virol. 1995 Jul;4(1):1-13.

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 07/05/2018 Page 1 of 2

NOTE:	Acknowledgement for publications should read "The following reagent was	
NOIL.	ACKIOWICUCEITETT TO DUDIICATIONS SHOULD FEAU THE TOHOWING FEAUCHT WAS	٠,

through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: HHV7 JI Infected Sup-T1 Cells from Dr. Dahram Ablashi." Also, include the reference cited above in any publications.

Last Updated July 05, 2018

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 07/05/2018 Page 2 of 2