

NIH AIDS Reagent Program

20301 Century Boulevard Building 6, Suite 200 Germantown, MD 20874 USA

Phone: 240 686 4740 Fax: 301 515 4015 aidsreagent.org

DATA SHEET

Reagent:	HTLV-II Cloning Vector (pGEMT-HTLV-2)
Catalog Number:	12496
Lot Number:	160169
Release Category:	E
Provided:	5 µg of dried purified DNA stabilized in DNAstable PLUS
Cloning Vector:	pGEMT-Easy
	Ampicillin resistant
Cloning Site:	TOPO TA cloned
	The size of the insert is 801 bp.
Description:	A HTLV-II cloning vector
Special Characteristics:	This construct is 3818 bp including the insert.
	The HTLV-II PCR amplified sequence was TOPO TA cloned between the EcoRI and SpeI sites of the pGEMT-Easy (Promega) vector.
	The sequence maps to partial LTR and gag sequences of the HTLV-II virus (isolate G12) between nucleotides 39 and 840 corresponding to partial LTR and gag sequences of the virus.
	This sequence can serve as a positive control in population screenings using the following primers: Forward: 5'-CGAGTCATCGACCCAAAAGGTC-3' Reverse: 5'-GGAGTTGGGGAAAGCCCGTGG-3'
	This construct is not for protein expression.
	GenBank Accession Number: L11456
	Plasmid map and sequence file lot 160169

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

	This reagent is currently being provided as dried purified DNA stabilized in DNAstable <i>PLUS</i> . Please see the notice for additional information and the protocol for reconstitution of dried DNA reagents. <u>Dried DNA Notice</u>
	Plasmids can be propagated in STBL2 cells and grown at 37°C. Larger plasmids may benefit from growth at 30°C. This construct may also be grown in other competent cells.
Recommended Storage:	Keep the reagent at room temperature in a dry storage cabinet or in a moisture barrier bag.
Contributor:	Dr. Elisa Oltra
References:	Oltra, E., Garcia-Escudero, M., Mena-Duran, A. V., Monsalve, V., & Cerda-Olmedo, G. (2013). Lack of evidence for retroviral infections formerly related to chronic fatigue in Spanish fibromyalgia patients. Virol J, 10, 332. doi: 10.1186/1743-422X-10-332 <u>PUBMED</u>
NOTE:	Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: HTLV-II Cloning Vector (pGEMT-HTLV-2) from Dr. Elisa Oltra." Also include the references cited above in any publications.
Last Updated	June 21, 2018

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.