

## NIH AIDS Reagent Program

20301 Century Boulevard Building 6, Suite 200 Germantown, MD 20874 USA

Phone: 240 686 4740 Fax: 301 515 4015 aidsreagent.org

## **DATA SHEET**

Reagent: ★ FeLV 61E Infected 3201B Cells

Catalog Number: 198

Lot Number: 99068

Release Category:

 $9 \times 10^6$  cells. Provided:

**Propagation** Medium:

RPMI 1640, 45%; Leibovitz's L-15, 45%; fetal bovine serum, 10%.

Freeze Medium: Propagation medium with 50% fetal bovine serum, 90%; DMSO, 10%.

Growth

3201B feline T cells grow as a clumpy suspension. Maintain at 0.5- $3.0 \times 10^6$  cells/ml. Characteristics: Passage every 4-6 days; the doubling time is 15-20 hours. These cells do not grow in

other media.

Sterility: Negative for aerobic and anaerobic bacteria, mycoplasma, fungi, and yeast.

**Description:** 3201B cells infected with FeLV 61E.

Special Characteristics:

Virus titer reaches 10<sup>5</sup> at best. FeLV 61/E induces persistent viremia and caused low

incidence, long latency lymphoma in domestic cats<sup>2</sup>. The plasmid clone p61E is

available as Catalog #109.

Recommended

Storage:

Liquid nitrogen.

**Contributor:** Dr. James I. Mullins.

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 07/02/2018 Page 1 of 2 **References:** 1Snyder HW Jr, Hardy WD Jr, Zuckerman EE, Fleissner E. Characterisation of a

tumour-specific antigen on the surface of feline lymphosarcoma cells. *Nature* 

**275**:656-658, 1978.

 $^2$ Overbaugh J, Donahue PR, Quackenbush SL, Hoover EA, Mullins JI. Molecular cloning of a feline leukemia virus that induces fatal immunodeficiency disease in cats. *Science* 

**239**:906-910, 1988.

NOTE: Acknowledgment for publications should read "The following reagent was obtained

through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: FeLV 61E Infected 3201B Cells from Dr. James Mullins." Also include the reference cited above

in any publications.

Last Updated July 02, 2018

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 07/02/2018 Page 2 of 2