



## NIH AIDS Reagent Program

20301 Century Boulevard  
Building 6, Suite 200  
Germantown, MD 20874  
USA

Phone: 240 686 4740  
Fax: 301 515 4015  
aidsreagent.org

### DATA SHEET

<b>Reagent:</b>	☒ HIV-1 HXB2 gp120 Secreting CHO Cells (CHO-SEC)
<b>Catalog Number:</b>	2237
<b>Lot Number:</b>	94003
<b>Release Category:</b>	D
<b>Provided:</b>	1 vial frozen cells.
<b>Cell Type:</b>	Derived from CHO-K1 cells (ATCC)
<b>Propagation Medium:</b>	See below. In addition to HIV-1 env genes, these cells have been stably transfected with a glutamine synthetase gene. Do not add glutamine to the culture medium, as this may select for cells that do not contain the desired env inserts.
<b>Freeze Medium:</b>	GMEM-S medium without MSX, 60%; fetal bovine serum, 30%; DMSO, 10%.
<b>Growth Characteristics:</b>	Split cells every 3-4 days at 1:12. Cells should just reach confluency on day of passage. Cells are heterogeneous and grow as a flat, adherent monolayer, singly or in clusters.
<b>Sterility:</b>	Negative for bacteria, fungi, and mycoplasma.
<b>Description:</b>	CHO cells that secrete cleaved and uncleaved HXB2 gp120
<b>Special Characteristics:</b>	CHO-K1 cells were cotransfected with HIV-1 env and rev expression vectors. The HXB2 env vector was mutated to include a stop codon just after the second amino acid residue before the predicted transmembrane domain, and was inserted into the vector pEE14 (Celltech), which expresses glutamine synthetase. The cells secrete equal amounts of uncleaved and cleaved envelope glycoprotein, on the order of 100 ng/ml supernatant from $5 \times 10^6$ cells. CHO-SEC expresses roughly five times more envelope glycoprotein than CHO-WT.

#### [Table 1. CHO-Cell Lines](#)

---

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

**Recommended Storage:** Liquid nitrogen.

**Contributor:** Dr. Carol Weiss and Dr. Judith White.

**References:** Weiss CD, White JM. Characterization of stable Chinese hamster ovary cells expressing wild-type, secreted, and glycosylphosphatidylinositol-anchored human immunodeficiency virus type 1 envelope glycoprotein. *J Virol* **67**:7060-7066, 1993.

**NOTE:** Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: HIV-1 HXB2 gp120 Secreting CHO Cells (CHO-SEC) from Dr. Carol Weiss and Dr. Judith White." Also include the reference cited above in any publications.

**Last Updated** July 03, 2018

---

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.