



NIH AIDS Reagent Program

20301 Century Boulevard
Building 6, Suite 200
Germantown, MD 20874
USA

Phone: 240 686 4740
Fax: 301 515 4015
aidsreagent.org

DATA SHEET

Reagent:	CHO-EE Cells
Catalog Number:	2238
Lot Number:	94004
Release Category:	D
Provided:	1 vial frozen cells. 1 x 10 ⁷ /vial, >98% viability
Cell Type:	Derived from CHO-K1 cells (ATCC).
Propagation Medium:	See below. These cells have been stably transfected with a glutamine synthetase gene. Do not add glutamine to the culture medium, as this may select for cells that do not contain the desired env inserts.
Freeze Medium:	GMEM-S medium without MSX, 60%; fetal bovine serum, 30%; DMSO, 10%.
Growth Characteristics:	Split cells every 3-4 days at 1:12. Cells should just reach confluency on day of passage. Cells are heterogeneous and grow as a flat, adherent monolayer, singly or in clusters.
Sterility:	Negative for mycoplasma, bacteria, and fungi
Description:	The cells can serve as a control for CHO-SEC and CHO-WT.
Special Characteristics:	CHO-K1 cells were transfected with pEE14 (Celltech), which expresses glutamine synthetase.

[Table 1. CHO-Cell Lines](#)

[Protocol: Culture and Syncytium Detection Using CHO-Env Cell Lines](#)

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

Recommended Storage: Keep the reagent in liquid nitrogen.

Contributor: Dr. Carol Weiss and Dr. Judith White.

References: Weiss CD, White JM. Characterization of stable Chinese hamster ovary cells expressing wild-type, secreted, and glycosylphosphatidylinositol-anchored human immunodeficiency virus type 1 envelope glycoprotein. *J Virol* **67**:7060-7066, 1993.

NOTE: Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: CHO-SEC from Dr. Carol Weiss and Dr. Judith White." Also include the reference cited above in any publications.

Last Updated August 20, 2020

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.