

## NIH AIDS Reagent Program

20301 Century Boulevard Building 6, Suite 200 Germantown, MD 20874 USA Phone: 240 686 4740 Fax: 301 515 4015 aidsreagent.org

## **DATA SHEET**

Reagent: FIV 34TF10 Pol Expression Vector

Catalog Number: 3716

**Lot Number:** 140029

Release Category: B

**Provided:** 5 μg of dried purified DNA stabilized in DNAstable *PLUS* 

Cloning Vector: PT7-7

Ampicillin resistant

**Description:** The FIV-34F10 pol gene was cloned into the PT7-7 vector as an EcoRI-BqlII fragment

(nt 1875-6460). An adaptor was added between NdeI and EcoRI to put pol in frame

for translation.

Special

**Characteristics:** 

This clone expresses all Pol proteins expressed by FIV PR, including PR, RT, DU, and IN. Expression is IPTG inducible. The polyprotein will facilitate structural, enzymatic, and drug efficacy studies. GenBank #M25381.1 for the complete FIV genome. The infectious molecular clone FIV-34F10 is also available.

Sequence file lot 140029

This reagent is currently being provided as dried purified DNA stabilized in DNAstable

PLUS. Please see the notice for additional information and the protocol for

reconstitution of dried DNA reagents. Dried DNA Notice

Recommended

Storage:

Keep the reagent at room temperature in a dry storage cabinet or in a moisture

barrier bag.

**Contributor:** Dr. John Elder

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 11/19/2018 Page 1 of 2

References: Talbott RL, et al. Proc Natl Acad Sci USA 86:5743, 1989.

> Elder JH, et al. J Virol 66:1791, 1992. Elder JH, et al. J Virol 67:1869, 1993.

NOTE:

Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: FIV 34TF10 Pol Expression Vector from Dr. John Elder (cat# 3716)." Also include the references cited

above in any publications.

Last Updated: November 19, 2018

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

REV: 11/19/2018 Page 2 of 2