



## NIH AIDS Reagent Program

20301 Century Boulevard  
Building 6, Suite 200  
Germantown, MD 20874  
USA

Phone: 240 686 4740  
Fax: 301 515 4015  
aidsreagent.org

### DATA SHEET

**Reagent:** Anti-HIV-1 IIIB Nef Monoclonal (AE6)

**Catalog Number:** 709

**Lot Number:** 080174

**Release Category:** D

**Provided:** 200 µg of purified antibody at 2 mg/mL

**Description:** A monoclonal antibody to HIV-1 IIIB Nef, specifically the C-terminus region (epitope: VARELHPEYFKNC)

**Host:** Mouse

**Titer:** The user should determine the optimal concentration for any application.

**Special Characteristics:** This antibody was produced using the mouse ascites method and purified by Protein A/G chromatography. This antibody originates from a hybridoma. The hybridoma was created by immunizing mice with a HIV-1 IIIB nef recombinant protein.

Please see the [LANL HIV Molecular Database](#) for more information.

The critical binding residue for this antibody is the "N" in the IIIB C-terminal "YFKNC" sequence. This antibody does not react with HIV-1 SF2 nef.

Applications: Western Blot, Immunoprecipitation, IFA

**Recommended Storage:** Keep the reagent at 4°C for short term storage and at -80°C for long term storage. Avoid freeze-thaw cycles as reagent degradation may result.

**Contributor:** Dr. James Hoxie

**Isotype:** IgG<sub>1</sub> κ

---

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

**References:**

Chang, A. H., Hoxie, J. A., Cassol, S., O'Shaughnessy, M., & Jirik, F. (1998). Construction of single-chain antibodies that bind an overlapping epitope of HIV-1 Nef. FEBS Lett, 441(2), 307-312. [PUBMED](#)

**NOTE:**

Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: Anti-HIV-1 IIIB Nef Monoclonal (AE6) from Dr. James Hoxie."

**Last Updated**

October 31, 2017

---

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.