



NIH AIDS Reagent Program
20301 Century Boulevard
Building 6, Suite 200
Germantown, MD 20874
USA

Phone: 240 686 4740
Fax: 301 515 4015
aidsreagent.org

DATA SHEET

Reagent: HIV-1 NL4-3 1499 Infectious Molecular Clone (p6463-13)

Catalog Number: 7401

Lot Number: 100211

Release Category: A

Provided: 5 µg of dried purified DNA stabilized in DNASTable *Plus*

Cloning Vector: pNLPFB
Ampicillin resistant

Cloning Site: MscI/PfIMI cloning site
The size of the insert is 871 bp.

GenBank: [AY351774](#)

Host Strain: Plasmids can be propagated in STBL2 cells and grown at 37°C. Larger plasmids may benefit from growth at 30°C. This construct may also be grown in other competent cells.

Description: A full length replication competent, infectious HIV-1 chimeric molecular clone. This plasmid contains the entire HIV-1 NL4-3 genome where the wild-type RT has been replaced with a mutant RT from patient that contains multiple drug resistance mutations.

Part of panel of 12 prototypical infectious multidrug resistant HIV-1 reverse transcriptase (RT) nucleoside/nucleotide RT inhibitor (NRTI) clones. The panel includes clones with each of the published nucleoside analog RT mutations in the combinations that occur most frequently in HIV-infected individuals.

The insert was amplified by RT-PCR from viral RT RNA isolated from patient plasma.

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

**Special
Characteristics:**

This construct is 14,825 bp including the insert.

The panel includes clones with each of the published nucleoside analog RT mutations in the combinations that occur most frequently in HIV-infected individuals. Nucleic acid sequence data of the RT as well as phenotype susceptibility results for each clone is available from the [Stanford HIV Drug Resistance Database](#).

[Table of Catalog Numbers for all members of the panel](#)

The panel can be used for the following purposes:

- 1) Screening and testing new compounds designed to be effective against the most commonly isolated multidrug resistant variants.
- 2) Biochemical and biophysical studies that require a representative set of drug-resistant variants.

[Plasmid map and sequence file lot 100211](#)

This reagent is currently being provided as dried purified DNA stabilized in DNASTable Plus. Please see the notice for additional information and the protocol for reconstitution of dried DNA reagents. [Dried DNA Notice](#)

**Recommended
Storage:**

Keep the reagent at room temperature in a dry storage cabinet or in a moisture barrier bag.

Contributor:

Dr. Robert W. Shafer

References:

Johnston E., Dupnik K.M., Gonzales M.J., Winters M.A., Rhee S.Y., Imamichi T. and Shafer R.W. Panel of prototypical infectious molecular HIV-1 clones containing multiple nucleoside reverse transcriptase inhibitor resistance mutations. *AIDS* **19** (7):731-3, 2005. [Abstract](#)

NOTE:

Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: HIV-1 NL4-3 1499 Infectious Molecular Clone (p6463-13) from Dr. Robert W. Shafer (cat# 7401)." Also include the reference cited above in any publications.

Last Updated:

August 02, 2018

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.