SUPPORTING INFECTIOUS DISEASE RESEARCH

Escherichia coli K-12, Strain IM08B

Catalog No. NR-49806

Product Description:

Escherichia coli (*E. coli*) K-12, strain IM08B is an insertion mutant produced in *E. coli* K-12, strain DC10B via recombination-mediated genetic engineering (recombineering), and contains the *hsdM* (methylase) and *hsdS* (specificity) genes from *Staphylococcus aureus* (*S. aureus*), strain NRS384 clonal complex 8 (CC8). NR-49806 was produced by inoculation of BEI Resources seed lot 63849806 into Tryptic Soy broth with 25 µg per mL streptomycin and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar with 25 µg per mL streptomycin kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70039948

Manufacturing Date: 11NOV2020

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology	Report results	Circular, low convex, entire, smooth and cream (Figure 1)
Motility (wet mount)	Report results	Motile
VITEK [®] MS (MALDI-TOF)	E. coli	<i>E. coli</i> (99.9%)
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene (~ 1480 base pairs)	≥ 99% sequence identity to <i>E. coli</i> K-12 strain (GenBank: NZ_CP014225.1)	99.5% sequence identity to <i>E. coli</i> K-12 strain (GenBank: NZ_CP014225.1) ¹
Analysis of <i>hsdMS</i> by PCR Assay ^{2,3}		
hsdM2 and hsdS2 (CC8)	~ 3400 base pair amplicon	~ 3400 base pair amplicon
hsdS1 (CC8)	~ 1770 base pair amplicon	~ 1770 base pair amplicon
Analysis of <i>hsdMS</i> by Sequence Analysis ^{2,3}		
hsdM2 and hsdS2 (CC8)	Consistent with depositor sequence	Consistent with depositor sequence
hsdS1 (CC8)	Consistent with depositor sequence	Consistent with depositor sequence
Confirmation of Streptomycin Resistance 1 day at 37°C in an aerobic atmosphere on Tryptic Soy agar with 25 µg per mL streptomycin	Growth	Growth
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Also consistent with other *Escherichia* and *Shigella* species


²Performed on BEI Resources NRS-49806 lot 63849806

³PCR primers used for amplification were IM434 forward primer 5'-ACTTTCTTTAAGGCTTAGAGTCAAGC-3', IM435 reverse primer 5'-TTTAACGCCACGTTCACTCTTTGC-3', 179 forward primer 5'-CGGCCATTTATACAGGAAAAGCCTA-3' and 180 reverse primer 5'-GTTACCTTCTCATAGAGAGTGGTG-3'. For additional information, refer to Monk, I., et al. "Complete Bypass of Restriction Systems for Major *Staphylococcus aureus* Lineages." <u>mBio</u> 26 (2015): e00308-15. PubMed: 26015493. biei resources

Certificate of Analysis for NR-49806

SUPPORTING INFECTIOUS DISEASE RESEARCH

Figure 1: Colony Morphology



/Heather Couch/

Heather Couch

16 MAR 2021

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

 $ATCC^{\circledast}$ is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.

