

Certificate of Analysis for NR-52156

Acinetobacter baumannii, Strain MRSN 1187

Catalog No. NR-52156

This reagent is the tangible property of the U.S. Government.

Product Description:

Acinetobacter baumannii (A. baumannii), strain MRSN 1187 was isolated in 2010 from a human wound in the United States as part of a global surveillance program. A. baumannii, strain MRSN 1187 was deposited as sensitive to colistin, imipenem, meropenem and tetracycline and resistant to amikacin, ceftazidime, ciprofloxacin, ceftriaxone, gentamicin, trimethoprim/sulfamethoxazole, levofloxacin, tobramycin and ampicillin/sulbactam, with intermediate resistance to cefepime. NR-52156 lot 70039386 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70039386 Manufacturing Date: 24SEP2020

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

| TEST | SPECIFICATIONS | RESULTS |
|--|--|--|
| Phenotypic Analysis | | |
| Cellular morphology | Gram-negative rods | Gram-negative rods |
| Colony morphology | Report results | Circular, convex, entire, smooth and cream (Figure 1) |
| Growth at 44°C ± 2°C¹ 1 day in an aerobic atmosphere on Tryptic Soy agar | Growth | Growth |
| Motility Remel™ Motility Test Medium w/TTC Indicator for 1 day at 37°C in an aerobic atmosphere | Report results | Non-motile |
| VITEK® MS (MALDI-TOF) | A. baumannii | A. baumannii (99.9%) |
| Antibiotic Susceptibility Profile ^{2,3} | | |
| Amikacin | Resistant | Resistant (128 µg/mL) |
| Ampicillin/sulbactam | Resistant | Sensitive (8 µg/mL) ⁴ |
| Cefepime | Intermediate | Intermediate (12 to 16 µg/mL) |
| Ceftriaxone | Resistant | Resistant (≥ 64 µg/mL) |
| Ceftazidime | Resistant | Resistant (≥ 64 µg/mL) |
| Ciprofloxacin | Resistant | Resistant (≥ 4 µg/mL) |
| Colistin | Sensitive | Sensitive (≤ 0.25 µg/mL) |
| Gentamicin | Resistant | Resistant (≥ 16 µg/mL) |
| Imipenem | Sensitive | Sensitive (≤ 1 µg/mL) |
| Levofloxacin | Resistant | Resistant (> 4 µg/mL) |
| Meropenem | Sensitive | Sensitive to Intermediate (2 to 3 µg/mL) |
| Trimethoprim/sulfamethoxazole | Resistant | Resistant (> 4 µg/mL) |
| Tobramycin | Resistant | Resistant (64 µg/mL) |
| Tetracycline | Sensitive | Resistant (> 256 µg/mL) ⁵ |
| Genotypic Analysis | | , i j |
| Sequencing of 16S ribosomal RNA gene (1480 base pairs) | ≥ 99% sequence identity to A. baumannii, strain MRSN 1187 (GenBank: VHHC01000094.1) | 99.9% sequence identity to A. baumannii, strain MRSN 1187 (GenBank: VHHC01000094.1) |
| Purity 12 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar | Growth consistent with expected colony morphology | Growth consistent with expected colony morphology |

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898



Certificate of Analysis for NR-52156

| TEST | SPECIFICATIONS | RESULTS |
|-----------|----------------|---------|
| Viability | Growth | Growth |

¹Growth at 44°C differentiates A. baumannii from A. calcoaceticus and A. pittii, which do not grow at 44°C.

Figure 1: Colony Morphology



/Heather Couch/ Heather Couch

23 JUN 2021

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

²Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

³Antibiotic susceptibility was tested using a combination of VITEK®2 GN82, Sensititre™ GNX2F AST and E-test strips.

⁴A. baumannii, strain MRSN 1187 was deposited as resistant to ampicillin/sulbactam but showed a MIC of 8 μg/mL (interpreted as sensitive) for ampicillin/sulbactam during QC testing. Testing was performed in duplicate.

⁵A. baumannii, strain MRSN 1187 was deposited as sensitive to tetracycline but showed a MIC of > 256 μg/mL (interpreted as resistant) for tetracycline during QC testing. Testing was performed in quadruplicate.