

# **Certificate of Analysis for NR-52199**

### Acinetobacter baumannii, Strain MRSN 15129

### Catalog No. NR-52199

This reagent is the tangible property of the U.S. Government.

### **Product Description:**

Acinetobacter baumannii (A. baumannii), strain MRSN 15129 was isolated in 2012 from a human respiratory sample in South America as part of a global surveillance program. A. baumannii, strain MRSN 15129 was deposited as sensitive to amikacin, ceftazidime, colistin, imipenem, meropenem and ampicillin/sulbactam, intermediately resistant to ceftriaxone, cefepime and tobramycin and resistant to ciprofloxacin, gentamicin, trimethoprim/sulfamethoxazole, levofloxacin and tetracycline. NR-52199 lot 70040796 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70040796 Manufacturing Date: 23DEC2020

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: <a href="mailto:Contact@BEIResources.org">Contact@BEIResources.org</a>. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology	Report results	Circular, convex, entire, smooth and
Growth at 44°C ± 2°C <sup>1</sup>	Cusually	cream (Figure 1) Growth
1 day in an aerobic atmosphere on Tryptic	Growth	Growin
Soy agar		
Motility	Report results	Motile
Remel™ Motility Test Medium w/TTC	Troport rocale	Motifo
Indicator for 1 day at 37°C in an aerobic		
atmosphere		
VITEK® MS (MALDI-TOF)	A. baumannii	A. baumannii (99.9%)
Antibiotic Susceptibility Profile <sup>2,3</sup>		
Amikacin	Sensitive	Sensitive (16 µg/mL)
Ampicillin/sulbactam	Sensitive	Sensitive (1.5 μg/mL)
Cefepime	Intermediate	Intermediate (12 µg/mL)
Ceftriaxone	Intermediate	Intermediate (16 µg/mL)
Ceftazidime	Sensitive	Sensitive (8 µg/mL)
Ciprofloxacin	Resistant	Resistant (≥ 4 µg/mL)
Colistin	Sensitive	Sensitive (≤ 0.5 μg/mL)
Gentamicin	Resistant	Resistant (16 µg/mL)
Imipenem	Sensitive	Sensitive (≤ 1 μg/mL)
Levofloxacin	Resistant	Resistant (> 32 μg/mL)
Meropenem	Sensitive	Sensitive (≤ 1 μg/mL)
Trimethoprim/sulfamethoxazole	Resistant	Resistant (≥ 4 µg/mL)
Tobramycin	Intermediate	Sensitive (4 μg/mL) <sup>4</sup>
Tetracycline	Resistant	Resistant (> 256 μg/mL)
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene	≥ 99% sequence identity to	100% sequence identity to
(~ 1480 base pairs)	A. baumannii, strain MRSN 15129	A. baumannii, strain MRSN 15129
B ''	(GenBank: VHGR01000048.1)	(GenBank: VHGR01000048.1)
Purity  7 days at 27°C in an earthic atmosphere with	Growth consistent with expected colony	Growth consistent with expected colony
7 days at 37°C in an aerobic atmosphere with and without 5% CO <sub>2</sub> on Tryptic Soy agar	morphology	morphology
and without 370 002 on Tryphic day agai		

BEI Resources www.beiresources.org E-mail: <a href="mailto:contact@beiresources.org">contact@beiresources.org</a>
Tel: 800-359-7370

Fax: 703-365-2898



# **Certificate of Analysis for NR-52199**

TEST	SPECIFICATIONS	RESULTS
Viability	Growth	Growth

<sup>&</sup>lt;sup>1</sup>Growth at 44°C differentiates *A. baumannii* from *A. calcoaceticus and A. pittii*, which do not grow at 44°C.

Figure 1: Colony Morphology



/Heather Couch/ Heather Couch

26 JUL 2021

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898

<sup>&</sup>lt;sup>2</sup>Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

<sup>&</sup>lt;sup>3</sup>Antibiotic susceptibility was tested using a combination of VITEK®2 GN82, Sensititre™ GNX2F AST and E-test strips.

<sup>&</sup>lt;sup>4</sup>The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.