

Certificate of Analysis for NR-2634

Genomic DNA from Vaccinia Virus, Modified Vaccinia Ankara (MVA)

Catalog No. NR-2634

Product Description:


Genomic DNA was isolated from a preparation of cell lysate and supernatant from chicken embryo fibroblast cells infected with vaccinia virus, Modified Vaccinia Ankara (MVA)¹

Lot: 62153956 Manufacturing Date: 06NOV2013

TEST	SPECIFICATIONS	RESULTS
Sequencing of Vaccinia Virus Specific Region (989 base pairs)	Consistent with vaccinia virus, MVA	99% identity with vaccinia virus, MVA (GenBank: U94848)
Total DNA Content by PicoGreen® Measurement (Viral and Cellular)²	Report results	9.4 ng per 100 μL
Functional Activity by RT-PCR Amplification ³	~ 1100 base pair amplicon	~ 1100 bp amplicon (Figure 1)
Virus Inactivation 10% of total yield incubated on Vero cells and evaluated by hemagglutination assay ^{1,4,5}	No viable virus detected	No viable virus detected

¹SL-29 cells: ATCC[®] CRL-1590™

Figure 1: Functional Activity of NR-2634 by RT-PCR



/Heather Couch/ Heather Couch

06 OCT 2021

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898

²Nucleic acid was extracted from a preparation of Vaccinia Virus, MVA (BEI Resources NR-1 lot 7704329) using a QIAamp[®] Viral RNA Mini Kit (Qiagen 52906).

³Amplified using iTag™ DNA Polymerase (BioRad 170-8870) with 5 µL of NR-2634 in a 50 µL reaction

⁴Vero cells: ATCC[®] CCL-81[™]; 11 days at 37°C and 5% CO₂ with media overlay

⁵This extraction procedure has been demonstrated to consistently inactivate 100% of orthopoxviruses.