

Certificate of Analysis for NR-55305

Vector pCMV/R Containing the SARS-Related Coronavirus 2, Spike Glycoprotein Gene, Lineage B.1.351, Beta Variant

Catalog No. NR-55305

This reagent is the tangible property of the U.S. Government.

Product Description:

NR-55305 expresses the full-length, Beta variant spike (S) glycoprotein, and is intended for producing pseudotyped particles/pseudovirions. NR-55305 is not intended for recombinant protein expression. The vector for the S glycoprotein gene from severe acute respiratory syndrome-related coronavirus 2 (SARS-CoV-2), Wuhan-Hu-1 (GenBank: MN908947) was designed by codon optimizing the full-length S sequence (residues 1 to 1273) for mammalian expression and introducing point mutations found in the B.1.351 lineage, resulting in a spike glycoprotein gene representative of the Beta variant. The spike gene was subcloned into the pCMV/R mammalian expression vector (also referred to as VRC8400). The protein encoded by NR-55305 contains the following point mutations: L18F, D80A, D215G, LAL242-244del, R246l, K417N, E484K, N501Y, D614G and A701V. The kanamycin resistance gene, *aph*, provides transformant selection through kanamycin resistance in *Escherichia coli* (E. coli). NR-55305 is also referred to as VRC7597. The deposited plasmid was transformed into One Shot™ TOP10 E. coli (Invitrogen™ C404003), grown in Luria-Bertani broth with kanamycin (50 µg per mL) for 1 day at 37°C in an aerobic atmosphere, extracted using a Plasmid *Plus* Maxi Kit (QIAGEN® 12963) and vialed in TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0).

Lot: 70046965 Manufacturing Date: 07SEP2021

TEST	SPECIFICATIONS	RESULTS
Next-Generation DNA Sequencing	~ 8240 base pairs	8240 base pairs ¹
Genotypic Analysis Sequencing of S glycoprotein insert (~ 3820 base pairs) Lineage B.1.351 mutations	≥ 99% sequence identity to depositor's sequence Mutations present	100% sequence identity to depositor's sequence ² Mutations present
Antibiotic Resistance Kanamycin (encoded by kanamycin gene aph)	aph sequence present	aph sequence present
Concentration by Qubit Fluorometer®	≥ 2 µg per mL	1.7 μg in 90 μL per vial (19.4 μg per mL)
Amount per Vial	Report results	1.7 μg per vial
OD ₂₆₀ /OD ₂₈₀ Ratio	1.7 to 2.1	2.1
Effective Bacterial Transformation Invitrogen™ One Shot™ TOP10 <i>E. coli</i>	≥ 50 colonies per ng	189 colonies per ng

¹The sequence was assembled pre-vial using the depositor's provided sequence as the reference sequence. The complete plasmid sequence and map are provided on the BEI Resources webpage.

/Heather Couch/ Heather Couch

04 NOV 2021

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by ATCC® and the contributor to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org Tel: 800-359-7370

Fax: 703-365-2898

²The NR-55305 insert was codon optimized and B.1.351 mutations were introduced, but otherwise is 100% identical to the SARS-CoV-2, Wuhan-Hu-1 S protein (GenPept: YP 009724390.1).