

Nebraska Transposon Mutant Library (NTML) Screening Array

Catalog No. NR-48501

For research use only. Not for use in humans.

Contributor:

Kenneth Bayles, Ph.D., Director, Center for Staphylococcal Research, Department of Pathology and Microbiology and Paul Fey, Ph.D., Medical Director, Nebraska Medical Center Microbiology Section, University of Nebraska Medical Center, Omaha, Nebraska, USA, Network on Antimicrobial Resistance in *Staphylococcus aureus* (NARSA), and University of Chicago, Chicago, Illinois, USA (contribution of transposon, *bursa aurealis*)

Manufacturer:

University of Nebraska Medical Center

Product Description:

Production in the 384-well format has increased risk of cross-contamination between adjacent wells. Individual clones should be purified (e.g. single colony isolation and purification using good microbiological practices) and sequence-verified prior to use. BEI Resources does not confirm or validate individual clone identities provided by the contributor.

The Nebraska Transposon Mutant Library (NTML) Screening Array consists of approximately 2000 *Staphylococcus aureus* (*S. aureus*) subsp. *aureus* USA300 JE2, transposon (Tn) mutants arrayed in five 384-well microtiter plates. The screening array was developed to be a resource for performing high-throughput phenotypic screens to identify candidate genes for future research. The mutants were derived from *S. aureus* subsp. *aureus*, strain JE2. Mutagenesis occurred through the use of the *mariner*-based transposon *bursa aurealis* resulting in erythromycin-resistant deletion strains of JE2.^{1,2} The insertion sites were identified by sequencing the junction fragments containing the end of the transposon and the flanking DNA. The gene names and descriptions associated with each of the Tn mutants were obtained from the National Center for Biotechnology Information.²

S. aureus subsp. *aureus*, strain JE2 is a plasmid-cured derivative of strain LAC that was isolated in 2002 from a skin and soft tissue infection of an inmate in the Los Angeles County Jail in California, USA.² Strain JE2 is a methicillin-resistant *S. aureus* (MRSA) strain and is a USA300 isolate. USA300 isolates have the same MLST profile (ST 8), SCCmec (subtype IV), *agr* group (I) and *spa* motif (MBQBLO) and typically carry the Pantone-Valentine leukocidin (PVL) toxin genes and the arginine catabolic mobile element (ACME).^{3,4} USA300 is the most common cause of community-associated MRSA infection and an increasing cause of hospital-acquired infections.⁴

This screening array is intended for high-throughput phenotypic screens. The mutant strains have not had their Tn-insertion site re-sequenced to confirm the identity of each mutation; therefore, it is recommended that array users confirm each mutant selected from a screen. To ensure the quality of this resource for the research community and to reduce the number of passages from user to user, this library is not to be reproduced and distributed to other investigators in any way.

The Nebraska Transposon Mutant Library (NTML) was constructed in the laboratories of Dr. Kenneth Bayles and Dr. Paul Fey at the University of Nebraska Medical Center. Additional information is available at the [NTML](#) website.

Material Provided:

Each well contains approximately 150 µL of bacterial culture in Tryptic Soy broth containing 5 µg per mL erythromycin supplemented with 25% glycerol.

Packaging/Storage:

NR-48501 was packaged in 384-well microtiter plates. The product is provided frozen and should be stored at -60°C or colder immediately upon arrival. For long-term storage, the vapor phase of a liquid nitrogen freezer is recommended. Freeze-thaw cycles should be avoided.

Growth Conditions:

Media:

Tryptic Soy broth containing 5 µg per mL erythromycin or equivalent

Tryptic Soy agar containing 5 µg per mL erythromycin or equivalent

Incubation:

Temperature: 37°C

Atmosphere: Aerobic

Propagation:

1. Scrape top of frozen well with a pipette tip and streak onto agar plate.
2. Incubate the plates at 37°C for 1 day.

Citation:

Acknowledgment for publications should read "The following reagent was provided by the Network on Antimicrobial Resistance in *Staphylococcus aureus* (NARSA) for distribution by BEI Resources, NIAID, NIH: Nebraska Transposon Mutant Library (NTML) Screening Array, NR-48501."

Biosafety Level: 2

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. [Biosafety in Microbiological and Biomedical Laboratories](#). 6th ed. Washington, DC: U.S. Government Printing Office, 2020; see www.cdc.gov/biosafety/publications/bmb15/index.htm.

Disclaimers:

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at www.beiresources.org.

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC® nor the U.S. Government makes any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC® nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. ATCC® and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, ATCC®, their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

Use Restrictions:

This material is distributed for internal research, non-commercial purposes only. This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale.

References:

1. Bae, T., et al. "Staphylococcus aureus Virulence Genes Identified by bursa aurealis Mutagenesis and Nematode Killing." Proc. Natl. Acad. Sci. USA 101 (2004): 12312-12317. PubMed: 15304642.
2. Fey, P. D., et al. "A Genetic Resource for Rapid and Comprehensive Phenotype Screening of Nonessential Staphylococcus aureus Genes." mBio. 4 (2013): e00537-12. PubMed: 23404398.
3. Diep, B. A., et al. "Roles of 34 Virulence Genes in the Evolution of Hospital- and Community-Associated Strains of Methicillin-Resistant Staphylococcus aureus." J. Infect. Dis. 193 (2006): 1495-1503. PubMed: 16652276.
4. Diekema, D. J., et al. "Continued Emergence of USA300 Methicillin-Resistant Staphylococcus aureus in the United States: Results from a Nationwide Surveillance Study." Infect. Control Hosp. Epidemiol. 35 (2014): 285-292. PubMed: 24521595.

ATCC® is a trademark of the American Type Culture Collection.

