

Klebsiella pneumoniae* MRSN Diversity Panel*Catalog No. NR-55604**

This reagent is the tangible property of the U.S. Government.

Product Description:

The *Klebsiella pneumoniae* (*K. pneumoniae*) Multidrug-Resistant Organism Repository and Surveillance Network (MRSN) strains that comprise NR-55604 were isolated between 2003 and 2020 as part of a surveillance program in the United States.

Lot: 70052443**Manufacturing Date: 2021 and 2022**

QC testing was performed, and the results are provided on the Certificate of Analysis for each isolate.

/Sonia Bjorum Brower/

Sonia Bjorum Brower

15 AUG 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Heather Couch/

Heather Couch

22 MAR 2022

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 4111**

Catalog No. NR-55505

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 4111 was isolated in 2011 from a human perianal sample in North America as part of a global surveillance program. NR-55505 was deposited as a susceptible strain, sensitive to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55505 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70049674

Manufacturing Date: 19JAN2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (87%) ¹
Antibiotic Susceptibility Profile^{2,3} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive	Sensitive (≤ 2 µg/mL) Sensitive (≤ 2 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.38 µg/mL) Sensitive (0.25 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.032 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.25 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 0.5 µg/mL) ⁴ Sensitive (≤ 1 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1480 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 4111 (GenBank: JAGYFB010000168.1)	99.6% sequence identity to <i>K. pneumoniae</i> , strain MRSN 4111 (GenBank: JAGYFB010000168.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Although the result is out of specification, performance of the product should not be affected.


²Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

³Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

09 AUG 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 4759**

Catalog No. NR-55506

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 4759 was isolated in 2011 from a human urine sample in North America as part of a global surveillance program. NR-55506 was deposited as a multidrug-resistant strain, sensitive to amikacin, cefepime, ceftazidime/avibactam, ceftolozane/tazobactam, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline, tobramycin and trimethoprim/sulfamethoxazole and resistant to ampicillin/sulbactam, aztreonam, ceftazidime and ceftriaxone. NR-55506 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70049672

Manufacturing Date: 20JAN2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card) VITEK® MS (MALDI-TOF)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%) <i>K. pneumoniae</i>	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%) <i>K. pneumoniae</i> (99.9%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Sensitive Resistant Sensitive Sensitive Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Sensitive (2 µg/mL) ³ Sensitive (2 µg/mL) Sensitive (1 µg/mL) ⁴ Sensitive (0.19 µg/mL) Sensitive (0.25 to 0.38 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.25 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 0.5 µg/mL) ⁵ Sensitive (≤ 1 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 4759 (GenBank: JAGYFA010000074.1)	99.3% sequence identity to <i>K. pneumoniae</i> , strain MRSN 4759 (GenBank: JAGYFA010000074.1) ⁶

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].


³*K. pneumoniae*, strain MRSN 4759 was deposited as resistant to aztreonam, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴*K. pneumoniae*, strain MRSN 4759 was deposited as resistant to ceftazidime, but showed a MIC of 1 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Heather Couch/

Heather Couch

18 APR 2022

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 4815**

Catalog No. NR-55507

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 4815 was isolated in 2011 from a human wound sample in North America as part of a global surveillance program. NR-55507 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, imipenem, meropenem, piperacillin/tazobactam, tobramycin, intermediately resistant to cefepime and resistant to ampicillin/sulbactam, aztreonam, ceftazidime, ceftriaxone, ciprofloxacin, gentamicin, levofloxacin, tetracycline, tigecycline and trimethoprim/sulfamethoxazole. NR-55507 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70049670

Manufacturing Date: 20JAN2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card) VITEK® MS (MALDI-TOF)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%) <i>K. pneumoniae</i>	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%) <i>K. pneumoniae</i> (99.9%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Intermediate Sensitive Sensitive Sensitive Sensitive Resistant Resistant Sensitive Resistant Sensitive Resistant Sensitive Resistant Resistant Sensitive Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Sensitive (2 to 4 µg/mL) ³ Intermediate (3 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.75 µg/mL) Sensitive (0.75 µg/mL) Resistant (32 µg/mL) Resistant (6 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.125 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (16 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 8 µg/mL) ⁴ Sensitive (4 µg/mL) Resistant (160 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 4815 (GenBank: JAGYEZ010000093.1)	99.7% sequence identity to <i>K. pneumoniae</i> , strain MRSN 4815 (GenBank: JAGYEZ010000093.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 4815 was deposited as resistant to aztreonam, but showed a MIC of 2 µg per mL and 4 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Heather Couch/

Heather Couch

Program Manager or designee, ATCC Federal Solutions

15 APR 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 5613**

Catalog No. NR-55508

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 5613 was isolated in 2011 from a human urine sample in North America as part of a global surveillance program. NR-55508 was deposited as an extensively drug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ertapenem, imipenem, meropenem and tigecycline, intermediately resistant to ciprofloxacin, levofloxacin and piperacillin/tazobactam and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, gentamicin, tetracycline, trimethoprim/sulfamethoxazole and tobramycin. NR-55508 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70049668

Manufacturing Date: 06JAN2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® MS (MALDI-TOF)	Gram-negative rods Report results Report results <i>K. pneumoniae</i>	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99.9%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Intermediate Sensitive Resistant Sensitive Intermediate Sensitive Intermediate Resistant Sensitive Resistant Resistant	Sensitive ($\leq 2 \mu\text{g/mL}$) Resistant ($\geq 32 \mu\text{g/mL}$) Resistant (32 $\mu\text{g/mL}$) Sensitive (2 $\mu\text{g/mL}$)³ Resistant (16 $\mu\text{g/mL}$) Sensitive (0.75 $\mu\text{g/mL}$) Sensitive (1 $\mu\text{g/mL}$)⁴ Resistant ($\geq 64 \mu\text{g/mL}$) Resistant (4 $\mu\text{g/mL}$) ⁵ Sensitive ($\leq 0.5 \mu\text{g/mL}$) Resistant ($\geq 16 \mu\text{g/mL}$) Sensitive (0.25 $\mu\text{g/mL}$) Sensitive (1.5 to 2 $\mu\text{g/mL}$)⁶ Sensitive ($\leq 0.25 \mu\text{g/mL}$) Intermediate (24 $\mu\text{g/mL}$) Resistant ($\geq 16 \mu\text{g/mL}$) Resistant (2 $\mu\text{g/mL}$) ^{5,7} Resistant ($\geq 16 \mu\text{g/mL}$) Resistant ($\geq 320 \mu\text{g/mL}$)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	$\geq 99\%$ sequence identity to <i>K. pneumoniae</i> , strain MRSN 5613 (GenBank: JAGY EY010000099.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 5613 (GenBank: JAGY EY010000099.1) ⁸

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 5613 was deposited as resistant to cefepime, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴*K. pneumoniae*, strain MRSN 5613 was deposited as resistant to ceftolozane/tazobactam, but showed a MIC of 1 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed twice in duplicate.


⁵The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁶*K. pneumoniae*, strain MRSN 5613 was deposited as intermediately resistant to levofloxacin, but showed a MIC of 1.5 to 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed twice in duplicate.

⁷MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁸Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

10 AUG 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 5741**

Catalog No. NR-55509

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 5741 was isolated in 2003 from a human respiratory sample in North America as part of a global surveillance program. NR-55509 was deposited as a multidrug-resistant strain, sensitive to amikacin, cefepime, ceftazidime/avibactam, ertapenem, imipenem, meropenem, tetracycline and tigecycline, intermediately resistant to gentamicin and tobramycin and resistant to ampicillin/sulbactam, aztreonam, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, levofloxacin, piperacillin/tazobactam and trimethoprim/sulfamethoxazole. NR-55509 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70049666

Manufacturing Date: 06JAN2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Sensitive Resistant Sensitive Resistant Resistant Resistant Sensitive Intermediate Sensitive Resistant Sensitive Resistant Sensitive Sensitive Intermediate Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Sensitive (4 µg/mL) ³ Sensitive (≤ 1 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (6 µg/mL) Inconclusive ⁴ Resistant (8 to 12 µg/mL) Resistant (≥ 32 µg/mL) Sensitive (≤ 0.5 µg/mL) Intermediate (8 µg/mL) Sensitive (0.19 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 0.25 µg/mL) Resistant (≥ 128 µg/mL) Sensitive (4 µg/mL) Sensitive (0.38 µg/mL) ⁵ Intermediate (8 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1480 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 5741 (GenBank: JAGYEX01000090.1)	99.8% sequence identity to <i>K. pneumoniae</i> , strain MRSN 5741 (GenBank: JAGYEX01000090.1) ⁶

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].


³*K. pneumoniae*, strain MRSN 5741 was deposited as resistant to aztreonam, but showed a MIC of 4 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴*K. pneumoniae*, strain MRSN 5741 was deposited as resistant to ceftolozane/tazobactam. Antibiotic susceptibility testing performed in quadruplicate determined that for strain MRSN 5741, the ceftolozane/tazobactam MICs are 3 and 4 µg per mL, which are interpreted as sensitive and intermediately resistant, respectively.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

21 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 5881**

Catalog No. NR-55510

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 5881 was isolated in 2005 from a human wound sample in North America as part of a global surveillance program. NR-55510 was deposited as an extensively drug-resistant strain, sensitive to ceftazidime/avibactam, ertapenem, imipenem, meropenem and tigecycline, intermediately resistant to levofloxacin and resistant to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, gentamicin, piperacillin/tazobactam, tetracycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55510 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70049663

Manufacturing Date: 20JAN2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Sensitive Resistant Sensitive Intermediate Sensitive Resistant Resistant Resistant Sensitive Resistant Resistant	Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Intermediate (4 µg/mL) ³ Resistant (16 to 32 µg/mL) Sensitive (0.5 to 0.75 µg/mL) Resistant (16 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.19 to 0.25 µg/mL) Intermediate (4 µg/mL) Sensitive (≤ 0.25 µg/mL) Resistant (≥ 128 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (1 µg/mL) ⁴ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 5881 (GenBank: JAGYEW010000112.1)	99.3% sequence identity to <i>K. pneumoniae</i> , strain MRSN 5881 (GenBank: JAGYEW010000112.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)


²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 5881 was deposited as resistant to cefepime, but showed a MIC of 4 µg/mL (interpreted as intermediately resistant) for cefepime during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

27 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 6031**

Catalog No. NR-55511

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 6031 was isolated in 2011 from a human wound sample in Europe as part of a global surveillance program. NR-55511 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam and tigecycline, intermediately resistant to aztreonam and ciprofloxacin and resistant to ampicillin/sulbactam, cefepime, ceftazidime, ceftriaxone, tetracycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55511 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70049661

Manufacturing Date: 19JAN2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Intermediate Resistant Resistant Sensitive Sensitive Resistant Intermediate Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Resistant Sensitive Resistant Resistant	Sensitive (4 µg/mL) Resistant (≥ 32 µg/mL) Resistant (16 µg/mL) ³ Sensitive (2 µg/mL) ⁴ Sensitive (4 µg/mL) ⁵ Sensitive (0.125 µg/mL) Sensitive (0.25 to 0.75 µg/mL) Resistant (≥ 64 µg/mL) Resistant (3 µg/mL) ³ Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.25 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (2 µg/mL) ^{3,6} Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1480 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 6031 (GenBank: JAGYEV010000085.1)	99.6% sequence identity to <i>K. pneumoniae</i> , strain MRSN 6031 (GenBank: JAGYEV010000085.1) ⁷

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.


⁴*K. pneumoniae*, strain MRSN 6031 was deposited as resistant to cefepime, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵*K. pneumoniae*, strain MRSN 6031 was deposited as resistant to ceftazidime, but showed a MIC of 4 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁶MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁷Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

08 AUG 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 6778**

Catalog No. NR-55512

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 6778 was isolated in 2011 from a human urine sample in North America as part of a global surveillance program. NR-55512 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ertapenem, imipenem, meropenem, tetracycline and tigecycline, intermediately resistant to levofloxacin and piperacillin/tazobactam and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, ceftolozane/tazobactam, ciprofloxacin, gentamicin, tobramycin and trimethoprim/sulfamethoxazole. NR-55512 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70049659

Manufacturing Date: 14JAN2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Sensitive Resistant Sensitive Sensitive Intermediate Sensitive Intermediate Sensitive Sensitive Resistant Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Sensitive (0.064 to 0.19 µg/mL) ³ Sensitive (0.25 µg/mL) ⁴ Sensitive (0.19 to 0.38 µg/mL) ⁵ Sensitive (0.25 µg/mL) Sensitive (0.25 µg/mL) ⁶ Sensitive (0.094 to 0.125 µg/mL) ⁷ Resistant (4 to 6 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (0.5 µg/mL) ⁸ Sensitive (0.25 µg/mL) Sensitive (0.75 to 1 µg/mL) ⁹ Sensitive (≤ 0.25 µg/mL) Intermediate (24 to 34 µg/mL) ¹⁰ Resistant (≥ 256 µg/mL) ¹¹ Sensitive (1 µg/mL) ¹² Intermediate (6 µg/mL) ¹³ Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 6778 (GenBank: JAGYEU010000114.1)	99.4% sequence identity to <i>K. pneumoniae</i> , strain MRSN 6778 (GenBank: JAGYEU010000114.1) ¹⁴

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 6778 was deposited as resistant to aztreonam, but showed a MIC of 0.064 µg per mL to 0.19 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴*K. pneumoniae*, strain MRSN 6778 was deposited as resistant to cefepime, but showed a MIC of 0.25 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵*K. pneumoniae*, strain MRSN 6778 was deposited as resistant to ceftazidime, but showed a MIC of 0.19 µg per mL to 0.38 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁶*K. pneumoniae*, strain MRSN 6778 was deposited as resistant to ceftolozane/tazobactam, but showed a MIC of 0.25 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁷*K. pneumoniae*, strain MRSN 6778 was deposited as resistant to ceftriaxone, but showed a MIC of 0.094 µg per mL to 0.125 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁸*K. pneumoniae*, strain MRSN 6778 was deposited as resistant to gentamicin, but showed a MIC of 0.25 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁹*K. pneumoniae*, strain MRSN 6778 was deposited as resistant to levofloxacin, but showed a MIC of 0.75 µg per mL to 1 µg per mL (interpreted as intermediately resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

¹⁰The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.


¹¹*K. pneumoniae*, strain MRSN 6778 was deposited as sensitive to tetracycline, but showed a MIC of ≥ 256 µg per mL (interpreted as resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

¹²MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

¹³*K. pneumoniae*, strain MRSN 6778 was deposited as resistant to tobramycin, but showed a MIC of 6 µg per mL (interpreted as intermediately resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

¹⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

04 AUG 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 7076**

Catalog No. NR-55513

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 7076 was isolated in 2003 from a human wound sample in North America as part of a global surveillance program. NR-55513 was deposited as an extensively drug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ciprofloxacin, ertapenem, imipenem, levofloxacin and meropenem, intermediately resistant to piperacillin/tazobactam and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, gentamicin, tetracycline, trimethoprim/sulfamethoxazole, tigecycline and tobramycin. NR-55513 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70049657

Manufacturing Date: 19JAN2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® MS (MALDI-TOF)	Gram-negative rods Report results Report results <i>K. pneumoniae</i>	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Intermediate Resistant Resistant Resistant Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Intermediate (8 µg/mL) ³ Resistant (16 µg/mL) Sensitive (1.5 µg/mL) Sensitive (1.5 µg/mL)⁴ Resistant (≥ 64 µg/mL) Sensitive (0.5 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.19 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.25 µg/mL) Intermediate (32 to 48 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 8 µg/mL) ⁵ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 7076 (GenBank: JAGYET010000100.1)	99.7% sequence identity to <i>K. pneumoniae</i> , strain MRSN 7076 (GenBank: JAGYET010000100.1) ⁶

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK®2 GN74 and ETEST®.


³The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁴*K. pneumoniae*, strain MRSN 7076 was deposited as resistant to ceftolozane/tazobactam, but showed a MIC of 1.5 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

01 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

14 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 13748**

Catalog No. NR-55515

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 13748 was isolated in 2011 from a human blood sample in Europe as part of a global surveillance program. NR-55515 was deposited as an extensively drug-resistant strain, sensitive to amikacin and ceftazidime/avibactam, intermediately resistant to ceftazidime and tobramycin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline and trimethoprim/sulfamethoxazole. NR-55515 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70049653

Manufacturing Date: 14JAN2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Intermediate Sensitive Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Intermediate Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (16 µg/mL) ³ Sensitive (1.5 µg/mL) Resistant (16 to 24 µg/mL) Resistant (≥ 64 µg/mL) Intermediate (2 to 3 µg/mL) ³ Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.5 to 0.75 µg/mL) ⁴ Resistant (≥ 8 µg/mL) Intermediate (2 µg/mL) ³ Resistant (≥ 128 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 8 µg/mL) ⁵ Intermediate (8 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1490 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 13748 (GenBank: JAGYER010000099.1)	99.7% sequence identity to <i>K. pneumoniae</i> , strain MRSN 13748 (GenBank: JAGYER010000099.1) ⁶

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK®2 GN74 and ETEST®.

³The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁴*K. pneumoniae*, strain MRSN 13748 was deposited as resistant to imipenem, but showed a MIC of 0.5 µg per mL and 0.75 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

08 AUG 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 13761 was deposited as intermediately resistant to aztreonam, but showed a MIC of 24 to 32 µg per mL (interpreted as resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴*K. pneumoniae*, strain MRSN 13761 was deposited as resistant to cefepime, but showed a MIC of 4 µg per mL (interpreted as intermediately resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁶*K. pneumoniae*, strain MRSN 13761 was deposited as intermediately resistant to levofloxacin, but showed a MIC of 1.5 to 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁷MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁸Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

27 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 13768**

Catalog No. NR-55517

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 13768 was isolated in 2011 from a human blood sample in Europe as part of a global surveillance program. NR-55517 was deposited as an extensively drug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, gentamicin, tigecycline and trimethoprim/sulfamethoxazole and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline and tobramycin. NR-55517 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70049649

Manufacturing Date: 06JAN2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® MS (MALDI-TOF)	Gram-negative rods Report results Report results <i>K. pneumoniae</i>	Gram-negative rods Circular, convex, entire, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99.9%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Sensitive	Sensitive (8 to 16 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (4 µg/mL) Resistant (> 256 µg/mL) Resistant (≥ 64 µg/mL) Resistant (> 32 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 1 µg/mL) Resistant (> 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Resistant (≥ 16 µg/mL) Inconclusive ^{3,4} Resistant (≥ 16 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 13768 (GenBank: JAGYEP010000113.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 13768 (GenBank: JAGYEP010000113.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)


²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK®2 GN74 and ETEST®.

³*K. pneumoniae*, strain MRSN 13768 was deposited as being sensitive to tigecycline. Antibiotic susceptibility testing performed in duplicate determined that for strain MRSN 13768, the tigecycline MICs are 1 µg per mL and 2 µg per mL, which are interpreted as sensitive and resistant, respectively.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

11 AUG 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

/Sonia Bjorum Brower/
Sonia Bjorum Brower

29 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 15219**

Catalog No. NR-55519

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 15219 was isolated in 2013 from a human urine sample in North America as part of a global surveillance program. NR-55519 was deposited as a multidrug-resistant strain, sensitive to amikacin, cefepime, ceftazidime/avibactam, ceftolozane/tazobactam, ciprofloxacin, ertapenem, imipenem, levofloxacin, meropenem, tigecycline and piperacillin/tazobactam and resistant to aztreonam, ampicillin/sulbactam, ceftazidime, ceftriaxone, gentamicin, tetracycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55519 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70049644

Manufacturing Date: 12JAN2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Sensitive Resistant Sensitive Sensitive Sensitive Resistant Sensitive Sensitive Resistant Sensitive Sensitive Sensitive Sensitive Resistant Sensitive Resistant Resistant	Sensitive (16 µg/mL) Sensitive (6 µg/mL) ³ Resistant (≥ 64 µg/mL) Sensitive (0.75 to 1.0 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.5 µg/mL) Sensitive (0.38 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.5 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.19 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (8 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (≤ 0.5 µg/mL) ⁴ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 15219 (GenBank: JAGYEN010000159.1)	99.6% sequence identity to <i>K. pneumoniae</i> , strain MRSN 15219 (GenBank: JAGYEN010000159.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 15219 was deposited as resistant to ampicillin/sulbactam, but showed a MIC of 6 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed twice in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

04 AUG 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 15687**

Catalog No. NR-55520

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 15687 was isolated in 2013 from a human urine sample in North America as part of a global surveillance program. NR-55520 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, trimethoprim/sulfamethoxazole, tetracycline and tigecycline, intermediately resistant to piperacillin/tazobactam and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone and tobramycin. NR-55520 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70049643

Manufacturing Date: 26JAN2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, mucoid and cream Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Intermediate Sensitive Sensitive Sensitive Resistant Sensitive	Sensitive (4 µg/mL) Resistant (≥ 32 µg/mL) Resistant (32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (16 µg/mL) Sensitive (0.5 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.125 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.25 to 0.5 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Intermediate (48 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 0.5 µg/mL) ³ Resistant (≥ 16 µg/mL) Resistant (> 32 µg/mL) ⁴
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1480 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 15687 (GenBank: JAGYEM010000080.1)	99.6% sequence identity to <i>K. pneumoniae</i> , strain MRSN 15687 (GenBank: JAGYEM010000080.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴*K. pneumoniae*, strain MRSN 15687 was deposited as sensitive to trimethoprim/sulfamethoxazole, but showed a MIC of > 32 µg per mL (interpreted as resistant) for this antibiotic during QC testing. Testing was performed twice in duplicate.

⁵Also consistent with other *Klebsiella* species

/Sonia Bjorum Brower/
Sonia Bjorum Brower

27 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 15882 was deposited as resistant to ampicillin/sulbactam, but showed a MIC of 8 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵*K. pneumoniae*, strain MRSN 15882 was deposited as sensitive to tigecycline, but showed a MIC of 4 µg per mL (interpreted as resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

Figure 1: Colony Morphology



/Heather Couch/
Heather Couch

28 APR 2022

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

09 JUN 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

28 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁴*K. pneumoniae*, strain MRSN 16233 was deposited as resistant to piperacillin/tazobactam, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶Also consistent with other *Klebsiella* species

/Sonia Bjorum Brower/

Sonia Bjorum Brower

09 AUG 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth


¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

29 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 8 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK®2 GN74 and ETEST®.


³*K. pneumoniae*, strain MRSN 19073 was deposited as resistant to cefepime, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorun Brower/

Sonia Bjorun Brower

08 AUG 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 20522**

Catalog No. NR-55527

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 20522 was isolated in 2013 from a human respiratory sample in North America as part of a global surveillance program. NR-55527 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ertapenem, imipenem, levofloxacin, meropenem, tigecycline and trimethoprim/sulfamethoxazole, intermediately resistant to ciprofloxacin and piperacillin/tazobactam and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, gentamicin, tetracycline and tobramycin. NR-55527 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050350

Manufacturing Date: 10FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Intermediate Sensitive Resistant Sensitive Sensitive Sensitive Intermediate Resistant Sensitive Resistant Sensitive	Sensitive (4 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (8 to 16 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.5 µg/mL) Sensitive (1 µg/mL) ³ Resistant (≥ 64 µg/mL) Sensitive (1 µg/mL) ⁴ Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.25 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.25 µg/mL) Intermediate (32 µg/mL) Resistant (≥ 16 µg/mL) Resistant (2 µg/mL) ^{4,5} Resistant (≥ 16 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 20522 (GenBank: JAGYEF010000101.1)	99.4% sequence identity to <i>K. pneumoniae</i> , strain MRSN 20522 (GenBank: JAGYEF010000101.1) ⁶

TEST	SPECIFICATIONS	RESULTS
Purity 8 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 20522 was deposited as resistant to ceftolozane/tazobactam, but showed a MIC of 1 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

08 AUG 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 8 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)


²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 21304 was deposited as being sensitive to piperacillin/tazobactam. Antibiotic susceptibility testing performed in duplicate determined that for strain MRSN 21304, the piperacillin/tazobactam MICs are 12 µg per mL and 24 µg per mL, which are interpreted as sensitive and resistant, respectively.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

08 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 21352**

Catalog No. NR-55529

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 21352 was isolated in 2014 from a human urine sample in North America as part of a global surveillance program. NR-55529 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, gentamicin, imipenem, meropenem, tetracycline and tigecycline, intermediately resistant to ceftazidime, piperacillin/tazobactam and tobramycin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftriaxone, ciprofloxacin, levofloxacin and trimethoprim/sulfamethoxazole. NR-55529 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050354

Manufacturing Date: 17FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Intermediate Sensitive Sensitive Sensitive Resistant Resistant Sensitive Sensitive Sensitive Resistant Sensitive Intermediate Sensitive Sensitive Intermediate Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (2 µg/mL) ³ Resistant (16 µg/mL) ⁴ Sensitive (0.5 µg/mL) Sensitive (0.75 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.125 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 0.25 µg/mL) Intermediate (32 µg/mL) Sensitive (2 µg/mL) Sensitive (1 µg/mL) ⁵ Intermediate (8 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 21352 (GenBank: JAGYED020000062.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 21352 (GenBank: JAGYED020000062.1) ⁶

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 21352 was deposited as resistant to cefepime, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

29 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 22232**

Catalog No. NR-55530

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 22232 was isolated in 2013 from a human respiratory sample in South America as part of a global surveillance program. NR-55530 was deposited as an extensively drug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ertapenem, imipenem, levofloxacin, meropenem and tigecycline, intermediately resistant to ciprofloxacin and piperacillin/tazobactam, and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, gentamicin, tetracycline, trimethoprim/sulfamethoxazole and tobramycin. NR-55530 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050356

Manufacturing Date: 17FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, mucoid, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Intermediate Sensitive Resistant Sensitive Sensitive Sensitive Intermediate Resistant Sensitive Resistant Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (16 to 64 µg/mL) Resistant (16 µg/mL) Sensitive (0.25 µg/mL) Sensitive (0.38 µg/mL)³ Resistant (≥ 64 µg/mL) Intermediate (2 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.25 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.25 µg/mL) Intermediate (64 to 96 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (1 µg/mL) ⁴ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 22232 (GenBank: JAGYEC01000097.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 22232 (GenBank: JAGYEC01000097.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 22232 was deposited as resistant to ceftolozane/tazobactam, but showed a MIC of 0.38 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

21 JUN 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 22265**

Catalog No. NR-55531

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 22265 was isolated in 2012 from a human respiratory sample in South America as part of a global surveillance program. NR-55531 was deposited as a multidrug-resistant strain, sensitive to amikacin, cefepime, ceftazidime/avibactam, ertapenem, imipenem, meropenem, tetracycline and tigecycline, intermediately resistant to piperacillin/tazobactam and resistant to ampicillin/sulbactam, aztreonam, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, gentamicin, levofloxacin, trimethoprim/sulfamethoxazole and tobramycin. NR-55531 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050358

Manufacturing Date: 18FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, mucoid, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Sensitive Resistant Sensitive Resistant Resistant Resistant Sensitive Resistant Sensitive Resistant Sensitive Intermediate Sensitive Sensitive Resistant Resistant	Sensitive (≤ 2 µg/mL) Sensitive (6 to 8 µg/mL) ³ Resistant (≥ 64 µg/mL) Sensitive (≤ 1.5 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.75 µg/mL) Sensitive (0.5 to 0.75 µg/mL) ⁴ Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.19 to 0.25 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 0.25 µg/mL) Intermediate (32 µg/mL) Sensitive (2 µg/mL) Sensitive (≤ 0.5 µg/mL) ⁵ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 22265 (GenBank: JAGYEB010000089.1)	99.4% sequence identity to <i>K. pneumoniae</i> , strain MRSN 22265 (GenBank: JAGYEB010000089.1)

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)


²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 22265 was deposited as resistant to ampicillin/sulbactam but showed a MIC of 6 to 8 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴*K. pneumoniae*, strain MRSN 22265 was deposited as resistant to ceftolozane/tazobactam but showed a MIC of 0.5 to 0.75 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

14 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth


¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

26 JUL 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 15219 was deposited as intermediately resistant to piperacillin/tazobactam, but showed a MIC of 3 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed twice in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

/Heather Couch/
Heather Couch

06 MAY 2022

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

/Sonia Bjorum Brower/

Sonia Bjorum Brower

16 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

/Heather Couch/

Heather Couch

25 APR 2022

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 27106**

Catalog No. NR-55536

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 27106 was isolated in 2014 from human urine sample in South America as part of a global surveillance program. NR-55536 was deposited as a susceptible strain, sensitive to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tigecycline, tobramycin, tetracycline and trimethoprim/sulfamethoxazole. NR-55536 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050369

Manufacturing Date: 16FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive	Sensitive (≤ 2 µg/mL) Sensitive (4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.19 µg/mL) Sensitive (0.19 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 0.016 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 0.5 µg/mL) ³ Sensitive (≤ 1 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1460 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 27106 (GenBank: JAGYDX010000088.1 and JAGYDX010000089.1)	99.9% sequence identity to <i>K. pneumoniae</i> , strain MRSN 27106 (GenBank: JAGYDX010000088.1 and JAGYDX010000089.1) ⁴

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

/Sonia Bjorum Brower/

Sonia Bjorum Brower

27 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 27778**

Catalog No. NR-55537

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 27778 was isolated in 2015 from a human sample in North America as part of a global surveillance program. NR-55537 was deposited as a multidrug-resistant strain, sensitive to ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, imipenem, meropenem, piperacillin/tazobactam and tigecycline, intermediately resistant to tetracycline and resistant to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, ciprofloxacin, gentamicin, levofloxacin, tobramycin and trimethoprim/sulfamethoxazole. NR-55537 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050371

Manufacturing Date: 11FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Resistant Resistant Resistant Resistant Resistant Sensitive Sensitive Resistant Resistant Sensitive Resistant Resistant Sensitive Sensitive Sensitive Intermediate Sensitive Resistant Resistant	Resistant (> 256 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (2 µg/mL) ³ Resistant (16 µg/mL) Sensitive (0.25 µg/mL) Sensitive (0.75 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.19 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (16 µg/mL) Intermediate (8 µg/mL) Resistant (≥ 2 µg/mL) ^{4,5} Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 27778 (GenBank: JAGYDW010000047.1)	99.6% sequence identity to <i>K. pneumoniae</i> , strain MRSN 27778 (GenBank: JAGYDW010000047.1) ⁶

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 27778 was deposited as resistant to cefepime, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁶Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Heather Couch/
Heather Couch

14 APR 2022

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 27989**

Catalog No. NR-55538

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 27989 was isolated in 2014 from a human wound sample in North America as part of a global surveillance program. NR-55538 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ertapenem, imipenem, meropenem, tigecycline and piperacillin/tazobactam, intermediately resistant to levofloxacin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, gentamicin, tetracycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55538 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050373

Manufacturing Date: 11FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Sensitive Resistant Sensitive Resistant Sensitive Resistant Sensitive Resistant Resistant Resistant Intermediate Sensitive Sensitive Resistant Sensitive Resistant Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.25 µg/mL) Sensitive (0.5 to 0.75 µg/mL) ³ Resistant (≥ 64 µg/mL) Resistant (4 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.19 to 0.25 µg/mL) Intermediate (4 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (16 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (≤ 0.5 µg/mL) ⁴ Resistant (≥ 256 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 27989 (GenBank: JAGYDV010000076.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 27989 (GenBank: JAGYDV010000076.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)


²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 27989 was deposited as resistant to ceftolozane/tazobactam, but showed a MIC of 0.5 µg per mL and 0.75 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

26 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 28183**

Catalog No. NR-55539

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 28183 was isolated in 2008 from a human respiratory sample in North America as part of a global surveillance program. NR-55539 was deposited as an extensively drug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ertapenem, gentamicin, imipenem, meropenem and tigecycline, intermediately resistant to piperacillin/tazobactam and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, levofloxacin, tetracycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55539 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050375

Manufacturing Date: 11FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Sensitive Sensitive Sensitive Sensitive Resistant Resistant Sensitive Intermediate Resistant Sensitive Resistant Resistant Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (≤ 1 µg/mL) ³ Resistant (≥ 64 µg/mL) Sensitive (0.25 µg/mL) Sensitive (0.19 µg/mL) ⁴ Resistant (8 µg/mL) Resistant (≥ 32 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.25 to 0.38 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 0.25 µg/mL) Intermediate (24 µg/mL) Resistant (≥ 16 µg/mL) Resistant (2 µg/mL) ^{5,6} Inconclusive (4 to 6 µg/mL) ⁷ Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 28183 (GenBank: JAGYDU010000101.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 28183 (GenBank: JAGYDU010000101.1) ⁸

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 28183 was deposited as resistant to cefepime, but showed a MIC of ≤ 1 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴*K. pneumoniae*, strain MRSN 28183 was deposited as resistant to ceftolozane/tazobactam, but showed a MIC of 0.19 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.


⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁷*K. pneumoniae*, strain MRSN 28183 was deposited as resistant to tobramycin, but showed a MIC of 4 µg per mL and 6 µg per mL (interpreted as sensitive and intermediate, respectively) for this antibiotic during QC testing. Testing was performed in duplicate.

⁸Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

30 JUL 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 28866**

Catalog No. NR-55540

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 28866 was isolated in 2015 from a human wound sample in North America as part of a global surveillance program. NR-55540 was deposited as a multidrug-resistant strain, sensitive to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ertapenem, imipenem, meropenem, piperacillin/tazobactam, tigecycline and trimethoprim/sulfamethoxazole, intermediately resistant to tobramycin and resistant to ciprofloxacin, gentamicin, levofloxacin and tetracycline. NR-55540 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050377

Manufacturing Date: 16FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, low convex, entire, mucoid, smooth and cream Non-motile <i>K. pneumoniae</i> (94%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Resistant Resistant Resistant Sensitive Sensitive Sensitive Sensitive Resistant Sensitive Sensitive Resistant Sensitive Intermediate Sensitive	Sensitive (≤ 2 µg/mL) Sensitive (≤ 2 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.023 µg/mL) Sensitive (0.032 µg/mL) Sensitive (≤ 1 µg/mL) Resistant (4 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.19 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (≤ 0.5 µg/mL) ³ Intermediate (8 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 28866 (GenBank: JAGYDT010000088.1)	99.6% sequence identity to <i>K. pneumoniae</i> , strain MRSN 28866 (GenBank: JAGYDT010000088.1) ⁴

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

/Sonia Bjorum Brower/

Sonia Bjorum Brower

30 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth


¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Heather Couch/

Heather Couch

Program Manager or designee, ATCC Federal Solutions

11 MAY 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 28887**

Catalog No. NR-55542

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 28887 was isolated in 2014 from a human urine sample in North America as part of a global surveillance program. NR-55542 was deposited as a susceptible strain, sensitive to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tigecycline, tobramycin, tetracycline and trimethoprim/sulfamethoxazole. NR-55542 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050381

Manufacturing Date: 25FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive	Sensitive (≤ 2 µg/mL) Sensitive (8 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.19 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.064 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.25 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Sensitive (2 µg/mL) Sensitive (1 to 2 µg/mL) ³ Sensitive (≤ 1 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 28887 (GenBank: JAGYDR010000138.1)	99.4% sequence identity to <i>K. pneumoniae</i> , strain MRSN 28887 (GenBank: JAGYDR010000138.1) ⁴

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Heather Couch/

Heather Couch

Program Manager or designee, ATCC Federal Solutions

25 APR 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Viability	Growth	Growth


¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

29 JUN 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 365679**

Catalog No. NR-55544

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 365679 was isolated in 2016 from a human urine sample in Asia as part of a global surveillance program. NR-55544 was deposited as an extensively drug-resistant strain, sensitive to tigecycline, intermediately resistant to amikacin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, ciprofloxacin, ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55544 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050648

Manufacturing Date: 02MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Intermediate Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Resistant	Intermediate (32 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Resistant (≥ 16 µg/mL) Resistant (2 to 3 µg/mL) ^{3,4} Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1110 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 365679 (GenBank: JAGYDP010000045.1)	99.6% sequence identity to <i>K. pneumoniae</i> , strain MRSN 365679 (GenBank: JAGYDP010000045.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴*K. pneumoniae*, strain MRSN 365679 was deposited as sensitive to tigecycline, but showed a MIC of 2 to 3 µg per mL (interpreted as resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

27 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 366562**

Catalog No. NR-55545

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 366562 was isolated in 2015 from a human urine sample in Africa as part of a global surveillance program. NR-55545 was deposited as an extensively drug-resistant strain, sensitive to ceftazidime/avibactam and resistant to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55545 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050651

Manufacturing Date: 25FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant	Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (1 µg/mL) Resistant (48 to 96 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (6 to 8 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 8 µg/mL) ³ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1460 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 366562 (GenBank: JAGYDO010000055.1)	99.6% sequence identity to <i>K. pneumoniae</i> , strain MRSN 366562 (GenBank: JAGYDO010000055.1) ⁴

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

/Sonia Bjorum Brower/

Sonia Bjorum Brower

18 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 368001**

Catalog No. NR-55546

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 368001 was isolated in 2015 from a human wound sample in Africa as part of a global surveillance program. NR-55546 was deposited as an extensively drug-resistant strain (XDR), sensitive to tetracycline and tigecycline and resistant to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, ceftazidime/avibactam, ceftolozane/tazobactam, ciprofloxacin, ertapenem, meropenem, gentamicin, imipenem, levofloxacin, piperacillin/tazobactam, tobramycin and trimethoprim/sulfamethoxazole. NR-55546 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050653

Manufacturing Date: 03MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (95%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Resistant Sensitive Sensitive Resistant Resistant	Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 0.5 µg/mL) ³ Resistant (≥ 16 µg/mL) Sensitive (≤ 32 µg/mL) ⁴
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1480 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 368001 (GenBank: JAGYDN010000074.1)	99.9% sequence identity to <i>K. pneumoniae</i> , strain MRSN 368001 (GenBank: JAGYDN010000074.1)

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴*K. pneumoniae*, strain MRSN 368001 was deposited as resistant to trimethoprim/sulfamethoxazole, but showed a MIC of ≤ 32 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

/Heather Couch/

Heather Couch

13 MAY 2022

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 368320**

Catalog No. NR-55547

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 368320 was isolated in 2015 from a human blood sample in Africa as part of a global surveillance program. NR-55547 was deposited as an extensively drug-resistant strain (XDR), sensitive to amikacin and ceftazidime/avibactam, intermediately resistant to imipenem and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55547 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050657

Manufacturing Date: 03MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Resistant Resistant Intermediate Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (1.5 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 32 µg/mL) ³ Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 8 µg/mL) ⁴ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 368320 (GenBank: JAGYDM010000069.1)	99.2% sequence identity to <i>K. pneumoniae</i> , strain MRSN 368320 (GenBank: JAGYDM010000069.1)

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 368320 was deposited as intermediately resistant to imipenem, but showed a MIC of ≥ 32 µg per mL (interpreted as resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

Figure 1: Colony Morphology



/Heather Couch/
Heather Couch

12 MAY 2022

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 371351**

Catalog No. NR-55548

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 371351 was isolated in 2003 from a human wound in Europe as part of a global surveillance program. MRSN 371351 was deposited as a multidrug-resistant strain, sensitive to amikacin, cefepime, ceftazidime/avibactam, ceftolozane/tazobactam, ciprofloxacin, ertapenem, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline and tigecycline, intermediately resistant to ampicillin/sulbactam, and resistant to aztreonam, ceftazidime, ceftriaxone, gentamicin, tobramycin and trimethoprim/sulfamethoxazole. NR-55548 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050659

Manufacturing Date: 04MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (93%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Intermediate Resistant Sensitive Resistant Sensitive Sensitive Sensitive Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Resistant Resistant Resistant	Sensitive (16 µg/mL) Sensitive (8 µg/mL) ³ Resistant (≥ 64 µg/mL) Sensitive (≤1 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.25 µg/mL) Sensitive (0.19 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.25 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.19 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤4 µg/mL) Sensitive (≤1 µg/mL) Sensitive (≤ 0.5 µg/mL) ⁴ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1480 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 371351 (GenBank: JAGYDL010000087.1)	99.9% sequence identity to <i>K. pneumoniae</i> , strain MRSN 371351 (GenBank: JAGYDL010000087.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)


²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Heather Couch/
Heather Couch

11 MAY 2022

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 374613**

Catalog No. NR-55549

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 374613 was isolated in 2016 from a human urine sample in North America as part of a global surveillance program. NR-55549 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tigecycline and tobramycin, intermediately resistant to cefepime and resistant to ampicillin/sulbactam, aztreonam, ceftazidime, ceftriaxone, tetracycline and trimethoprim/sulfamethoxazole. NR-55549 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050664

Manufacturing Date: 04MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided. E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Intermediate Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Resistant Sensitive Sensitive Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Sensitive (2 µg/mL) ³ Sensitive (1 to 1.5 µg/mL) ⁴ Sensitive (≤ 1 µg/mL) ⁵ Sensitive (0.19 µg/mL) Sensitive (0.25 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.75 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.125 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (16 µg/mL) Resistant (≥ 16 µg/mL) Resistant (2 µg/mL) ^{6,7} Sensitive (≤ 1 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 374613 (GenBank: JAGYDK010000086.1)	99.7% sequence identity to <i>K. pneumoniae</i> , strain MRSN 374613 (GenBank: JAGYDK010000086.1) ⁸

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 374613 was deposited as resistant to aztreonam, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴*K. pneumoniae*, strain MRSN 374613 was deposited as intermediately resistant to cefepime, but showed a MIC of 1 to 1.5 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵*K. pneumoniae*, strain MRSN 374613 was deposited as resistant to ceftazidime, but showed a MIC of ≤ 1 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁶MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁷The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁸Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

17 JUN 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 375436**

Catalog No. NR-55550

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 375436 was isolated in 2016 from a human urine sample in North America as part of a global surveillance program. NR-55550 was deposited as an extensively drug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, tigecycline, tetracycline and tobramycin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam and trimethoprim/sulfamethoxazole. NR-55550 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050666

Manufacturing Date: 04MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Sensitive Sensitive Sensitive Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (24 to 32 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (4 µg/mL) Resistant (48 µg/mL) Resistant (32 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Sensitive (4 µg/mL) Sensitive (1 µg/mL) ³ Sensitive (2 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 375436 (GenBank: JAGYDJ010000110.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 375436 (GenBank: JAGYDJ010000110.1) ⁴

TEST	SPECIFICATIONS	RESULTS
Purity 11 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

/Sonia Bjorum Brower/

Sonia Bjorum Brower

26 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

01 AUG 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 401050**

Catalog No. NR-55552

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 401050 was isolated in 2016 from a human fluid sample in North America as part of a global surveillance program. NR-55552 was deposited as a multidrug-resistant strain, sensitive to amikacin, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ertapenem, gentamicin, imipenem, meropenem, tetracycline, tobramycin and tigecycline, intermediately resistant to levofloxacin, and resistant to ampicillin/sulbactam, ciprofloxacin, piperacillin/tazobactam and trimethoprim/sulfamethoxazole. NR-55552 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050670

Manufacturing Date: 04MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Resistant Sensitive Sensitive Sensitive Intermediate Sensitive Resistant Sensitive Sensitive Sensitive Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.50 µg/mL) Sensitive (0.50 to 1 µg/mL) Sensitive (≤ 1 µg/mL) Intermediate (3 µg/mL) ³ Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.125 to 0.19 µg/mL) Intermediate (4 µg/mL) Sensitive (≤ 0.25 µg/mL) Resistant (≥ 128 µg/mL) Sensitive (4 µg/mL) Sensitive (0.75 to 1 µg/mL) ⁴ Sensitive (≤ 1 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 401050 (GenBank: JAGYDH010000094.1)	99.8% sequence identity to <i>K. pneumoniae</i> , strain MRSN 401050 (GenBank: JAGYDH010000094.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)


²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Heather Couch/
Heather Couch

06 MAY 2022

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 410359**

Catalog No. NR-55553

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 410359 was isolated in 2016 from a human urine sample in North America as part of a global surveillance program. NR-55553 was deposited as an extensively drug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, imipenem and meropenem, intermediately resistant to ciprofloxacin, levofloxacin, piperacillin/tazobactam and tobramycin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, gentamicin, tetracycline, tigecycline and trimethoprim/sulfamethoxazole. NR-55553 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050672

Manufacturing Date: 25FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Sensitive Sensitive Resistant Intermediate Sensitive Resistant Sensitive Intermediate Sensitive Intermediate Resistant Resistant Intermediate Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (2 µg/mL) ³ Resistant (16 µg/mL) Sensitive (0.5 µg/mL) Sensitive (0.38 µg/mL) Resistant (≥ 64 µg/mL) Intermediate (2 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.19 µg/mL) Intermediate (4 µg/mL) Sensitive (≤ 0.25 µg/mL) Intermediate (32 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 8 µg/mL) ⁴ Intermediate (8 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1480 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 410359 (GenBank: JAGYDG010000093.1)	99.6% sequence identity to <i>K. pneumoniae</i> , strain MRSN 410359 (GenBank: JAGYDG010000093.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 410359 was deposited as resistant to cefepime, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

15 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 414780**

Catalog No. NR-55554

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 414780 was isolated in 2016 from a human urine sample in North America as part of a global surveillance program. NR-55554 was deposited as an extensively drug-resistant strain (XDR), sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, imipenem and meropenem, intermediately resistant to piperacillin/tazobactam, tigecycline and tobramycin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, ciprofloxacin, gentamicin, levofloxacin, tetracycline and trimethoprim/sulfamethoxazole. NR-55554 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050674

Manufacturing Date: 25FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card) VITEK® MS (MALDI-TOF)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%) <i>K. pneumoniae</i>	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%) <i>K. pneumoniae</i> (99.9%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Sensitive Resistant Sensitive Sensitive Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Intermediate Resistant Intermediate Intermediate Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (32 µg/mL) Sensitive (2 µg/mL) Resistant (16 µg/mL) Sensitive (1.0 to 1.5 µg/mL) Sensitive (1.5 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.25 to 0.75 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (0.5 µg/mL) Inconclusive ³ Sensitive (0.5 to 0.75 µg/mL) Sensitive (≤ 0.25 µg/mL) Intermediate (32 to 64 µg/mL) Resistant (≥ 16 µg/mL) Inconclusive ^{4,5} Intermediate (8 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1450 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 414780 (GenBank: JAGYDF010000080.1 and JAGYDF010000081.1)	99.8% sequence identity to <i>K. pneumoniae</i> , strain MRSN 414780 (GenBank: JAGYDF010000080.1 and JAGYDF010000081.1) ⁶

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK®2 GN74 and ETEST®.


³*K. pneumoniae*, strain MRSN 414780 was deposited as being resistant to imipenem. Repeated antibiotic susceptibility testing determined that for strain MRSN 414780, the imipenem MICs are 0.75 µg per mL and 1 µg per mL (interpreted as sensitive), 1.5 µg per mL (interpreted as intermediately resistant) and 2 µg per mL and 3 µg per mL (interpreted as resistant).

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵*K. pneumoniae*, strain MRSN 414780 was deposited as being intermediately resistant to tigecycline. Repeated antibiotic susceptibility testing determined that for strain MRSN 414780, the tigecycline MICs are 0.75 µg per mL and 2 µg per mL, which are interpreted as sensitive and resistant, respectively.

⁶Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

01 AUG 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 430405**

Catalog No. NR-55555

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 430405 was isolated in 2013 from a human blood sample in the Middle East as part of a global surveillance program. NR-55555 was deposited as an extensively drug-resistant strain, sensitive to ceftazidime/avibactam and tigecycline, intermediately resistant to gentamicin and tetracycline and resistant to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tobramycin and trimethoprim/sulfamethoxazole. NR-55555 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050676

Manufacturing Date: 02MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Resistant Intermediate Resistant Resistant Resistant Resistant Resistant Intermediate Sensitive Resistant Resistant	Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (2 µg/mL) Resistant (64 to 96 µg/mL) Resistant (≥ 64 µg/mL) Resistant (> 32 µg/mL) Resistant (≥ 8 µg/mL) Intermediate (8 µg/mL) Resistant (> 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Intermediate (8 µg/mL) Resistant (1.5 µg/mL) ^{3,4} Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1460 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 430405 (GenBank: JAGYDE010000058.1)	99.4% sequence identity to <i>K. pneumoniae</i> , strain MRSN 430405 (GenBank: JAGYDE010000058.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)


²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK®2 GN74 and ETEST®.

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

26 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 430414**

Catalog No. NR-55556

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 430414 was isolated in 2013 from a human blood sample in the Middle East as part of a global surveillance program. NR-55556 was deposited as an extensively drug-resistant strain, sensitive to amikacin, ceftazidime/avibactam and imipenem, intermediately resistant to tigecycline and tobramycin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline and trimethoprim/sulfamethoxazole. NR-55556 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050678

Manufacturing Date: 02MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology	Report results	Circular, convex, entire, smooth and cream
Motility (wet mount)	Report results	Non-motile
VITEK® 2 (GN card)	<i>K. pneumoniae</i> (≥ 89%)	<i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2}		
Amikacin	Sensitive	Sensitive (16 µg/mL)
Ampicillin/sulbactam	Resistant	Resistant (≥ 32 µg/mL)
Aztreonam	Resistant	Resistant (≥ 64 µg/mL)
Cefepime	Resistant	Resistant (≥ 64 µg/mL)
Ceftazidime	Resistant	Resistant (≥ 64 µg/mL)
Ceftazidime/avibactam	Sensitive	Sensitive (1 to 1.5 µg/mL)
Ceftolozane/tazobactam	Resistant	Resistant (24 to 48 µg/mL)
Ceftriaxone	Resistant	Resistant (≥ 64 µg/mL)
Ciprofloxacin	Resistant	Resistant (≥ 32 µg/mL)
Ertapenem	Resistant	Resistant (≥ 8 µg/mL)
Gentamicin	Resistant	Resistant (≥ 16 µg/mL)
Imipenem	Sensitive	Resistant (3 µg/mL) ³
Levofloxacin	Resistant	Resistant (≥ 8 µg/mL)
Meropenem	Resistant	Resistant (≥ 16 µg/mL)
Piperacillin/tazobactam	Resistant	Resistant (≥ 128 µg/mL)
Tetracycline	Resistant	Resistant (≥ 16 µg/mL)
Tigecycline	Intermediate	Resistant (2 to 3 µg/mL) ^{4,5}
Tobramycin	Intermediate	Intermediate (8 µg/mL)
Trimethoprim/sulfamethoxazole	Resistant	Resistant (≥ 320 µg/mL)
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 430414 (GenBank: JAGYDD010000079.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 430414 (GenBank: JAGYDD010000079.1) ⁶

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].


³*K. pneumoniae*, strain MRSN 430414 was deposited as sensitive to imipenem, but showed a MIC of 3 µg per mL (interpreted as resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵*K. pneumoniae*, strain MRSN 430414 was deposited as intermediately resistant to tigecycline, but showed a MIC of 2 to 3 µg per mL (interpreted as resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

⁶Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

15 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

30 JUL 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 468268**

Catalog No. NR-55558

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 468268 was isolated in 2017 from a human urine sample in North America as part of a global surveillance program. NR-55558 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ciprofloxacin, ertapenem, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tigecycline and tobramycin, intermediately resistant to ampicillin/sulbactam and resistant to aztreonam, cefepime, ceftazidime, ceftriaxone, gentamicin, tetracycline and trimethoprim/sulfamethoxazole. NR-55558 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050682

Manufacturing Date: 02MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Intermediate Resistant Resistant Resistant Sensitive Sensitive Resistant Sensitive Sensitive Resistant Sensitive Sensitive Sensitive Sensitive Resistant Sensitive Sensitive Resistant	Sensitive (≤ 2 µg/mL) Intermediate (16 µg/mL) Sensitive (2 µg/mL) ³ Sensitive (2 µg/mL) ⁴ Sensitive (4 µg/mL) ⁵ Sensitive (0.125 µg/mL) Sensitive (0.19 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.38 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.25 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Resistant (≥ 16 µg/mL) Resistant (2 µg/mL) ^{6,7} Sensitive (2 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 468268 (GenBank: JAGYDB01000076.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 468268 (GenBank: JAGYDB01000076.1) ⁸

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 468268 was deposited as resistant to aztreonam, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴*K. pneumoniae*, strain MRSN 468268 was deposited as resistant to cefepime, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.


⁵*K. pneumoniae*, strain MRSN 468268 was deposited as resistant to ceftazidime, but showed a MIC of 4 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁶MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁷The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁸Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

31 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 479404**

Catalog No. NR-55559

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 479404 was isolated in 2017 from a human wound sample in Asia as part of a global surveillance program. NR-55559 was deposited as an extensively drug-resistant strain, sensitive to amikacin, gentamicin and tigecycline and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55559 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050684

Manufacturing Date: 02MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Resistant	Sensitive (16 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 1 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Resistant (≥ 16 µg/mL) Resistant (2 µg/mL) ^{3,4} Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (1480 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 479404 (GenBank: JAGYDA010000110.1)	99.8% sequence identity to <i>K. pneumoniae</i> , strain MRSN 479404 (GenBank: JAGYDA010000110.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)


²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

26 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 10 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁴*K. pneumoniae*, strain MRSN 499958 was deposited as resistant to cefepime, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶Also consistent with other *Klebsiella* species

/Sonia Bjorum Brower/

Sonia Bjorum Brower

08 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 511348**

Catalog No. NR-55561

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 511348 was isolated in 2014 from a human sample in the Middle East as part of a global surveillance program. NR-55561 was deposited as an extensively drug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, levofloxacin, tetracycline and tigecycline, intermediately resistant to ciprofloxacin and imipenem and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ertapenem, gentamicin, meropenem, piperacillin/tazobactam, tobramycin and trimethoprim/sulfamethoxazole. NR-55561 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050688

Manufacturing Date: 02MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Intermediate Resistant Resistant Intermediate Sensitive Resistant Resistant Sensitive Sensitive Resistant Resistant	Sensitive (4 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 64 µg/mL) Inconclusive ³ Inconclusive ⁴ Sensitive (4 µg/mL) Resistant (192 µg/mL) Resistant (≥ 256 µg/mL) Intermediate (1.5 µg/mL) Resistant (1.5 to 2 µg/mL) Resistant (≥ 16 µg/mL) Intermediate (1.5 to 2 µg/mL) Sensitive (1 µg/mL) Sensitive (0.75 to 1 µg/mL) ⁵ Resistant (≥ 128 µg/mL) Sensitive (4 µg/mL) Sensitive (≤ 0.5 µg/mL) ⁶ Resistant (≥ 16 µg/mL) Resistant (160 to ≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 511348 (GenBank: JAGYCY010000118.1)	99.9% sequence identity to <i>K. pneumoniae</i> , strain MRSN 511348 (GenBank: JAGYCY010000118.1) ⁷

TEST	SPECIFICATIONS	RESULTS
Purity 10 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*Klebsiella pneumoniae*, strain MRSN 511348 was deposited as being resistant to cefepime. Repeated antibiotic susceptibility testing determined that for strain MRSN 511348, the cefepime MICs are 2 µg per mL and ≥ 64 µg per mL, which are interpreted as sensitive and resistant, respectively.

⁴*Klebsiella pneumoniae*, strain MRSN 511348 was deposited as being resistant to ceftazidime. Repeated antibiotic susceptibility testing determined that for strain MRSN 511348, the ceftazidime MICs are 4 µg per mL and ≥ 64 µg per mL, which are interpreted as sensitive and resistant, respectively.

⁵*K. pneumoniae*, strain MRSN 511348 was deposited as resistant to meropenem, but showed a MIC of 0.75 to 1 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁶MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁷Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

11 JUL 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 513382**

Catalog No. NR-55562

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 513382 was isolated in 2014 from a human sample in the Middle East as part of a global surveillance program. NR-55562 was deposited as a susceptible strain, sensitive to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55562 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050690

Manufacturing Date: 25FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream Non-motile <i>K. pneumoniae</i> (94%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive	Sensitive (≤ 2 µg/mL) Sensitive (4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.125 µg/mL) Sensitive (0.25 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.094 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.25 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Sensitive (1 µg/mL) Sensitive (0.5 µg/mL) ³ Sensitive (≤ 1 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1460 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 513382 (GenBank: JAGYCX010000073.1)	99.4% sequence identity to <i>K. pneumoniae</i> , strain MRSN 513382 (GenBank: JAGYCX010000073.1)

TEST	SPECIFICATIONS	RESULTS
Purity 11 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

/Heather Couch/

Heather Couch

19 MAY 2022

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Viability	Growth	Growth


¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

06 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 515432**

Catalog No. NR-55564

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 515432 was isolated in 2014 from an unknown human sample in the Middle East as part of a global surveillance program. NR-55564 was deposited as a multidrug-resistant strain (MDR), sensitive to amikacin, ceftazidime/avibactam, ertapenem, imipenem, meropenem, tetracycline and tigecycline and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, gentamicin, levofloxacin, piperacillin/tazobactam, tobramycin and trimethoprim/sulfamethoxazole. NR-55564 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051098

Manufacturing Date: 23MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (95%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Sensitive Resistant Sensitive Resistant Resistant Sensitive Resistant Resistant Sensitive Resistant	Sensitive (16 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Intermediate (8 µg/mL) ³ Resistant (≥ 64 µg/mL) Sensitive (3 µg/mL) Resistant (48 µg/mL) Resistant (≥ 64 µg/mL) Resistant (8 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.25 to 0.38 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 0.25 µg/mL) Resistant (≥ 128 µg/mL) Sensitive (2 µg/mL) Sensitive (1 µg/mL) ⁴ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 515432 (GenBank: JAGYCV010000163.1)	99.7% sequence identity to <i>K. pneumoniae</i> , strain MRSN 515432 (GenBank: JAGYCV010000163.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®]

³The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

18 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 516635**

Catalog No. NR-55565

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 516635 was isolated in 2016 from an unknown human sample in the Middle East as part of a global surveillance program. NR-55565 was deposited as a multidrug-resistant strain (MDR), sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ciprofloxacin, ertapenem, imipenem, levofloxacin, meropenem, piperacillin/tazobactam and tigecycline, intermediately resistant to tobramycin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, gentamicin, tetracycline and trimethoprim/sulfamethoxazole. NR-55565 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051100

Manufacturing Date: 30MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Sensitive Sensitive Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Resistant Sensitive Intermediate Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (16 µg/mL) Inconclusive ³ Resistant (16 µg/mL) Sensitive (0.25 µg/mL) Sensitive (0.38 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.38 to 0.5 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.25 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (≤ 0.5 µg/mL) ⁴ Intermediate (8 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1430 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 516635 (GenBank: JAGYCU010000101.1)	99.2% sequence identity to <i>K. pneumoniae</i> , strain MRSN 516635 (GenBank: JAGYCU010000101.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®]

³*K. pneumoniae*, strain MRSN 516635 was deposited as being resistant to cefepime. Antibiotic susceptibility testing performed in duplicate determined that for strain MRSN 516635, the cefepime MICs are 2 µg per mL and 4 µg per mL, which are interpreted as sensitive and intermediately resistant, respectively.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

/Sonia Bjorum Brower/

Sonia Bjorum Brower

11 AUG 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth


¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

15 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 518712**

Catalog No. NR-55567

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 518712 was isolated in 2016 from a human sample in the Middle East as part of a global surveillance program. NR-55567 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline and tobramycin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone and trimethoprim/sulfamethoxazole. NR-55567 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051104

Manufacturing Date: 23MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (16 µg/mL) Sensitive (2 µg/mL) ³ Resistant (16 µg/mL) Sensitive (0.19 µg/mL) Sensitive (0.19 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.25 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.25 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 0.5 µg/mL) ⁴ Sensitive (≤ 1 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1460 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 518712 (GenBank: JAGYCS010000120.1)	99.4% sequence identity to <i>K. pneumoniae</i> , strain MRSN 518712 (GenBank: JAGYCS010000120.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 518712 was deposited as resistant to cefepime, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

16 JUL 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 526410**

Catalog No. NR-55568

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 526410 was isolated in 2017 from a human urine sample in Europe as part of a global surveillance program. NR-55568 was deposited as a multidrug-resistant strain, sensitive to amikacin and ceftazidime/avibactam, intermediately resistant to imipenem and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55568 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051106

Manufacturing Date: 23MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Resistant Resistant Intermediate Resistant Resistant Resistant Resistant Resistant Resistant Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.75 µg/mL) Resistant (8 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Intermediate (2 µg/mL) Resistant (≥ 8 µg/mL) Resistant (4 µg/mL) Resistant (≥ 128 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 8 µg/mL) ³ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 526410 (GenBank: JAGYCR010000086.1)	99.6% sequence identity to <i>K. pneumoniae</i> , strain MRSN 526410 (GenBank: JAGYCR010000086.1)


TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

09 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 539414**

Catalog No. NR-55569

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 539414 was isolated in 2017 from a human blood sample in North America as part of a global surveillance program. NR-55569 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam and tigecycline, intermediately resistant to aztreonam and ceftazidime and resistant to ampicillin/sulbactam, cefepime, ceftriaxone, tetracycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55569 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051108

Manufacturing Date: 24MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (95%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Intermediate Resistant Intermediate Sensitive Sensitive Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Resistant Sensitive Resistant Resistant	Sensitive (4 µg/mL) Resistant (≥ 32 µg/mL) Resistant (48 µg/mL) ³ Resistant (≥ 64 µg/mL) Intermediate (8 to 12 µg/mL) Sensitive (0.19 µg/mL) Sensitive (0.19 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.064 to 0.094 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.19 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (8 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (1 µg/mL) ⁴ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 539414 (GenBank: JAGYQCQ010000070.1)	99.8% sequence identity to <i>K. pneumoniae</i> , strain MRSN 539414 (GenBank: JAGYQCQ010000070.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)


²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 539414 was deposited as intermediate to aztreonam, but showed a MIC of 48 µg per mL (interpreted as resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

13 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 546733**

Catalog No. NR-55570

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 546733 was isolated in 2018 from a human urine sample in North America as part of a global surveillance program. NR-55570 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, imipenem, levofloxacin, meropenem, piperacillin/tazobactam and tigecycline, intermediately resistant to ciprofloxacin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, gentamicin, tetracycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55570 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051110

Manufacturing Date: 24MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Sensitive Sensitive Resistant Intermediate Sensitive Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Resistant Sensitive Resistant Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (16 µg/mL) Sensitive (2 µg/mL) ³ Intermediate (8 µg/mL) ⁴ Sensitive (0.125 µg/mL) Sensitive (0.19 µg/mL) Resistant (≥ 64 µg/mL) Intermediate (1.5 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.25 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (16 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (≤ 0.5 µg/mL) ⁵ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 546733 (GenBank: JAGYCP010000085.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 546733 (GenBank: JAGYCP010000085.1) ⁶

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].


³*K. pneumoniae*, strain MRSN 546733 was deposited as resistant to cefepime, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

27 JUL 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 560539**

Catalog No. NR-55571

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 560539 was isolated in 2018 from a human urine sample in North America as part of a global surveillance program. NR-55571 was deposited as an extensively drug-resistant strain (XDR), sensitive to amikacin, ceftazidime/avibactam, and tigecycline, intermediate to ceftazidime and tobramycin, resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, and trimethoprim/sulfamethoxazole. NR-55571 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70051112

Manufacturing Date: 23MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Intermediate Sensitive Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Sensitive Intermediate Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (2 µg/mL) ³ Resistant (24 µg/mL) ⁴ Sensitive (1.5 µg/mL) Resistant (24 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (4 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Resistant (≥ 16 µg/mL) Resistant (4.0 µg/mL) ^{5,6} Intermediate (8 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 560539 (GenBank: JAGYCO010000175.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 560539 (GenBank: JAGYCO010000175.1)

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].


³*K. pneumoniae*, strain MRSN 560539 was deposited as resistant to cefepime but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴*K. pneumoniae*, strain MRSN 560539 was deposited as intermediate to ceftazidime but showed a MIC of 24 µg per mL (interpreted as resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵*K. pneumoniae*, strain MRSN 560539 was deposited as immediate to tigecycline but showed a MIC of 4 µg per mL (interpreted as resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

⁶MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

02 JUN 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 562722**

Catalog No. NR-55572

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 562722 was isolated in 2018 from a human urine sample in North America as part of a global surveillance program. NR-55572 was deposited as a susceptible strain, sensitive to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55572 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70051114

Manufacturing Date: 23MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive	Sensitive (≤ 2 µg/mL) Sensitive (4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.19 µg/mL) Sensitive (0.19 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.032 to 0.047 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.25 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (1 µg/mL) ³ Sensitive (≤ 1 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 562722 (GenBank: JAGYCN010000106.1)	99.6% sequence identity to <i>K. pneumoniae</i> , strain MRSN 562722 (GenBank: JAGYCN010000106.1) ⁴

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Heather Couch/

Heather Couch

26 MAY 2022

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 564304 was deposited as resistant to cefepime, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶*K. pneumoniae*, strain MRSN 564304 was deposited as resistant to tobramycin, but showed a MIC of 6 µg per mL (interpreted as intermediate) for this antibiotic during QC testing. Testing was performed in duplicate.

/Sonia Bjorum Brower/

Sonia Bjorum Brower

02 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth


¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

29 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 581745**

Catalog No. NR-55575

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 581745 was isolated in 2018 from a human urine sample in Asia as part of a global surveillance program. NR-55575 was deposited as an extensively drug-resistant strain (XDR), sensitive to amikacin and gentamicin, intermediately resistant to levofloxacin and meropenem, and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, imipenem, piperacillin/tazobactam, tetracycline, tigecycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55575 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70051120

Manufacturing Date: 18MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card) VITEK® MS (MALDI-TOF)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%) <i>K. pneumoniae</i>	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (95%) <i>K. pneumoniae</i> (99.9%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Intermediate Intermediate Resistant Resistant Resistant Resistant Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (> 256 µg/mL) Resistant (> 256 µg/mL) Resistant (≥ 64 µg/mL) Resistant (> 32 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 1 µg/mL) Resistant (> 32 µg/mL) Intermediate (4 µg/mL) Resistant (> 32 µg/mL)³ Resistant ≥ 128 µg/mL Resistant (≥ 16 µg/mL) Resistant (≥ 8 µg/mL) ⁴ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 581745 (GenBank: JAGYCK010000159.1)	99.6% sequence identity to <i>K. pneumoniae</i> , strain MRSN 581745 (GenBank: JAGYCK010000159.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 581745 was deposited as intermediately resistant to meropenem, but showed a MIC of >32 µg per mL (interpreted as resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Heather Couch/
Heather Couch

31 MAY 2022

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 582610**

Catalog No. NR-55576

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 582610 was isolated in 2017 from a human respiratory sample in Asia as part of a global surveillance program. NR-55576 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem and tigecycline, intermediately resistant to ceftazidime, ceftolozane/tazobactam, piperacillin/tazobactam and tobramycin, and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftriaxone, tetracycline and trimethoprim/sulfamethoxazole. NR-55576 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70051122

Manufacturing Date: 18MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Intermediate Sensitive Intermediate Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Intermediate Resistant Sensitive Intermediate Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (16 µg/mL) Intermediate (4 µg/mL) ³ Resistant (12 to 16 µg/mL) ⁴ Sensitive (0.38 to 0.50 µg/mL) Sensitive (0.75 µg/mL) ⁵ Resistant (≥ 64 µg/mL) Sensitive (0.094 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.19 to 0.25 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Intermediate (32 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (2 µg/mL) ⁶ Intermediate (6 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 582610 (GenBank: JAGYCJ010000113.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 582610 (GenBank: JAGYCJ010000113.1) ⁷

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 582610 was deposited as resistant to cefepime, but showed a MIC of 4 µg per mL (interpreted as intermediate) for this antibiotic during QC testing. Testing was performed in duplicate.


⁴The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁵*K. pneumoniae*, strain MRSN 582610 was deposited as intermediately resistant to ceftolozane/tazobactam, but showed a MIC of 0.75 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁶MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁷Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

13 JUN 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].


³*K. pneumoniae*, strain MRSN 583114 was deposited as resistant to aztreonam, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴*K. pneumoniae*, strain MRSN 583114 was deposited as resistant to ceftazidime, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

29 JUN 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 583141**

Catalog No. NR-55578

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 583141 was isolated in 2018 from a human urine sample in North America as part of a global surveillance program. NR-55578 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, imipenem, levofloxacin, meropenem, piperacillin/tazobactam and tigecycline, intermediately resistant to ciprofloxacin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, gentamicin, tetracycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55578 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70051126

Manufacturing Date: 18MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card) VITEK® MS (MALDI-TOF)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%) <i>K. pneumoniae</i>	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (95%) <i>K. pneumoniae</i> (99.9%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Sensitive Resistant Intermediate Sensitive Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Resistant Sensitive Resistant Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (16 µg/mL) Sensitive (2 µg/mL) ³ Sensitive (4 µg/mL) ⁴ Sensitive (0.19 µg/mL) Sensitive (0.25 to 0.38 µg/mL) Resistant (≥ 64 µg/mL) Resistant (3 µg/mL) ⁵ Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.38 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (16 µg/mL) Resistant (≥ 16 µg/mL) Inconclusive ^{6,7} Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1460 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 583141 (GenBank: JAGYCH010000084.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 583141 (GenBank: JAGYCH010000084.1) ⁸

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK®2 GN74 and ETEST®.

³*K. pneumoniae*, strain MRSN 583141 was deposited as resistant to cefepime, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴*K. pneumoniae*, strain MRSN 583141 was deposited as resistant to ceftazidime, but showed a MIC of 4 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.


⁵The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁶MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁷*K. pneumoniae*, strain MRSN 583141 was deposited as being sensitive to tigecycline. Antibiotic susceptibility testing performed in duplicate determined that for strain MRSN 583141, the tigecycline MICs are 1 µg per mL and 2 µg per mL, which are interpreted as sensitive and resistant, respectively.

⁸Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

31 JUL 2022

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth


¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK®2 GN74 and ETEST®.

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
 Sonia Bjorum Brower

07 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.
 You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 599975**

Catalog No. NR-55580

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 599975 was isolated in 2018 from a human urine sample in North America as part of a global surveillance program. NR-55580 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, imipenem, levofloxacin, meropenem, piperacillin/tazobactam and tigecycline, intermediately resistant to ciprofloxacin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, gentamicin, tetracycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55580 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70051130

Manufacturing Date: 23MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Sensitive Resistant Intermediate Sensitive Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Resistant Sensitive Resistant Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (16 µg/mL) Sensitive (2 µg/mL) ³ Inconclusive ⁴ Sensitive (0.5 µg/mL) Sensitive (0.5 µg/mL) Resistant (≥ 64 µg/mL) Resistant (4 µg/mL) ⁵ Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.38 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (16 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (1 µg/mL) ⁶ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 599975 (GenBank: JAGYCF010000102.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 599975 (GenBank: JAGYCF010000102.1) ⁷

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO2 on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 599975 was deposited as resistant to cefepime, but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.


⁴*K. pneumoniae*, strain MRSN 599975 was deposited as being resistant to ceftazidime. Antibiotic susceptibility testing performed in duplicate determined that for strain MRSN 599975, the ceftazidime MICs are 8 µg per mL and 16 µg per mL, which are interpreted as intermediately resistant and resistant, respectively.

⁵The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁶MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁷Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

31 JUL 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 607210**

Catalog No. NR-55581

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 607210 was isolated in 2018 from a human urine sample in North America as part of a global surveillance program. NR-55581 was deposited as a multidrug-resistant strain (MDR), sensitive to ceftazidime/avibactam, ertapenem, imipenem, meropenem, tetracycline, tigecycline and trimethoprim/sulfamethoxazole, intermediately resistant to piperacillin/tazobactam and resistant to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, gentamicin, levofloxacin and tobramycin. NR-55581 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051132

Manufacturing Date: 30MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (95%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Sensitive Resistant Sensitive Sensitive Resistant Sensitive Intermediate Sensitive Sensitive Resistant Sensitive	Resistant (≥ 256 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (1.5 µg/mL) Resistant (12 to 48 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 256 µg/mL) Sensitive (0.38 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 0.25 µg/mL) Intermediate (64 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 0.5 µg/mL) ³ Resistant (≥ 16 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 607210 (GenBank: JAGYCE010000063.1)	99.8% sequence identity to <i>K. pneumoniae</i> , strain MRSN 607210 (GenBank: JAGYCE010000063.1) ⁴

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth


¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
 Sonia Bjorum Brower

21 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 613682**

Catalog No. NR-55582

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 613682 was isolated in 2018 from an environmental sample in Africa as part of a global surveillance program. NR-55582 was deposited as an extensively drug-resistant strain, sensitive to tigecycline and resistant to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55582 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was passaged once in Tryptic Soy broth for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed in Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051134

Manufacturing Date: 30MAR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (96%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Resistant Sensitive Resistant Resistant	Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (1 µg/mL) ³ Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 613682 (GenBank: JAGYCD010000073.1)	99.7% sequence identity to <i>K. pneumoniae</i> , strain MRSN 613682 (GenBank: JAGYCD010000073.1)

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

14 JUL 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 614201**

Catalog No. NR-55583

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 614201 was isolated in 2018 from an environmental sample in Africa as part of a global surveillance program. NR-55583 was deposited as a multidrug-resistant (MDR) strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline, tobramycin and trimethoprim/sulfamethoxazole and resistant to ampicillin/sulbactam, aztreonam, ceftazidime, ceftriaxone and cefepime. NR-55583 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051603

Manufacturing Date: 06APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Sensitive Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive	Sensitive (≤ 2 µg/mL) Sensitive (4 µg/mL) ³ Sensitive (2 µg/mL) ⁴ Sensitive (≤ 1 µg/mL) ⁵ Sensitive (4 µg/mL) ⁶ Sensitive (0.094 µg/mL) Sensitive (0.125 µg/mL) Resistant (16 µg/mL) Sensitive (0.016 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.25 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 0.5 µg/mL) ⁷ Sensitive (≤ 1 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 614201 (GenBank: JAGYCC010000052.1)	99.9% sequence identity to <i>K. pneumoniae</i> , strain MRSN 614201 (GenBank: JAGYCC010000052.1)

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 614201 was deposited as resistant to ampicillin/sulbactam but showed a MIC of 4 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴*K. pneumoniae*, strain MRSN 614201 was deposited as resistant to aztreonam but showed a MIC of 2 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵*K. pneumoniae*, strain MRSN 614201 was deposited as resistant to cefepime but showed a MIC of ≤ 1 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁶*K. pneumoniae*, strain MRSN 614201 was deposited as resistant to ceftazidime but showed a MIC of 4 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁷MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

23 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 669448**

Catalog No. NR-55584

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 669448 was isolated in 2017 from a blood sample in Europe as part of a global surveillance program. NR-55584 was deposited as an extensively drug-resistant (XDR) strain, sensitive to amikacin and ceftazidime/avibactam, intermediately resistant to tetracycline and tigecycline, and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tobramycin and trimethoprim/sulfamethoxazole. NR-55584 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051605

Manufacturing Date: 06APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Intermediate Intermediate Resistant Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (1.5 µg/mL) Resistant (48 to 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Intermediate (6 µg/mL) Resistant (4 µg/mL) ^{3,4} Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 669448 (GenBank: JAGYCB010000094.1)	99.8% sequence identity to <i>K. pneumoniae</i> , strain MRSN 669448 (GenBank: JAGYCB010000094.1)

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth


¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴*K. pneumoniae*, strain MRSN 669448 was deposited as intermediately resistant to tigecycline but showed a MIC of 4 µg per mL (interpreted as resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

20 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 669510**

Catalog No. NR-55585

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 669510 was isolated in 2018 from a human respiratory sample in Europe as part of a global surveillance program. NR-55585 was deposited as a multidrug-resistant (MDR) strain, sensitive to amikacin, ciprofloxacin, gentamicin, levofloxacin, tigecycline and tetracycline, intermediately resistant to aztreonam, ertapenem and tobramycin, and resistant to ampicillin/sulbactam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, imipenem, meropenem, piperacillin/tazobactam and trimethoprim/sulfamethoxazole. NR-55585 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70051607

Manufacturing Date: 07APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Intermediate Resistant Resistant Resistant Resistant Resistant Resistant Sensitive Intermediate Sensitive Resistant Sensitive Resistant Resistant Sensitive Sensitive Intermediate Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Intermediate (8 to 12 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.064 to 0.094 µg/mL) Sensitive (0.19 to 0.25 µg/mL) ³ Sensitive (≤ 1 µg/mL) Resistant (4 µg/mL) Sensitive (≤ 0.12 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (1 µg/mL) ⁴ Intermediate (8 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 669510 (GenBank: JAGYCA010000117.1)	99.7% sequence identity to <i>K. pneumoniae</i> , strain MRSN 669510 (GenBank: JAGYCA010000117.1)

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 669510 was deposited as intermediately resistant toertapenem, but showed a MIC of 0.19 to 0.25 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

/Sonia Bjorum Brower/
Sonia Bjorum Brower

14 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

/Sonia Bjorum Brower/

Sonia Bjorum Brower

20 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)


²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

30 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth


¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

20 JUN 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 681054**

Catalog No. NR-55589

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 681054 was isolated in 2015 from a urine sample in Asia as part of a global surveillance program. NR-55589 was deposited as an extensively drug-resistant (XDR) strain, sensitive to amikacin and tigecycline, intermediately resistant to tobramycin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline and trimethoprim/sulfamethoxazole. NR-55589 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051615

Manufacturing Date: 08APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Sensitive Intermediate Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 256 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (16 to 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Resistant (≥ 16 µg/mL) Resistant (2 µg/mL) ^{3,4} Intermediate (8 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 681054 (GenBank: JAGYBW010000061.1)	99.4% sequence identity to <i>K. pneumoniae</i> , strain MRSN 681054 (GenBank: JAGYBW010000061.1)

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth


¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

20 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 699478**

Catalog No. NR-55590

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 699478 was isolated in 2019 from a human urine sample in North America as part of a global surveillance program. NR-55590 was deposited as a susceptible strain, sensitive to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, ciprofloxacin, ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tobramycin, tetracycline, tigecycline and trimethoprim/sulfamethoxazole. NR-55590 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051617

Manufacturing Date: 06APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive	Sensitive (≤ 2 µg/mL) Sensitive (4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.064 µg/mL) Sensitive (0.064 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.023 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.25 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 0.5 µg/mL) ³ Sensitive (≤ 1 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 699478 (GenBank: JAGYBV010000073.1)	99.9% sequence identity to <i>K. pneumoniae</i> , strain MRSN 699478 (GenBank: JAGYBV010000073.1) ⁴

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

30 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 699654**

Catalog No. NR-55591

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 699654 was isolated in 2019 from a human urine sample in North America as part of a global surveillance program. NR-55591 was deposited as a susceptible strain, sensitive to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, ciprofloxacin, ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tobramycin, tetracycline, tigecycline and trimethoprim/sulfamethoxazole. NR-55591 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70051620

Manufacturing Date: 14APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive	Sensitive (≤ 2 µg/mL) Sensitive (4 to 8 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.19 µg/mL) Sensitive (0.125 to 0.19 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.032 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.25 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 0.5 µg/mL) ³ Sensitive (≤ 1 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 699654 (GenBank: JAGYBU01000091.1)	99.3% sequence identity to <i>K. pneumoniae</i> , strain MRSN 699654 (GenBank: JAGYBU01000091.1) ⁴

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK®2 GN74 and ETEST®.

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

01 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 702261**

Catalog No. NR-55592

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 702261 was isolated in 2019 from a human wound sample in Africa as part of a global surveillance program. NR-55592 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam and tigecycline, intermediately resistant to ceftazidime and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftriaxone, tetracycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55592 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70051622

Manufacturing Date: 14APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Intermediate Sensitive Sensitive Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Resistant Sensitive Resistant Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (16 µg/mL) Intermediate (4 µg/mL) ³ Intermediate (6 to 8 µg/mL) Sensitive (0.25 µg/mL) Sensitive (0.25 to 0.38 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.094 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.19 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (16 µg/mL) Resistant (≥ 16 µg/mL) Inconclusive ^{4,5} Intermediate (8 µg/mL) ⁶ Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 702261 (GenBank: JAGYBT010000080.1)	99.7% sequence identity to <i>K. pneumoniae</i> , strain MRSN 702261 (GenBank: JAGYBT010000080.1) ⁷

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 702261 was deposited as resistant to cefepime, but showed a MIC of 4 µg per mL (interpreted as intermediately resistant) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵*K. pneumoniae*, strain MRSN 702261 was deposited as being sensitive to tigecycline. Antibiotic susceptibility testing performed in duplicate determined that for strain MRSN 702261, the tigecycline MICs are 1 µg per mL and 2 µg per mL, which are interpreted as sensitive and resistant, respectively.

⁶The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁷Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

11 AUG 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 702325**

Catalog No. NR-55593

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 702325 was isolated in 2019 from a human wound sample in Africa as part of a global surveillance program. NR-55593 was deposited as a multidrug-resistant strain, sensitive to amikacin, aztreonam, cefepime, ceftazidime, ceftriaxone, ciprofloxacin, ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, tobramycin, tigecycline and trimethoprim/sulfamethoxazole, intermediately resistant to piperacillin/tazobactam and resistant to ampicillin/sulbactam and tetracycline. NR-55593 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70051624

Manufacturing Date: 14APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, low convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Intermediate Resistant Sensitive Sensitive Sensitive	Sensitive (≤ 2 µg/mL) Sensitive (6 µg/mL) ³ Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.19 µg/mL) Sensitive (0.125 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.032 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.19 to 0.25 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (1.5 µg/mL) ⁴ Resistant (≥ 16 µg/mL) Sensitive (1 µg/mL) ⁵ Sensitive (≤ 1 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 702325 (GenBank: JAGYBS010000047.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 702325 (GenBank: JAGYBS010000047.1) ⁶

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 702325 was deposited as resistant to ampicillin/sulbactam, but showed a MIC of 6 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴*K. pneumoniae*, strain MRSN 702325 was deposited as intermediately resistant to piperacillin/tazobactam, but showed a MIC of 1.5 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

30 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 728987**

Catalog No. NR-55594

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 728987 was isolated in 2018 from a human wound sample in Asia as part of a global surveillance program. NR-55594 was deposited as a multidrug-resistant strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ciprofloxacin, ertapenem, imipenem, levofloxacin, meropenem, piperacillin/tazobactam and tigecycline, intermediately resistant to tobramycin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, gentamicin, tetracycline and trimethoprim/sulfamethoxazole. NR-55594 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051626

Manufacturing Date: 01APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® MS (MALDI-TOF)	Gram-negative rods Report results Report results <i>K. pneumoniae</i>	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99.9%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Sensitive Sensitive Resistant Sensitive Sensitive Sensitive Resistant Sensitive Sensitive Resistant Sensitive Intermediate Resistant	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Resistant (16 µg/mL) Resistant (≥ 64 µg/mL) Intermediate (8 µg/mL) ³ Sensitive (0.25 µg/mL) Sensitive (0.38 to 0.5 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.38 to 0.5 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.25 µg/mL) Sensitive (1 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Resistant (≥ 16 µg/mL) Resistant (2 µg/mL) ^{3,4} Intermediate (8 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	$\geq 99\%$ sequence identity to <i>K. pneumoniae</i> , strain MRSN 728987 (GenBank: JAGYBR010000097.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 728987 (GenBank: JAGYBR010000097.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

12 AUG 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.




TEST	SPECIFICATIONS	RESULTS
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

20 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 731029**

Catalog No. NR-55596

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 731029 was isolated in 2019 from a human sample in Asia as part of a global surveillance program. NR-55596 was deposited as a multidrug-resistant strain (MDR) sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, imipenem, meropenem, piperacillin/tazobactam, tigecycline, tobramycin, and trimethoprim/sulfamethoxazole, intermediately resistant to ciprofloxacin and levofloxacin and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, gentamicin, and tetracycline. NR-55596 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051630

Manufacturing Date: 06APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Sensitive Sensitive Resistant Intermediate Sensitive Resistant Sensitive Intermediate Sensitive Sensitive Resistant Sensitive Sensitive Sensitive	Sensitive (≤ 2 µg/mL) Resistant (≥ 32 µg/mL) Sensitive (≤ 1 µg/mL) ³ Sensitive (≤ 1 µg/mL) ³ Sensitive (≤ 1 µg/mL) ³ Sensitive (0.125 µg/mL) Sensitive (0.19 µg/mL) Resistant (32 to 64 µg/mL) Intermediate (1.5 µg/mL) Sensitive (≤ 0.5 µg/mL) Resistant (≥ 16 µg/mL) Sensitive (0.19 to 0.25 µg/mL) Intermediate (4 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (2 µg/mL) ^{4,5} Sensitive (2 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1460 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 731029 (GenBank: JAGYBP010000103.1)	99.3% sequence identity to <i>K. pneumoniae</i> , strain MRSN 731029 (GenBank: JAGYBP010000103.1)

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 731029 was deposited as resistant to aztreonam, cefepime, and ceftazidime but showed a MIC of ≤ 1 µg per mL (interpreted as sensitive) for these antibiotics during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

27 JUN 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 736213**

Catalog No. NR-55597

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 736213 was isolated in 2019 from a human urine sample in North America as part of a global surveillance program. NR-55597 was deposited as a susceptible strain, sensitive to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55597 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051632

Manufacturing Date: 01APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive	Sensitive (≤ 2 µg/mL) Sensitive (4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.19 µg/mL) Sensitive (0.125 to 0.25 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.032 to 0.047 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.25 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 0.5 µg/mL) ³ Sensitive (≤ 1 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1460 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 736213 (GenBank: JAGYBO010000109.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 736213 (GenBank: JAGYBO010000109.1)

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Lead Technical Writer or designee, ATCC Federal Solutions

24 JUN 2022

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 740795**

Catalog No. NR-55598

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 740795 was isolated in 2019 from a human respiratory sample in Europe as part of a global surveillance program. NR-55598 was deposited as a multidrug-resistant strain, sensitive to amikacin, gentamicin, ceftazidime/avibactam, trimethoprim/sulfamethoxazole, tetracycline and tigecycline and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftriaxone, ceftolozane/tazobactam, ciprofloxacin, ertapenem, imipenem, levofloxacin, meropenem, piperacillin/tazobactam and tobramycin. NR-55598 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051634

Manufacturing Date: 01APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Sensitive Sensitive Resistant Sensitive	Sensitive (8 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (0.5 µg/mL) Resistant (12 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 1 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Sensitive (4 µg/mL) Resistant (2 µg/mL) ^{3,4} Resistant (≥ 16 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 740795 (GenBank: JAGYBN01000093.1)	99.3% sequence identity to <i>K. pneumoniae</i> , strain MRSN 740795 (GenBank: JAGYBN01000093.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

15 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

28 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 750877**

Catalog No. NR-55600

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 750877 was isolated in 2020 from a human urine sample in North America as part of a global surveillance program. NR-55600 was deposited as a susceptible strain, sensitive to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline, tobramycin and trimethoprim/sulfamethoxazole. NR-55600 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051638

Manufacturing Date: 01APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream Non-motile <i>K. pneumoniae</i> (89%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive	Sensitive (≤ 2 µg/mL) Sensitive (4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.19 µg/mL) Sensitive (0.19 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.094 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.38 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (1 µg/mL) ³ Sensitive (≤ 1 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 750877 (GenBank: JAGYBL010000154.1)	99.7% sequence identity to <i>K. pneumoniae</i> , strain MRSN 750877 (GenBank: JAGYBL010000154.1) ⁴

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁴Also consistent with other *Klebsiella* species

/Sonia Bjorum Brower/

Sonia Bjorum Brower

12 AUG 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 750999**

Catalog No. NR-55601

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 750999 was isolated in 2020 from a human urine sample in North America as part of a global surveillance program. NR-55601 was deposited as a susceptible strain, sensitive to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline, tobramycin, and trimethoprim/sulfamethoxazole. NR-55601 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051640

Manufacturing Date: 06APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream Non-motile <i>K. pneumoniae</i> (98%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive Sensitive	Sensitive (≤ 2 µg/mL) Sensitive (4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.125 µg/mL) Sensitive (0.125 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.047 µg/mL) Sensitive (≤ 0.5 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (0.25 to 0.38 µg/mL) Sensitive (≤ 0.12 µg/mL) Sensitive (≤ 0.25 µg/mL) Sensitive (≤ 4 µg/mL) Sensitive (≤ 1 µg/mL) Sensitive (≤ 0.5 µg/mL) ³ Sensitive (≤ 1 µg/mL) Sensitive (≤ 20 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 750999 (GenBank: JAGYBK010000105.1)	99.6% sequence identity to <i>K. pneumoniae</i> , strain MRSN 750999 (GenBank: JAGYBK010000105.1)

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
 Sonia Bjorum Brower

20 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 752729**

Catalog No. NR-55602

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 752729 was isolated in 2019 from a human perianal sample in Europe as part of a global surveillance program. NR-55602 was deposited as an extensively drug-resistant strain, sensitive to amikacin, gentamicin, tetracycline and tigecycline and resistant to ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tobramycin and trimethoprim/sulfamethoxazole. NR-55602 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051642

Manufacturing Date: 06APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Sensitive Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Resistant Resistant Sensitive Sensitive Resistant Resistant	Sensitive (16 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 64 µg/mL) Resistant (> 256 µg/mL) Resistant (> 256 µg/mL) Resistant (≥ 64 µg/mL) Resistant (> 32 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 1 µg/mL) Resistant (> 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Sensitive (4 µg/mL) Resistant (2 µg/mL) ^{3,4} Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 752729 (GenBank: JAGYBJ010000100.1)	99.5% sequence identity to <i>K. pneumoniae</i> , strain MRSN 752729 (GenBank: JAGYBJ010000100.1) ⁵

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)


²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species

Figure 1: Colony Morphology



/Sonia Bjorum Brower/
Sonia Bjorum Brower

12 AUG 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.



***Klebsiella pneumoniae*, Strain MRSN 761403**

Catalog No. NR-55603

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (*K. pneumoniae*), strain MRSN 761403 was isolated in 2020 from a blood sample in North America as part of a global surveillance program. NR-55603 was deposited as an extensively drug-resistant strain, sensitive to ceftazidime/avibactam, gentamicin and tigecycline, intermediately resistant to tetracycline and resistant to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tobramycin and trimethoprim/sulfamethoxazole. NR-55603 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was passaged once in Tryptic Soy broth for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051644

Manufacturing Date: 06APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® 2 (GN card)	Gram-negative rods Report results Report results <i>K. pneumoniae</i> (≥ 89%)	Gram-negative rods Circular, convex, entire, smooth, mucoid and cream Non-motile <i>K. pneumoniae</i> (99%)
Antibiotic Susceptibility Profile^{1,2} Amikacin Ampicillin/sulbactam Aztreonam Cefepime Ceftazidime Ceftazidime/avibactam Ceftolozane/tazobactam Ceftriaxone Ciprofloxacin Ertapenem Gentamicin Imipenem Levofloxacin Meropenem Piperacillin/tazobactam Tetracycline Tigecycline Tobramycin Trimethoprim/sulfamethoxazole	Resistant Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Resistant Sensitive Resistant Resistant Resistant Resistant Resistant Intermediate Sensitive Resistant Resistant	Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 64 µg/mL) Resistant (16 µg/mL) Resistant (≥ 64 µg/mL) Sensitive (1 to 1.5 µg/mL) Resistant (16 µg/mL) Resistant (≥ 64 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Sensitive (≤ 1 µg/mL) Resistant (≥ 32 µg/mL) Resistant (≥ 8 µg/mL) Resistant (≥ 16 µg/mL) Resistant (≥ 128 µg/mL) Sensitive (3 to 4 µg/mL) ³ Resistant (1.5 to 2 µg/mL) ^{4,5} Resistant (≥ 16 µg/mL) Resistant (≥ 320 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1460 base pairs)	≥ 99% sequence identity to <i>K. pneumoniae</i> , strain MRSN 761403 (GenBank: JAGYBI010000093.1)	99.4% sequence identity to <i>K. pneumoniae</i> , strain MRSN 761403 (GenBank: JAGYBI010000093.1) ⁶

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK[®]2 GN74 and ETEST[®].

³*K. pneumoniae*, strain MRSN 761403 was deposited as intermediately resistant to tetracycline, but showed a MIC of 3 to 4 µg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

⁶Also consistent with other *Klebsiella* species

/Sonia Bjorum Brower/

Sonia Bjorum Brower

18 JUL 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

