SUPPORTING INFECTIOUS DISEASE RESEARCH

## Trypanosoma brucei subsp. gambiense, Strain STIB 386 (in vitro)

#### Catalog No. NR-44389

#### **Product Description:**

*Trypanosoma brucei (T. brucei)* subsp. *brucei*, strain STIB 386 *(in vitro)* was harvested from the blood of infected BALB/c mice and adapted to cell culture by BEI Resources. The parent strain STIB 386 (BEI Resources NR-36198) was derived from strain TH 114/78E (020), which was isolated in 1978 from a male patient in Koudougou, Ivory Coast, West Africa. NR-44389 lot 70022602 was produced by cultivation of BEI Resources seed lot 62069847 in SDM-79 medium (Life Technologies, custom order part number ME090164 P1) supplemented with 10% (v/v) heat-inactivated fetal bovine serum (HIFBS) and 7.5 µg/mL hemin for 5 days at 27°C in an aerobic atmosphere.

#### Lot: 70022602

### Manufacturing Date: 30JAN2019

SPECIFICATIONS	RESULTS
Report results	Elongated and refractive; rosettes visible
≥ 98% sequence identity to <i>T. brucei</i> subsp. <i>gambiense</i> , strain DA1972 (GenBank: AF306774.1)	98.2% sequence identity to <i>T. brucei</i> subsp. <i>gambiense</i> , strain DA1972 (GenBank: AF306774.1) <sup>3</sup>
No amplicon	No amplicon
> 10 <sup>6</sup> cells/mL	1 × 10 <sup>8</sup> parasites/mL
Growth	Growth
No growth No growth No growth No growth No growth	No growth No growth No growth No growth No growth No growth No growth
	Report results ≥ 98% sequence identity to <i>T. brucei</i> subsp. <i>gambiense</i> , strain DA1972 (GenBank: AF306774.1) No amplicon > 10 <sup>6</sup> cells/mL Growth No growth No growth

<sup>1</sup>Testing completed on vialed, post-freeze material.

<sup>2</sup>Testing completed on bulk material prior to vialing and freezing.

<sup>3</sup>Also consistent with *T. evansi* and/or *T. equiperdum*, which are putative subspecies of *T. brucei* (Lun, Z. R., et al. "*Trypanosoma brucei*: Two Steps to Spread Out from Africa." <u>Trends Parasitol.</u> 26 (2010): 424-427. PubMed: 20561822.)

<sup>4</sup>T. brucei subsp. gambiense is differentiated from *T. brucei* subsp. rhodesiense by lack of the SRA gene. (Radwanska, M., et al. "The Serum Resistance-Associated Gene as a Diagnostic Tool for the Detection of *Trypanosoma brucei rhodesiense*." <u>Am. J. Trop. Med. Hyg.</u> 67 (2002): 684-690. PubMed: 12518862.).

<sup>5</sup>Atlas, Ronald M. <u>Handbook of Microbiological Media</u>. 3rd ed. Ed. Lawrence C. Parks. Boca Raton: CRC Press, 2004, p. 798.

# /Heather Couch/

## Heather Couch

Program Manager or designee, ATCC Federal Solutions

ATCC<sup>®</sup>, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC<sup>®</sup>'s knowledge.

ATCC<sup>®</sup> is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.



19 APR 2022

BEI Resources www.beiresources.org E-mail: <u>contact@beiresources.org</u> Tel: 800-359-7370 Fax: 703-365-2898