

## **Certificate of Analysis for NR-50345**

### Genomic RNA from Chikungunya Virus, 181/25

#### Catalog No. NR-50345

#### **Product Description:**


Genomic RNA was isolated from a preparation of clarified supernatant from *Chlorocebus* (previously *Cercopithecus*) *aethiops* kidney epithelial cells (Vero; ATCC® CCL-81™) infected with chikungunya virus (CHIKV), 181/25 (BEI Resources lot 70050553) using QIAamp® Viral RNA Mini Kit (Qiagen® 52906). The viral genomic RNA is in a background of cellular nucleic acid and carrier RNA.

Lot: 70058080 Manufacturing Date: 18JAN2023

TEST	SPECIFICATIONS	RESULTS
Genotypic Analysis Sequencing of species-specific region (~ 1060 nucleotides)	≥ 98% identity with CHIKV, 181/25 (GenBank: MW473668.1)	100% identity with CHIKV, 181/25 (GenBank: MW473668.1)
Functional Activity by RT-PCR Amplification <sup>1</sup> NSP gene	~ 1180 base pair amplicon	~ 1180 base pair amplicon (Figure 1)
Estimated Concentration (post-dilution) by RiboGreen® Measurement (Viral, Cellular and Carrier)²	Report results	24.5 ng per 100 μL (0.25 μg/mL)
Estimated Amount per Vial <sup>2</sup>	Report results	24.5 ng
Genome Copy Number Using BioRad QX200 Droplet Digital PCR (ddPCR™) System (Post vial; 18 replicates)	Report results	1.4 × 10 <sup>8</sup> NDU/mL <sup>3</sup>
Virus Inactivation 10% of total yield inoculated on Vero cells and evaluated for cytopathic effect and by RT-PCR after serial passage <sup>4</sup>	No viable virus detected	No viable virus detected

<sup>&</sup>lt;sup>1</sup>Amplified using iTaq<sup>™</sup> Universal SYBR Green One-Step Kit (Bio-Rad<sup>®</sup> 172-5151) with 5 μL of NR-50345 in a 50 μL reaction

Figure 1: Functional Activity of NR-50345 by RT-PCR Amplification



Lane 1: Invitrogen™ TrackIt™ 1 Kb Plus DNA Ladder

Lane 2: PCR product from 1 µL of NR-50345

BEI Resources www.beiresources.org E-mail: contact@beiresources.org Tel: 800-359-7370

Tel: 800-359-7370 Fax: 703-365-2898

<sup>&</sup>lt;sup>2</sup>Measurement is determined pre-vial prior to dilution due to the limit of detection of the quantification method

<sup>&</sup>lt;sup>3</sup>NDU; NAAT(Nucleic acid amplification testing) Detectable Unit

<sup>&</sup>lt;sup>4</sup>Use of the QIAamp® Viral RNA Mini Kit has been demonstrated to consistently inactivate 100% of CHIKV as shown by the absence of cytopathic effect (CPE) and viral RNA expression by RT-PCR after plating the entire extract on virus-susceptible cells for two passages.



# Certificate of Analysis for NR-50345

/Sonia Bjorum Brower/ Sonia Bjorum Brower

16 MAY 2023

Technical Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898