

Certificate of Analysis for NR-52241

Acinetobacter baumannii, Strain MRSN 351162

Catalog No. NR-52241

This reagent is the tangible property of the U.S. Government.

Product Description:

Acinetobacter baumannii (A. baumannii), strain MRSN 351162 was isolated in 2011 from a human in Europe as part of a global surveillance program. A. baumannii, strain MRSN 351162 was deposited as sensitive to amikacin, ampicillin/sulbactam, cefepime, ceftazidime, ceftriaxone, ciprofloxacin, colistin, imipenem, levofloxacin, meropenem, tetracycline, tobramycin and trimethoprim/sulfamethoxazole and intermediately resistant to gentamicin. NR-52241 was produced by inoculation of BEI Resources seed lot 70042863 into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70061815 Manufacturing Date: 29JUN2023

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology	Report results	Circular, convex, entire, smooth and cream
Growth at 44°C ± 2°C ¹ 1 day in an aerobic atmosphere on Tryptic Soy agar	Growth	Growth
Motility Hardy Diagnostics™ Motility Test Medium w/TTC Indicator for 1 day at 37°C in an aerobic atmosphere	Report results	Non-motile
VITEK® MS (MALDI-TOF)	A. baumannii	A. baumannii (99.9%)
Antibiotic Susceptibility Profile ^{2,3}		
Amikacin	Sensitive	Sensitive (3 µg/mL)
Ampicillin/sulbactam	Sensitive	Sensitive (1.5 µg/mL)
Cefepime	Sensitive	Sensitive (2 µg/mL)
Ceftazidime	Sensitive	Sensitive (4 µg/mL)
Ceftriaxone	Intermediate	Intermediate (16 µg/mL) ⁴
Ciprofloxacin	Sensitive	Sensitive (≤ 0.25 µg/mL)
Gentamicin	Sensitive	Sensitive (≤ 1 µg/mL) ⁵
Imipenem	Sensitive	Sensitive (0.25 µg/mL)
Levofloxacin	Sensitive	Sensitive (≤ 0.12 µg/mL)
Meropenem	Sensitive	Sensitive (≤ 0.25 µg/mL)
Tetracycline	Sensitive	Sensitive (≤ 1 µg/mL)
Tobramycin	Sensitive	Sensitive (≤ 1 µg/mL)
Trimethoprim/sulfamethoxazole	Sensitive	Sensitive (≤ 20 µg/mL)
Genotypic Analysis		, ,
Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to A. baumannii, strain MRSN 351162 (GenBank: VHEX01000060.1)	100% sequence identity to A. baumannii, strain MRSN 351162 (GenBank: VHEX01000060.1)
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898



Certificate of Analysis for NR-52241

TEST	SPECIFICATIONS	RESULTS
Viability	Growth	Growth

¹Growth at 44°C differentiates A. baumannii from A. calcoaceticus and A. pittii, which do not grow at 44°C.

/Sonia Bjorum Brower/ Sonia Bjorum Brower

07 SEP 2023

Technical Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org Tel: 800-359-7370

Fax: 703-365-2898

²Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

³Antibiotic susceptibility was tested using a combination of VITEK®2 GN81 and E-test strips.

⁴A. baumannii, strain MRSN 351162 was deposited as being sensitive to ceftriaxone, but showed a MIC of 12 to 16 μg/mL (interpreted as intermediately resistant) for lot 70042862 during QC testing. Testing was performed in duplicate.

⁵A. baumannii, strain MRSN 351162 was deposited as being intermediately resistant to gentamicin, but showed a MIC of ≤ 2 μg/mL (interpreted as sensitive) for lot 70042862 during QC testing. Testing was performed in duplicate.