

Certificate of Analysis for NR-52225

Acinetobacter baumannii, Strain MRSN 31937

Catalog No. NR-52225

This reagent is the tangible property of the U.S. Government.

Product Description:

Acinetobacter baumannii (A. baumannii), strain MRSN 31937 was isolated in 2004 from a human wound in Europe as part of a global surveillance program. A. baumannii, strain MRSN 31937 was deposited as sensitive to amikacin, ampicillin/sulbactam, cefepime, ceftazidime, ceftriaxone, ciprofloxacin, colistin, gentamicin, imipenem, levofloxacin, meropenem, tetracycline, tobramycin and trimethoprim/sulfamethoxazole. NR-52225 was produced by inoculation of BEI Resources seed lot 70042847 into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70061534 Manufacturing Date: 21JUN2023

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology	Report results	Circular, convex, entire, smooth and cream (Figure 1)
Growth at 44°C ± 2°C¹ 1 day in an aerobic atmosphere on Tryptic Soy agar	Growth	Growth
Motility Hardy Diagnostics™ Motility Test Medium with TTC Indicator for 1 day at 37°C in an aerobic atmosphere	Report results	Non-motile
VITEK® MS (MALDI-TOF)	A. baumannii	A. baumannii (99.9%)
Antibiotic Susceptibility Profile ^{2,3}		
Amikacin	Sensitive	Sensitive (1.2 µg/mL)
Ampicillin/sulbactam	Sensitive	Sensitive (1 µg/mL)
Cefepime	Sensitive	Sensitive (1 µg/mL)
Ceftazidime	Sensitive	Sensitive (4 µg/mL)
Ceftriaxone	Sensitive	Sensitive (8 µg/mL)
Ciprofloxacin	Sensitive	Sensitive (≤ 0.25 μg/mL)
Gentamicin	Sensitive	Sensitive (≤ 1 µg/mL)
Imipenem	Sensitive	Sensitive (0.19 µg/mL)
Levofloxacin	Sensitive	Sensitive (≤ 0.12 µg/mL)
Meropenem	Sensitive	Sensitive (≤ 0.25 µg/mL)
Tetracycline	Sensitive	Sensitive (≤ 1 µg/mL)
Tobramycin	Sensitive	Sensitive (≤ 1 µg/mL)
Trimethoprim/sulfamethoxazole	Sensitive	Sensitive (≤ 20 µg/mL)
Genotypic Analysis		, , ,
Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to A. baumannii, strain MRSN 31937 (GenBank: VHFP01000152.1)	99.7% sequence identity to A. baumannii, strain MRSN 31937 (GenBank: VHFP01000152.1)
Purity	Growth consistent with expected	Growth consistent with expected
7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	colony morphology	colony morphology

BEI Resources

www.beiresources.org

E-mail: contact@beiresources.org
Tel: 800-359-7370
Fax: 703-365-2898



Certificate of Analysis for NR-52225

TEST	SPECIFICATIONS	RESULTS
Viability	Growth	Growth

¹Growth at 44°C differentiates *A. baumannii* from *A. calcoaceticus and A. pittii*, which do not grow at 44°C.

Figure 1: Colony Morphology



/Sonia Bjorum Brower/ Sonia Bjorum Brower

06 OCT 2023

Technical Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898

²Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

³Antibiotic susceptibility was tested using a combination of VITEK[®]2 GN81 and E-test strips.