

**Pan-*Rickettsia* Quantitative PCR (qPCR) Detection Assay**

**Catalog No. NR-35520**

This reagent is the tangible property of the U.S. Government.

**Product Description:**

NR-35220 is a Pan-*Rickettsia* Quantitative PCR (qPCR) Detection Assay to detect and quantitate *Rickettsia* species. The assay consists of the following components:

Catalog Number	Component	Description
NR-35354	Forward Primer	The forward primer pairs with NR-35355 to amplify a sequence of 135 nucleotides from the 16S ribosomal RNA gene of <i>Rickettsia</i> species.
NR-35355	Reverse Primer	The reverse primer pairs with NR-35354 to amplify a sequence of 135 nucleotides from the 16S ribosomal RNA gene of <i>Rickettsia</i> species.
NR-35356	Probe	The probe is designed with 6-carboxyfluorescein (6-FAM) at the 5' end and both a minor groove binding moiety (MGB) and a non-fluorescent quenching dye (NFQ) at the 3' end.
NR-35519	Plasmid Standard	The plasmid standard contains <i>Rickettsia prowazekii</i> ribosomal RNA (rRNA) gene sequences.


**Lot: 70059015**

**Manufacturing Date: 24APR2023**

TEST	SPECIFICATIONS	RESULTS
<b>Genome Copy Number Using BioRad QX200 Droplet Digital PCR (ddPCR™) System</b> (Pre-vial; 9 replicates)	Report results	6.4 × 10 <sup>12</sup> genome copies/mL
<b>Quantitative PCR – Standard Curve (Figure 1)</b> Correlation coefficient PCR efficiency Dilution separations (C <sub>t</sub> values) Quantitative sensitivity	≥ 0.98 85 to 110% ~ 3.33 cycles Report results	0.996 88.8% ~ 3.62 cycles 10 molecules/reaction
<b>Concentration (µM)</b> NR-35354 Lot 70059016 (forward primer) NR-35355 Lot 70059017 (reverse primer) NR-35356 Lot 70059018 (probe)	Report results Report results Report results	10 10 5
<b>Plasmid-Based Standard (NR-35519 Lot 70059019)<sup>1</sup></b> Agarose Gel Electrophoresis of Linearized Plasmid DNA <sup>1</sup> Concentration of DNA Molecules	Single band at ~ 4,100 base pairs Report results	Single band at ~ 4,100 base pairs 2.0 × 10 <sup>11</sup> molecules/mL (1.0 × 10 <sup>9</sup> molecules/5 µL)

<sup>1</sup>Plasmid DNA was extracted using a Plasmid Plus Maxi Kit (QIAGEN® 12963). Purified plasmid DNA was linearized with EcoRI-HF (New England BioLabs, Inc. R3101).

Figure 1: Quantitative PCR Cycle Graph (A) and Associated Standard Curve (B)



/Sonia Bjorum Brower/

Sonia Bjorum Brower

Technical Manager or designee, ATCC Federal Solutions

26 JUL 2024

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.  
 You are authorized to use this product for research use only. It is not intended for human use.

