

Certificate of Analysis for NR-10383

SUPPORTING INFECTIOUS DISEASE RESEARCH

Genomic RNA from Kilbourne F161: A/Philippines/2/82 (HA, NA) x A/Puerto Rico/8/34 (H3N2), Reassortant X-79 (Mouse-adapted)

Catalog No. NR-10383

Product Description: Genomic RNA was isolated from a preparation of pooled allantoic fluid from specific-pathogen free embryonated chicken eggs¹ infected with mouse-adapted reassortant influenza A virus, A/Philippines/2/82 (HA, NA) x A/Puerto Rico/8/34 (H3N2) (Kilbourne F161; X-79).

Lot²: 58776486 Manufacturing Date: 15MAR2010

TEST	SPECIFICATIONS	RESULTS
Sequencing of Influenza A Specific Region Matrix gene (915 nucleotides)	Influenza A virus	Influenza A virus
Total RNA Content by RiboGreen [®] Measurement (Viral, Cellular and Carrier)	Report results	12 ng per 100 μL
Functional Activity by RT-PCR Amplification ^{3,4}	~ 1030 bp amplicon	~ 1030 bp amplicon (See Figure 1)
Virus Inactivation 10% of total yield incubated in SPF embryonated chicken eggs ¹ and evaluated by hemagglutination assay ⁵	No virus detected	No virus detected
Sodium Azide Content	Report results	0.0004%

¹⁰ to 11-day-old SPF Fertile Chicken Eggs acquired from B&E Eggs, York Springs, Pennsylvania

Date: 02 FEB 2011 Signature: Ocrothy C. Young

Title: Technical Manager, BEI Authentication or designee

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

Biodefense and Emerging Infections Research Resources Repository www.beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

E-mail: contact@beiresources.org


²Nucleic acid was extracted from a preparation of mouse-adapted reassortant influenza A virus, A/Philippines/2/82 (HA, NA) x A/Puerto Rico/8/34 (H3N2) (Kilbourne F161; X-79) (BEI Resources NR-3668; Lot 58312186) using a QIAamp[®] Viral RNA Mini kit (Qiagen 52906).

³Amplified using 1:10 dilution of NR-10376 and One-Step RT-PCR Kit (Qiagen 210212)

⁴The primers are described in Hoffmann, E., et al. "Universal Primer Set for the Full-Length Amplification of All Influenza A Viruses." <u>Arch. Virol.</u> 146 (2001): 2275-2289. PubMed: 11811679.

⁵This extraction procedure has been shown to consistently inactivate 100% of influenza viruses.

Certificate of Analysis for NR-10383



Lane 1: 10-fold dilution of NR-10383

Lane 2: Invitrogen™ TrackIt™ 1 Kb Plus DNA Ladder

NR-10383_58776486_02FEB2011

E-mail: contact@beiresources.org