

Certificate of Analysis for NR-4344

Salmonella enterica subsp. enterica, Strain LS1029

Catalog No. NR-4344


Product Description: Salmonella enterica (S. enterica) subsp. enterica, strain LS1029 is a derivative of strain 14028. Strain 14028 (also referred to as CDC 6516-60) was isolated from chickens and deposited to the ATCC[®] in 1960. Strain LS1029 expresses STM3117, a putative glutathione-dependent glyoxalase which may have a role in neutralizing methylglyoxal, a compound lethal to Salmonella.

Lot¹: 59586819 Manufacturing Date: 24NOV2010

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rod	Gram-negative rod
Colony morphology ²	Report results	Circular, convex, entire, and translucent (Figure 1)
Analytical profile index (API [®] 20 E) ³	Consistent with Salmonella spp.	Consistent with Salmonella spp.
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene (~ 1470 bp)	Consistent with S. enterica	Consistent with S. enterica
Serogrouping	Report results	Serogroup B ⁴
Viability (post-freeze) ²	Growth	Growth

¹NR-4344 was produced by inoculation of the deposited material into Tryptic Soy Broth and grown 24 hours in an aerobic atmosphere. Broth inoculum was added to Kolles which was grown 24 hours at 37°C and aerobic atmosphere to produce this lot.

Figure 1



Date: 13 SEP 2011 **Signature:**

www.beiresources.org

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

²24 hours at 37°C and aerobic atmosphere on Tryptic Soy Agar

³Additional serological tests are required to determine Salmonella species using API[®] 20 E

⁴Serogroup B contains serotype Typhimurium in addition to other serotypes