

Certificate of Analysis for HM-210D

Genomic DNA from Bacteroides eggerthii, Strain 1_2_48FAA

Catalog No. HM-210D

Product Description: Genomic DNA was obtained from a preparation of *Bacteroides eggerthii* (*B. eggerthii*), strain 1_2_48FAA.

Lot^{1,2}: 59925940 Manufacturing Date: 02MAR2012

TEST	SPECIFICATIONS	RESULTS
Sequencing of 16S Ribosomal RNA Gene (~ 1470 base pairs)	≥ 99% identical to GenBank: ACWG01000065 (<i>B. eggerthii</i> , strain 1_2_48FAA)	≥ 99% identical to GenBank: ACWG01000065 (<i>B. eggerthii</i> , strain 1_2_48FAA)
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by PicoGreen® Measurement	0.7 to 1.5 μg in 25 to 100 μL per vial	1.0 μg in 37 μL per vial (28 μg/mL)
Functional Activity by PCR Amplification 16S ribosomal RNA gene	~ 1500 bp amplicon	~ 1500 bp amplicon
OD ₂₆₀ /OD ₂₈₀ Ratio	1.7 to 2.0	1.9
Bacterial Inactivation 10% of total yield plated on Tryptic Soy Agar with 5% defibrinated sheep blood ^{3,4}	No viable bacteria detected	No viable bacteria detected

¹Quality control of HMP organisms used for DNA extraction is only performed to demonstrate that the material produced by BEI Resources is identical to the deposited material. It should not be considered a complete characterization of the deposited organism.

Date: 23 JUL 2012 Signature:

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

BEI Resources

www.beiresources.org

You are authorized to use this product for research use only. It is not intended for human use.

E-mail: contact@beiresources.org Tel: 800-359-7370

Fax: 703-365-2898


²The bacterial preparation used for extraction of genomic DNA was produced by culture of BEI Resources HM-210 (Lot 60685355). Genomic DNA was extracted using proprietary technology.

³7 days at 37°C and anaerobic atmosphere (80% N₂:10% CO₂:10% H₂)

⁴An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-positive and Gram-negative bacteria.



Certificate of Analysis for HM-210D



Lane 1: Invitrogen™ TrackIt 1 Kb Plus DNA Ladder™

Lane 2: 200 ng of HM-210D

E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898