

## **Certificate of Analysis for HM-242D**

## Genomic DNA from Neisseria mucosa, Strain C102

Catalog No. HM-242D

**Product Description:** Genomic DNA was obtained from a preparation of *Neisseria mucosa* (*N. mucosa*), strain C102.

Lot<sup>1,2</sup>: 61012590 Manufacturing Date: 25MAY2012

TEST	SPECIFICATIONS	RESULTS
Sequencing of 16S Ribosomal RNA Gene (~ 1490 base pairs)	≥ 99% identical to GenBank: ACRG01000007 ( <i>N. mucosa</i> , strain C102)	≥ 99% identical to GenBank: ACRG01000007 ( <i>N. mucosa</i> , strain C102)
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by PicoGreen® Measurement	0.7 to 1.5 μg in 25 to 100 μL per vial	0.8 μg in 43 μL per vial (19 μg/mL)
Functional Activity by PCR Amplification 16S ribosomal RNA gene	~ 1500 bp amplicon	~ 1500 bp amplicon
OD <sub>260</sub> /OD <sub>280</sub> Ratio	1.7 to 2.0	2.0
Bacterial Inactivation 10% of total yield plated on Chocolate Agar (GC Medium) <sup>3,4</sup>	No viable bacteria detected	No viable bacteria detected

<sup>&</sup>lt;sup>1</sup>Quality control of HMP organisms used for DNA extraction is only performed to demonstrate that the material produced by BEI Resources is identical to the deposited material. It should not be considered a complete characterization of the deposited organism.

Date: 27 SEP 2012 Signature:

**Title:** Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

**BEI Resources** 


<sup>&</sup>lt;sup>2</sup>The bacterial preparation used for extraction of genomic DNA was produced by culture of BEI Resources HM-242 (Lot 60609341). Genomic DNA was extracted using proprietary technology.

<sup>&</sup>lt;sup>3</sup>7 days at 37°C in an aerobic atmosphere with 5% CO<sub>2</sub>

<sup>&</sup>lt;sup>4</sup>An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-positive and Gram-negative bacteria.



## **Certificate of Analysis for HM-242D**



Lane 1: Invitrogen™ TrackIt 1 Kb Plus DNA Ladder™ Lane 2: 200 ng of HM-242D