

## Certificate of Analysis for NR-28784

## Salmonella enterica subsp. enterica, Strain SARA23 (CDC B1722) (Serovar Saint Paul)

## Catalog No. NR-28784

**Product Description**: Salmonella enterica (S. enterica) subsp. enterica, strain SARA23 (also referred to as strain CDC B1722) was isolated from human stool in Pennsylvania, USA. This strain is reported to be antibiotic-susceptible.

Lot<sup>1</sup>: 61664483 Manufacturing Date: 24APR2013

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology <sup>2</sup>	Report results	Circular, convex, entire, smooth and cream (Figure 1)
Analytical profile index (API <sup>®</sup> 20 E) <sup>3</sup>	Consistent with Salmonella spp.	Consistent with Salmonella spp.
Genotypic Analysis Sequencing of 16S ribosomal RNA gene		
(~ 1460 base pairs)	Consistent with S. enterica	Consistent with S. <i>enterica</i> subsp. <i>enterica</i> <sup>4</sup>
Riboprinter <sup>®</sup> Microbial Characterization System	Consistent with S. enterica	Consistent with S. enterica subsp. enterica serovar Saint Paul
Serogroup Analysis	Report results	Serogroup B <sup>5</sup>
Viability (post-freeze) <sup>2</sup>	Growth	Growth

<sup>&</sup>lt;sup>1</sup>NR-28784 was produced by inoculation of the deposited material into Nutrient broth and incubated 24 hours at 37°C in an aerobic atmosphere. Broth inoculum was added to kolles, which were incubated for 24 hours at 37°C in an aerobic atmosphere to produce this lot. Purity of this lot was assessed for 7 days under propagation conditions.

Figure 1



Date: 06 NOV 2013 Signature:

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

**BEI Resources** 

www.beiresources.org

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

<sup>&</sup>lt;sup>2</sup>24 hours at 37°C in an aerobic atmosphere on Nutrient agar

<sup>&</sup>lt;sup>3</sup>Additional serological tests are required to determine *Salmonella* species using API<sup>®</sup> 20 E.

<sup>&</sup>lt;sup>4</sup>≥ 99% identical to GenBank: ABAM02000001 (*S. enterica* subsp. *enterica*, serovar Saintpaul, strain SARA23)

<sup>&</sup>lt;sup>5</sup>Serogroup B contains serovar Saint Paul in addition to other serovars [Ng, S. P., et al. "Detection and Serogroup Differentiation of Salmonella spp. in Food within 30 Hours by Enrichment-Immunoassay with a T6 Monoclonal Antibody Capture Enzyme-Linked Immunosorbent Assay." <u>Appl.</u> Environ. Microbiol. 62 (1996): 2794-2302. PubMed: 8779567.].