

Certificate of Analysis for NR-43012

Genomic RNA from Influenza A virus, A/duck/Czechoslovakia/1956 (H4N6)

Catalog No. NR-43012

Product Description: Genomic RNA was isolated from a preparation of pooled allantoic fluid from specific pathogen free (SPF) embryonated chicken eggs¹ infected with influenza A virus, A/duck/Czechoslovakia/1956 (H4N6).

Lot²: 62244658 Manufacturing Date: 08JAN2014

TEST	SPECIFICATIONS	RESULTS
Sequencing of Influenza A Specific Region Matrix gene (942 nucleotides)	Consistent with A/duck/Czechoslovakia/1956 (H4N6) Consistent with NR-21652	100% identity with A/duck/Czechoslovakia/1956 (H4N6) (GenBank: CY130023) Consistent with NR-21652
Functional Activity by RT-PCR Amplification ^{3,4}	~ 1025 bp amplicon	~ 1025 bp amplicon (See Figure 1)
Total RNA Content by RiboGreen [®] Measurement (Viral, Cellular and Carrier)	Report results	259 ng per 100 μL
Virus Inactivation 10% of total yield incubated in SPF embryonated chicken eggs ¹ and evaluated by hemagglutination assay ⁵	No virus detected	No virus detected

⁹⁻ to 11-day-old SPF Embryonated Chicken Eggs acquired from B&E Eggs, York Springs, Pennsylvania

Signature: Milhard G. Com la **Date:** 21 MAR 2014

> Technical Manager, BEI Authentication or designee Title:

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org Tel: 800-359-7370

Fax: 703-365-2898

²Nucleic acid was extracted from a preparation of influenza A virus, A/duck/Czechoslovakia/1956 (H4N6) (BEI Resources NR-21652, Lot 59772071), using a QlAamp[®] Viral RNA Mini Kit (Qiagen 52906).


³Amplified using iScript™ One-Step RT-PCR Kit with SYBR[®] Green (Bio-Rad 170-8893) with 5μL of NR-43012 in a 50 μL reaction

⁴The primers are described in Hoffmann, E., et al. "Universal Primer Set for the Full-Length Amplification of All Influenza A Viruses." Arch. Virol. 146 (2001): 2275-2289. PubMed: 11811679.

⁵Use of the QIAamp[®] Viral RNA Mini Kit has been demonstrated to consistently inactivate 100% of influenza viruses.



Certificate of Analysis for NR-43012



Lane 1: Invitrogen™ TrackIt™ 1 Kb Plus DNA Ladder Lane 2: Negative control

Lane 3: NR-43012

E-mail: contact@beiresources.org Tel: 800-359-7370

Fax: 703-365-2898