

***Leptospira interrogans*, Strain HAI0156 (Serovar Copenhageni)**

**Catalog No. NR-19891**

**Product Description:** *Leptospira interrogans* (*L. interrogans*), strain HAI0156 (serovar Copenhageni) is a human isolate obtained from a patient at Hospital de Apoyo in Iquitos, Peru.

**Lot<sup>1</sup>: 62710998**

**Manufacturing Date: 15SEP2014**

TEST	SPECIFICATIONS	RESULTS
<b>Phenotypic Analysis</b> Cellular morphology Colony morphology  Motility (wet mount)	Spirochete Growth below the soft agar surface (Dinger's disk) Motile	Spirochete Growth below the soft agar surface (Dinger's disk) <sup>2</sup> Motile
<b>Genotypic Analysis</b> Sequencing of 16S ribosomal RNA gene (~ 850 base pairs)	Consistent with <i>L. interrogans</i>	Consistent with <i>L. interrogans</i> <sup>3</sup>
<b>Purity (post-freeze)<sup>4</sup></b>	No growth	Growth observed <sup>5</sup>
<b>Viability (post-vialing)</b> Visual observation LIVE/DEAD <sup>®</sup> BacLight <sup>™</sup> Bacterial Viability	Growth Green fluorescence visible	Growth <sup>2</sup> Green fluorescence visible (Figure 1) <sup>6</sup>

<sup>1</sup>NR-19891 was produced by inoculation of the deposited material into Ellinghausen-McCullough-Johnson-Harrison (EMJH) semisolid agar (0.15%) for 2 days at 30°C in an aerobic atmosphere, and the resulting subculture was viald and frozen. NR-19891 was produced by inoculation of the frozen subculture into EMJH semisolid agar (0.15%) and grown 32 days at 30°C in an aerobic atmosphere. The material from the initial growth was passaged once in EMJH semisolid agar (0.15%) for 38 days under propagation conditions to produce this lot.

<sup>2</sup>Disk of dense growth below the soft agar surface (Dinger's disk) (Czekalowski, J. W., J. W. McLeod and J. Rodican. "The Growth and Respiration of *Leptospira* in Solid or Semi-Solid Media with Special Reference to Dinger's Phenomenon." *Br. J. Exp. Pathol.* 34 (1953): 588-595.) was evident after 12 days at 30°C in EMJH semisolid agar (0.15%).


<sup>3</sup>100% identical to *L. interrogans*, strain HAI0156 (GenBank: AFLP01000082.1)

<sup>4</sup>Purity of this lot was assessed for 7 days at 37°C in an aerobic atmosphere on Tryptic Soy agar with 5% defibrinated sheep blood.

<sup>5</sup>Unexpected growth was observed. Sequencing of the 16S ribosomal RNA gene from the growth confirmed that it was consistent with *L. interrogans* [100% identical to *L. interrogans*, strain HAI0156 (GenBank: AFLP01000082.1)]

<sup>6</sup>Determined after 13 days incubation under cultivation conditions with LIVE/DEAD<sup>®</sup> BacLight<sup>™</sup> Bacterial Viability Kit, 100x magnification (Invitrogen<sup>™</sup> L34856). Cells with a compromised membrane that are dead or dying will stain red, while cells with an intact membrane will stain green.

**Figure 1**



**Date:** 15 DEC 2014

**Signature:** 

**Title:** Technical Manager, BEI Authentication or designee

ATCC<sup>®</sup>, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC<sup>®</sup>'s knowledge.

ATCC<sup>®</sup> is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

