

## **Certificate of Analysis for MRA-1237G**

## Genomic DNA from *Plasmodium falciparum*, Strain IPC 3663

Catalog No. MRA-1237G

**Product Description:** Genomic DNA was extracted from a preparation of *Plasmodium falciparum* (P. falciparum), strain IPC 3663.

Lot<sup>1</sup>: 62401471 Manufacturing Date: 21MAR2014

TEST	SPECIFICATIONS	RESULTS
Sequencing of Merozoite Surface Protein 2 (MSP2) gene (~ 740 base pairs)	Consistent with P. falciparum	Consistent with P. falciparum
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by NanoDrop™ Measurement	Report results	10 μg/mL
Amount	Report results	~ 500 ng per vial
Functional Activity by PCR Amplification MSP2 locus <sup>2</sup>	~ 600-900 base pair amplicon	~ 800 base pair amplicon
OD <sub>260</sub> /OD <sub>280</sub> Ratio	1.6 to 2.1	1.9
Protozoan Inactivation Human erythrocytes exposed to 10% of total yield of MRA-1237G³	No parasitemia observed	No parasitemia observed
Mycoplasma Contamination DNA Detection by PCR	None detected	None detected

<sup>&</sup>lt;sup>1</sup>MRA-1237G was produced from a cell culture of the deposited material. Genomic DNA was extracted using proprietary technology.

Date: 01 JUN 2016 Signature:

**BEI Resources Authentication** 

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.


<sup>&</sup>lt;sup>2</sup>Primer sequences and conditions for PCR are available upon request.

<sup>319</sup> days in complete RPMI culture medium at 37°C in sealed flasks outgassed with blood-gas atmosphere (90% N2, 5% CO2, 5% O2). Complete RPMI culture medium was changed and parasitemia checked every 1 to 4 days.



## **Certificate of Analysis for MRA-1237G**

Figure 1: Agarose Gel Electrophoresis



Lane 1: Invitrogen™ 1 Kb Plus DNA Ladder

Lane 2: ~ 511 ng of MRA-1237G

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

MRA-1237G\_62401471\_01JUN2016