

Certificate of Analysis for NR-50331

Genomic RNA from Zika Virus, H/PAN/2015/CDC-259364

Catalog No. NR-50331

Product Description: Genomic RNA was isolated from a preparation of cell lysate and supernatant from *Cercopithecus aethiops* kidney epithelial (Vero)¹ cells infected with Zika virus (ZIKV), H/PAN/2015/CDC-259364.

Lot²: 64237877 Manufacturing Date: 29APR2016

TEST	SPECIFICATIONS	RESULTS
Sequencing of Species-Specific Region Polyprotein gene NS3–NS4 region (975 nucleotides)	Consistent with ZIKV, H/PAN/2015/CDC-259364 Consistent with NR-50221	100% identity with ZIKV, H/PAN/2015/CDC-259364 ³ (GenBank: KX156776) Consistent with NR-50221
Functional Activity by RT-PCR Amplification ⁴	~ 1100 bp amplicon	~ 1100 bp amplicon (Figure 1)
Total RNA Content by RiboGreen® Measurement (Viral, Cellular and Carrier)	Report results	2.02 ng per 100 µL
Viral RNA Content by Droplet Digital RT-PCR	Report results	3.2 x 10 ³ copies per µL
Virus Inactivation 10% of total yield inoculated on Vero cells and evaluated for cytopathic effect ⁵	No viable virus detected	No viable virus detected

¹Vero 76, clone E6; ATCC[®] CRL-1586™

Date: 27 JUN 2016

Signature:

BEI Resources Authentication

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

Support Provided by NIAID

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898

²Nucleic acid was extracted from a preparation of ZIKV, H/PAN/2015/CDC-259364 (BEI Resources NR-50221, Lot 64158298), using a QIAamp[®] Viral RNA Mini Kit (Qiagen 52906).


³Published in NCBI as ZIKV/Homo sapiens/PAN/CDC-259364_V1-V2/2015

⁴Amplified using iTaq™ DNA Polymerase Kit (Bio-Rad 170-8870) with 2 μL of NR-50331 in a 50 μL reaction

⁵Use of the QlAamp[®] Viral RNA Mini Kit has been demonstrated to consistently inactivate 100% of Zika virus as shown by the absence of cytopathic effect (CPE) after plating the entire extract on virus-susceptible cells.



Certificate of Analysis for NR-50331



Lane 1: Invitrogen™ TrackIt™ 1 Kb Plus DNA Ladder

Lane 2: Negative control Lane 3: NR-50331

E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898