

## **Certificate of Analysis for NR-43737**

## Helicobacter pylori, Strain R056a

Catalog No. NR-43737

**Product Description:** Helicobacter pylori (H. pylori), strain R056a was isolated from gastric biopsy homogenate from an asymptomatic post-menopausal female patient in Alberta, Canada.

Lot<sup>1</sup>: 64136584 Manufacturing Date: 11MAY2016

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology <sup>2</sup>	Report results	Circular, low convex, entire,
Colony morphology	Report results	smooth and gray (Figure 1)
Motility (wet mount)	Report results	Motile
Biochemical characterization	'	
Catalase	Positive	Positive
Oxidase	Positive	Positive
Urease	Positive	Positive
Nitrate reduction	Negative	Negative
H2S (lead acetate paper)	Report results	Positive
Hippurate hydrolysis	Negative	Negative
Growth with 5% CO <sub>2</sub> Growth at 25°C	Growth No growth	Growth No growth
Growth at 37°C	Growth	Growth
Growth at 42°C	Report results	Growth
Brucella albimi + 0.16% agar (growth control)	Growth	Growth
Brucella albimi + 0.16% agar with 1% glycine	No growth	No growth
Brucella albimi + 0.16% agar with 3.5% NaCl	No growth	No growth
Antibiotic Susceptibility Profile  BD BBL™ Sensi-Disc™ susceptibility test discs  Metronidazole (80 µg)³  Nalidixic acid (30 µg)⁴	Report results Report results	54 mm 6 mm
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene (~ 1400 base pairs)	≥ 99% sequence identity to <i>H. pylori</i> , strain R056a (GenBank: AMOY01000003)	100% sequence identity to  H. pylori, strain R056a (GenBank: AMOY01000003.1)
Confirmation of <i>H. pylori</i> by PCR Amplification		
of Extracted DNA		
Positive control (16S ribosomal RNA gene)	Amplicon present	Amplicon present
Negative control ( <i>H. acinonychis</i> )	No amplicon present	No amplicon present
ureA	Amplican present	Amplican present
ssaA	Amplicon present	Amplicon present
Purity (post-freeze)	Consistent with synapted colors:	Consistent with expected estate
Microaerophilic growth <sup>5</sup>	Consistent with expected colony morphology	Consistent with expected colony morphology
Aerobic growth <sup>6</sup>	Consistent with expected colony	Consistent with expected colony
, totobio giowiii	morphology	morphology
Viability (post-freeze) <sup>2</sup>	Growth	Growth

NR-43737 was produced by inoculation of the deposited material into Brucella broth. Broth inoculum was added to Tryptic Soy agar with 5% defibrinated sheep blood. The inoculated broth was grown for 3 days at 37°C in a microaerophilic atmosphere (~ 6-16% O<sub>2</sub> and 2-10% CO<sub>2</sub>) and the inoculated agar was grown for 3 days at 37°C in an aerobic atmosphere with 5% CO<sub>2</sub>. Colonies from the Tryptic Soy agar culture were

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898



## Certificate of Analysis for NR-43737

suspended into the Brucella broth growth, and this biphasic culture was added to Tryptic Soy agar with 5% defibrinated sheep blood kolles, which were grown for 2 days at 37°C in an aerobic atmosphere with 5% CO<sub>2</sub> to produce this lot.

<sup>25</sup> days on Columbia agar with 7% defibrinated horse blood, 5 μg/mL trimethoprim, 5 μg/mL vancomycin, 10 μg/mL cefsulodin and 2.5 μg/mL amphotericin B at 37°C in a microaerophilic atmosphere

³Test performed using metronidazole 80 µg (MET-80, BBL™ catalog no. 231605)

<sup>4</sup>Test performed using nalidixic acid 30 μg (NA-30, BBL™ catalog no. 230874)

Figure 1: Colony Morphology



**Date:** 05 AUG 2016

Signature:

**BEI** Resources Authentication

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898

<sup>&</sup>lt;sup>5</sup>Purity of this lot was assessed for 7 days on Tryptic Soy agar with 5% defibrinated sheep blood at 37°C in a microaerophilic atmosphere.

<sup>&</sup>lt;sup>6</sup>Purity of this lot was assessed for 7 days on Tryptic Soy agar with 5% defibrinated sheep blood at 37°C in an aerobic atmosphere with 5% CO<sub>2</sub>.