

## Certificate of Analysis for NR-49291

## Clostridium difficile, Isolate 20120236

## Catalog No. NR-49291

**Product Description:** Clostridium difficile (C. difficile; also referred to as Peptoclostridium difficile), isolate 20120236 was obtained from the stool of an older female patient with a community-associated (CA) C. difficile infection in midwestern USA in 2011. Isolate 20120236 was deposited as PCR ribotype 027, North American pulsed-field gel electrophoresis type 1 (NAP1), containing tcdA, tcdB and tcdC (with 18 base pair deletion) of the PaLoc operon as well as the C. difficile binary toxin (CDT).

Lot<sup>1</sup>: 63719806 Manufacturing Date: 23SEP2015

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-positive rods	Gram-positive rods
Colony morphology <sup>2</sup>	Report results	Irregular, flat, lobate and gray (Figure 1)
Hemolysis <sup>2</sup>	Report results	Non-hemolytic
Motility (wet mount)	Report results	Motile
Biochemical tests:	·	
Esculin hydrolysis <sup>3</sup>	Positive	Positive
Gelatin hydrolysis <sup>3</sup>	Positive	Positive
VITEK® MS (MALDI-TOF)	C. difficile	C. difficile (99.9%)
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene (~ 1440 base pairs)	≥ 99% sequence identity to  C. difficile type strain  (GenBank: CP011968.1)	99.9% sequence identity to C. difficile type strain (GenBank: CP011968.1)
PCR Assay of Extracted DNA Presence of <i>C. difficile</i> -specific genes <sup>4</sup> Triose phosphate isomerase ( <i>tpi</i> )	~ 230 base pair amplicon	~ 230 base pair amplicon
Presence of toxin genes <sup>4,5</sup>		
cdtB	~ 510 base pair amplicon	~ 510 base pair amplicon
tcdA (wild type)	~ 370 base pair amplicon	~ 370 base pair amplicon
tcdA (partial deletion)	No amplicon	No amplicon
tcdB	~ 160 base pair amplicon	~ 160 base pair amplicon
Purity (post-freeze)		
Aerobic growth <sup>6</sup>	Consistent with expected colony morphology	Consistent with expected colony morphology
Aerobic growth <sup>7</sup>	No growth	No growth
Viability (post-freeze) <sup>2</sup>	Growth	Growth

<sup>&</sup>lt;sup>1</sup>NR-49291 was produced by inoculation of the deposited material into Modified Reinforced Clostridial medium and incubated for 1 day at 37°C in an anaerobic atmosphere (< 0.5% O₂; Remel™ Anaero Pack-Anaero™). The material from the initial growth was passaged once in Modified Reinforced Clostridial medium for 1 day at 37°C in an anaerobic atmosphere to produce this lot.

**BEI Resources** 

www.beiresources.org

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

<sup>&</sup>lt;sup>2</sup>2 days at 37°C in an anaerobic atmosphere on Tryptic Soy agar with 5% defibrinated sheep blood

<sup>&</sup>lt;sup>3</sup>Tests were assessed after 7 days at 37°C in an anaerobic atmosphere. The gelatin tube was placed at 4°C for one hour prior to result determination.

<sup>&</sup>lt;sup>4</sup>Lemee, L., et al. "Multiplex PCR Targeting *tpi* (Triose Phosphate Isomerase), *tcdA* (Toxin A), and *tcdB* (Toxin B) Genes for Toxigenic Culture of *Clostridium difficile*." J. Clin. Microbiol. 42 (2004): 5710-5714. PubMed: 15583303.

<sup>&</sup>lt;sup>5</sup>Antikainen, J., et al. "Detection of Virulence Genes of Clostridium difficile by Multiplex PCR." APMIS 117 (2009): 607-613. PubMed: 19664132.

<sup>&</sup>lt;sup>6</sup>Purity of this lot was assessed for 2 days at 37°C in an anaerobic atmosphere on Tryptic Soy agar with 5% defibrinated sheep blood.

Purity of this lot was assessed for 9 days at 37°C in an aerobic atmosphere with 5% CO<sub>2</sub> on Tryptic Soy agar with 5% defibrinated sheep blood.



## Certificate of Analysis for NR-49291

Figure 1: Colony Morphology



**Date:** 12 OCT 2016

Signature:

**BEI Resources Authentication** 

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources
www.beiresources.org

E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898