

Certificate of Analysis for NR-3047

Genomic DNA from Escherichia coli, Strain 94-G7771

Catalog No. NR-3047

Product Description: Genomic DNA was isolated from a preparation of *Escherichia coli*, strain 94-G7771, serotype O157:NM.

Lot¹: 7559645 Manufacturing Date: 05AUG2006

TEST	SPECIFICATIONS	RESULTS
Sequencing of 16S Ribosomal RNA Gene (~ 450 bp)	Consistent with Escherichia coli	Consistent with Escherichia coli ²
Presence of Plasmid Confirmed by PCR Amplification pO157	Confirmed	Confirmed
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by PicoGreen [®] Measurement	4 to 6 μg in 25 to 100 μL per vial	5.7 μg in 39 μL per vial (148 μg/mL)
Functional Activity by PCR Amplification 16S ribosomal RNA gene Gene-length DNA from pO157 (hlyA virulence marker)	~ 1500 bp amplicon ~ 3200 bp amplicon	~ 1500 bp amplicon ~ 3200 bp amplicon
OD ₂₆₀ /OD ₂₈₀ Ratio	1.7 to 1.9	1.9
Bacterial Inactivation 10% of total yield plated on Trypticase Soy Agar ^{3,4}	No viable bacteria detected	No viable bacteria detected

The bacterial preparation used for extraction of genomic DNA was produced by inoculation of BEI Resources NR-9 (Lot: 3592832) onto Trypticase Soy Agar. After incubation for 20 hours at 37°C, genomic DNA was extracted using proprietary technology.

www.beiresources.org

Date: 06 NOV 2007 **Signature:** Signature on File

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

800-359-7370

NR-3047 7559645 06NOV2007


²Also consistent with *Shigella* species.

³7 days at 37°C in an aerobic atmosphere.

⁴An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-negative bacteria.



Certificate of Analysis for NR-3047



Lane 1: Invitrogen™ 1Kb Plus DNA Ladder™

Lane 2: 200 ng of NR-3047

Page 2 of 2

Fax: 703-365-2898

E-mail: contact@beiresources.org

800-359-7370