SUPPORTING INFECTIOUS DISEASE RESEARCH

Escherichia coli, Strain JJ1886

Catalog No. NR-51620

Product Description:

Escherichia coli (*E. coli*), strain JJ1886 was isolated in 2007 from a woman with fatal urosepsis in the United States. *E. coli*, strain JJ1886 was deposited as sensitive to trimethoprim/sulfamethoxazole, doxycycline, meropenem, nitrofurantoin and polymyxin B and resistant to ciprofloxacin and levofloxacin.

Lot: 70022776¹

Manufacturing Date: 06FEB2019

| TEST | SPECIFICATIONS | RESULTS |
|---|--|--|
| Phenotypic Analysis | | |
| Cellular morphology | Gram-negative rods | Gram-negative rods |
| Colony morphology ² | Report results | Circular, low convex, entire, smooth and cream (Figure 1) |
| Motility (wet mount) | Report results | Motile |
| VITEK® 2 (GN card) | <i>E. coli</i> (≥ 89%) | E. coli (98%) |
| VITEK [®] MS (MALDÍ-TOF) | E. coli | E. coli (99.9%) |
| Antibiotic Susceptibility Profile ³ | | |
| VITEK [®] (AST-GN69 card) | | |
| ESBL ⁴ | Report results | Positive |
| Ampicillin | Report results | Resistant (≥ 32 µg/mL) |
| Amoxicillin/Clavulanic Acid | Report results | Resistant (≥ 32 µg/mL) |
| Ampicillin/Sulbactam | Report results | Resistant (≥ 32 µg/mL) |
| Piperacillin/Tazobactam | Report results | Sensitive (8-16 µg/mL) |
| Cefazolin | Report results | Resistant (≥ 64 µg/mL) |
| Ceftazidime | Report results | Sensitive (4 µg/mL)⁵ |
| Ceftriaxone | Report results | Resistant (≥ 64 µg/mL) |
| Cefepime | Report results | Sensitive (2 µg/mL) ⁶ |
| Ertapenem | Report results | Sensitive (≤ 0.5 µg/mL) |
| Imipenem | Report results | Sensitive (≤ 0.25 µg/mL) |
| Gentamicin | Report results | Sensitive (≤ 1µg/mL) |
| Tobramycin | Report results | Inconclusive ⁷ |
| Ciprofloxacin | Resistant | Resistant (≥ 4 µg/mL) |
| Levofloxacin | Resistant | Resistant (≥ 8 µg/mL) |
| Nitrofurantoin | Sensitive | Sensitive (≤ 16 µg/mL) |
| Trimethoprim/Sulfamethoxazole | Sensitive | Sensitive (≤ 20 µg/mL) |
| Etest [®] antibiotic test strips ⁸ | | |
| Doxycycline | Sensitive | Sensitive (2 µg/mL) |
| Meropenem | Sensitive | Sensitive (0.023 µg/mL) |
| Polymyxin B | Sensitive | 0.38 µg/mL ⁹ |
| Genotypic Analysis | | |
| Sequencing of 16S ribosomal RNA gene (~ 1480 base pairs) | ≥ 99% sequence identity with <i>E. coli,</i> strain JJ1886 (GenBank: CP006784.1) | 99.8% sequence identity with <i>E. coli,</i> strain JJ1886 (GenBank: CP006784.1) ¹⁰ |
| Purity (post-freeze) ^{11,12} | Growth consistent with expected colony morphology | Growth consistent with expected colony morphology |
| Viability (post-freeze) ² | Growth | Growth |

¹The deposited material was inoculated into Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere and preserved in 10% glycerol. NR-51620 was produced by the inoculation of the preserved material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot.

²1 day at 37°C in an aerobic atmosphere on Tryptic Soy agar

BEI Resources www.beiresources.org E-mail: <u>contact@beiresources.org</u> Tel: 800-359-7370 Fax: 703-365-2898 **b**|**e**|**i** resources

Certificate of Analysis for NR-51620

SUPPORTING INFECTIOUS DISEASE RESEARCH

³Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

⁴The VITEK[®]2 ESBL Test is a confirmatory test for Extended-Spectrum Beta-Lactamases (ESBLs) inhibited by clavulanic acid and utilizes cefepime, cefotaxime and ceftazidime, with and without clavulanic acid, to determine a positive or negative result.

⁵Because *E. coli*, strain JJ1886 is a confirmed ESBL-producer, CSLI recommendations are to modify the interpretation of ceftazidime based on the susceptibilities of other antibiotics in the same class, suggesting that this strain is resistant to ceftazidime.

⁶Because *E. coli*, strain JJ1886 is a confirmed ESBL-producer, CSLI recommendations are to modify the interpretation of cefepime based on the susceptibilities of other antibiotics in the same class, suggesting that this strain is resistant to cefepime.

⁷Antibiotic susceptibility testing performed in duplicate determined that for strain JJ1886, the tobramycin MICs are 8 µg/mL and 16 µg/mL, which are interpreted as intermediate and resistant, respectively.

⁸1 day at 37°C in an aerobic atmosphere on Mueller Hinton agar

⁹ CLŚI does not have published polymyxin B MIC breakpoints for *E. coli*. Isolates are defined only as wild type or non-wild type. For more information, please refer to Chew, K. L., et al. "Colistin and Polymyxin B Susceptibility Testing for Carbapenem-Resistant and *mcr*-Positive Enterobacteriaceae: Comparison of Sensititre, MicroScan, Vitek 2, and Etest with Broth Microdilution." <u>J. Clin. Microbiol.</u> 55 (2017): 2609-2616. PubMed: 28592552.
¹⁰Also consistent with *Shigella* and other *Escherichia* species

Figure 1: Colony Morphology

¹¹Purity of this lot was assessed for 8 days at 37°C in an aerobic atmosphere with 5% CO₂ on Tryptic Soy agar.

¹²Purity of this lot was assessed for 8 days at 37°C in an aerobic atmosphere on Tryptic Soy agar.

/Heather Couch/ Heather Couch

20 SEP 2019

Program Manager or designee, ATCC Federal Solutions

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.

