

Certificate of Analysis for NR-51590

Pseudomonas aeruginosa, Strain MRSN 18754

Catalog No. NR-51590

This reagent is the tangible property of the U.S. Government.

Product Description:

Pseudomonas aeruginosa (P. aeruginosa), strain MRSN 18754 was isolated in 2013 from human tissue as part of a surveillance program in the United States. *P. aeruginosa*, strain MRSN 18754 was deposited as sensitive to amikacin, aztreonam, cefepime, ceftazidime, ciprofloxacin, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam and tobramycin.

Lot: 70025086¹ Manufacturing Date: 01AUG2019

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology ²	Report results	Irregular, low convex, undulate,
· •	·	smooth and green (Figure 1)
Motility (wet mount)	Report results	Motile
VITEK [®] 2 (GN card)	P. aeruginosa (≥ 89%)	P. aeruginosa (99%)
Antibiotic Susceptibility Profile ³		
VITEK® (AST-GN81 Card)		
Ampicillin	Report results	Resistant (≥ 32 µg/mL)
Amoxicillin/clavulanic acid	Report results	Resistant (≥ 32 µg/mL)
Piperacillin/tazobactam	Sensitive	Sensitive (8 µg/mL)
Cefazolin	Report results	Resistant (≥ 64 µg/mL)
Cefoxitin	Report results	Resistant (≥ 64 µg/mL)
Ceftazidime	Sensitive	Sensitive (4 µg/mL)
Ceftriaxone	Report results	Resistant (≥ 64 µg/mL)
Cefepime	Sensitive	Sensitive (4 µg/mL)
Meropenem	Sensitive	Sensitive (1 µg/mL)
Amikacin	Sensitive	Sensitive (≤ 2 µg/mL)
Gentamicin	Sensitive	Sensitive (≤ 1 μg/mL)
Tobramycin	Sensitive	Sensitive (≤ 1 μg/mL)
Ciprofloxacin	Sensitive	Sensitive (≤ 0.25 μg/mL)
Levofloxacin	Sensitive	Sensitive (≤ 1 µg/mL)
Tetracycline	Report results	Resistant (≥ 16 µg/mL)
Nitrofurantoin	Report results	Resistant (≥ 256 µg/mL)
Trimethoprim/sulfamethoxazole	Report results	160 μg/mL ⁴
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene (~ 1410 base pairs)	≥ 99% sequence identity to P. aeruginosa, strain MRSN 18754	99.9% sequence identity to P. aeruginosa, strain MRSN 18754
(~ 1410 base pairs)	(GenBank: RXVH01000074.1)	(GenBank: RXVH01000074.1)
Purity (post-freeze) ⁵	Growth consistent with expected	Growth consistent with expected
	colony morphology	colony morphology
Viability (post-freeze) ²	Growth	Growth

¹NR-51590 was produced by inoculation of the depositor material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Tel: 800-359-7370 Fax: 703-365-2898

²1 day at 37°C in an aerobic atmosphere on Tryptic Soy agar

³Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

⁴Trimethoprim/sulfamethoxazole MIC interpretive standards are not available for *P. aeruginosa*, however most clinical isolates are resistant to trimethoprim/sulfamethoxazole. For more information, please refer to Köhler, T., et al. "Multidrug Efflux in Intrinsic Resistance to Trimethoprim and Sulfamethoxazole in *Pseudomonas aeruginosa.*" <u>Antimicrob. Agents Chemother.</u> 40 (1996): 2288-2290. PubMed: 9036831.

⁵Purity of this lot was assessed for 7 days at 37°C in an aerobic atmosphere with and without 5% CO₂ on Tryptic Soy agar.



Certificate of Analysis for NR-51590

Figure 1: Colony Morphology



/Heather Couch/ Heather Couch

04 NOV 2019

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

 $\label{eq:attack} \mbox{ATCC}^{\mbox{\tiny{\$}}} \mbox{ is a trademark of the American Type Culture Collection.}$ You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898