SUPPORTING INFECTIOUS DISEASE RESEARCH

## Borrelia afzelii, Strain Pko

#### Catalog No. NR-51676

#### **Product Description:**

Borrelia afzelii (B. afzelii), strain Pko was isolated in 1984 from the skin of a human with erythema migrans (Lyme borreliosis) in Germany. NR-51676 lot 70027343 was produced by inoculation of the deposited material into Revised Barbour-Stoenner-Kelly broth and grown for 3 days at 33°C in a microaerophilic atmosphere (6-16% O<sub>2</sub> and 2-10% CO<sub>2</sub>; BD GasPak<sup>™</sup> EZ Campy). Broth inoculum was added to Revised Barbour-Stoenner-Kelly broth and grown for 7 days at 33°C in a microaerophilic atmosphere to produce this lot.

#### Lot: 70027343


## Manufacturing Date: 28OCT2019

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Spirochete	Spirochete
6 days at 33°C in a microaerophilic atmosphere in		
Revised Barbour-Stoenner-Kelly broth		
Motility (wet mount)	Report results	Motile
Genotypic Analysis		
Sequencing of 16S ribosomal RNA (rRNA) gene	≥ 99% sequence identity to	100% sequence identity to
(~ 1350 base pairs)	<i>B. afzelii</i> , strain Pko	<i>B. afzelii</i> , strain Pko
	(GenBank: CP000395.1)	(GenBank: CP000395.1) <sup>1</sup>
Purity (post-freeze)		
7 days at 33°C in a microaerophilic atmosphere in	Growth consistent with colony	No growth
Tryptic Soy agar with 5% defibrinated sheep blood	morphology or no growth	
7 days at $37^{\circ}$ C in an aerobic atmosphere with $5\%$ CO <sub>2</sub>	No growth	No growth
in Tryptic Soy agar with 5% defibrinated sheep blood		
Viability (post-freeze)		
Visual observation	Growth	Growth
6 days at 33°C in a microaerophilic atmosphere in		
Revised Barbour-Stoenner-Kelly broth		
LIVE/DEAD <sup>®</sup> <i>Bac</i> Light™ Bacterial Viability	Green fluorescence visible	Green fluorescence visible (Figure 1) <sup>2</sup>

<sup>1</sup>Also consistent with other *Borrelia* species.

<sup>2</sup>Determined after 6 days at 33°C in a microaerophilic atmosphere in Revised Barbour-Stoenner-Kelly broth with LIVE/DEAD<sup>®</sup> BacLight<sup>™</sup> Bacterial Viability Kit, 1000× magnification (Invitrogen<sup>™</sup> L34856). Cells with a compromised membrane that are dead or dying will stain red, while cells with an intact membrane will stain green.





E-mail: <u>contact@beiresources.org</u> Tel: 800-359-7370 Fax: 703-365-2898 bei resources

# **Certificate of Analysis for NR-51676**

SUPPORTING INFECTIOUS DISEASE RESEARCH

# /Heather Couch/

Heather Couch

Program Manager or designee, ATCC Federal Solutions

ATCC<sup>®</sup>, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC<sup>®</sup>'s knowledge.

ATCC<sup>®</sup> is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.



22 APR 2020