

Vector pMCSG53 Containing the SARS-Related Coronavirus 2, Wuhan-Hu-1 Spike Glycoprotein N-Terminal Domain

Catalog No. NR-52428

This reagent is the tangible property of the U.S. Government.

Product Description:

The vector for the N-terminal domain (NTD) of the spike (S) glycoprotein gene from severe acute respiratory syndrome-related coronavirus 2 (SARS-CoV-2), Wuhan-Hu-1 (GenBank: [MN908947](#)) was designed by subcloning the codon-optimized S protein NTD (amino acids 26 to 305) into the pMCSG53 *Escherichia coli* (*E. coli*) expression vector. pMCSG53 is a ligation-independent cloning (LIC) vector containing an N-terminal hexa-histidine tag and tobacco etch virus (TEV) protease recognition site prior to the NTD. In addition, the vector includes tRNA genes covering rare codons for arginine (AGG/AGA) and isoleucine (AUA) to improve expression in the host, *E. coli*. NR-52428 contains the beta-lactamase gene, *bla*, to provide transformant selection through ampicillin resistance in *E. coli*. The deposited plasmid was transformed into One Shot™ TOP10 *E. coli* (Invitrogen™ C404003), grown in Luria-Bertani broth with ampicillin (50 µg per mL) for 1 day at 37°C in an aerobic atmosphere, extracted using a Plasmid Plus Maxi Kit (QIAGEN® 12963) and vialled in TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0).

Lot: 70035126

Manufacturing Date: 20APR2020

TEST	SPECIFICATIONS	RESULTS
Next-Generation DNA Sequencing	~ 5650 base pairs	5651 base pairs ¹
Genotypic Analysis Sequencing of Spike NTD insert (~ 850 base pairs) Sequencing of pMCSG53 vector	100% sequence identity to depositor's sequence His ₆ tag sequence confirmed TEV protease site sequence confirmed	100% sequence identity to depositor's sequence ² His ₆ tag sequence confirmed TEV protease site sequence confirmed
Antibiotic Resistance Ampicillin (encoded by beta-lactamase gene <i>bla</i>) ³	<i>bla</i> sequence present	<i>bla</i> sequence present
Concentration by PicoGreen® Measurement	≥ 2 µg/mL	0.04 µg in 20 µL per vial (2.2 µg/mL)
Amount per Vial	Report results	0.04 µg per vial
OD₂₆₀/OD₂₈₀ Ratio	1.7 to 2.1	1.9
Effective Bacterial Transformation Invitrogen™ One Shot™ TOP10 <i>E. coli</i>	≥ 50 colonies per ng	88 colonies per ng

¹The sequence was assembled pre-vial using the depositor's predicted sequence as the reference sequence. The complete plasmid sequence and map are provided on the BEI Resources webpage.

²The NR-52428 insert was codon optimized but is consistent with the SARS-CoV-2, Wuhan-Hu-1 S protein (GenPept: QHD43416).

³The antibiotic ampicillin degrades quickly during growth. Bacterial stationary phase should be minimized during plasmid expansion to avoid plasmid loss and increased antibiotic concentrations may be necessary.

/Heather Couch/

Heather Couch

01 MAY 2020

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by ATCC® and the contributor to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

