

## **Certificate of Analysis for NR-53911**

Whole Cell Antigen, SARS-Related Coronavirus 2, Isolate USA-WA1/2020 with Uninfected Control Vero E6 Cells with High Expression of Angiotensin-Converting Enzyme 2, Gamma-Irradiated

Catalog No. NR-53911

## **Product Description:**

NR-53911 consists of crude preparation of cell lysate and supernatant from uninfected *Cercopithecus aethiops* kidney epithelial cells with high expression of angiotensin-converting enzyme 2 (Vero E6-heACE2) and cells infected with severe acute respiratory syndrome-related coronavirus 2 (SARS-CoV-2), isolate USA-WA1/2020. These lysates were gamma-irradiated (5 × 10<sup>6</sup> RADs) on dry ice, followed by sonication.

The starting material for NR-53910 was produced by infecting VE6-heACE2 in Eagle's Minimum Essential Medium (ATCC® 30-2003™) supplemented with 2% fetal bovine serum (ATCC® 30-2020™) for 5 days at 37°C with 5% CO₂. Titer of the pre-inactivation material was 1.2 × 10<sup>6</sup> TCID<sub>50</sub> per mL by TCID<sub>50</sub> assay in Vero E6-heACE2 cells, determined by cytopathic effect in 5 days at 37°C with 5% CO₂, and sterility of the preparation was confirmed.<sup>1,2</sup>

Table 1: SARS-CoV-2 Whole Cell Antigen Kit Components

| Tuble 1: OARS GOV 2 Whole Gen Antigen Ric Components |   |            |                     |  |
|--|---|------------|---------------------|--|
| COMPONENT NUMBER                                     | DESCRIPTION   | LOT NUMBER | DATE OF MANUFACTURE |  |
| NR-53909   | Gamma-irradiated uninfected Vero E6-heACE2 cell lysate control                        | 70041014   | 07JAN2021           |  |
| NR-53910   | Gamma-irradiated SARS-CoV-2, isolate USA-WA1/2020 infected Vero E6-heACE2 cell lysate | 70041015   | 07JAN2021           |  |

## Table 2: Gamma-Irradiated SARS-CoV-2, isolate USA-WA1/2020 infected Vero E6-heACE2 Cell Lysate (NR-53910)

| 1444-000 (0)   |                          |                          |  |
|--|--------------------------|--------------------------|--|
| TEST   | SPECIFICATIONS           | RESULTS                  |  |
| Virus Inactivation   |                          |                          |  |
| 10% of total bulk gamma-irradiated preparation inoculated on Vero E6 cells and evaluated for cytopathic effect and expression of viral antigen by indirect immunofluorescence assay <sup>2,3</sup> | No viable virus detected | No viable virus detected |  |

¹The Tissue Culture Infectious Dose 50% (TCID₅₀) is the 50% infectious endpoint in cell culture. The TCID₅₀ is the dilution of virus that under the conditions of the assay can be expected to infect 50% of the culture vessels inoculated, just as a Lethal Dose 50% (LD₅₀) is expected to kill half of the animals exposed. A reciprocal of the dilution required to yield the TCID₅₀ provides a measure of the titer (or infectivity) of a virus preparation.

/Heather Couch/

Heather Couch 18 JUN 2021

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by ATCC® or the contractor to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: <a href="mailto:contact@beiresources.org">contact@beiresources.org</a>
Tel: 800-359-7370

Fax: 703-365-2898

<sup>&</sup>lt;sup>2</sup>Performed at University of Texas Medical Branch, Galveston, Texas, USA

<sup>&</sup>lt;sup>3</sup>The gamma-irradiated virus preparation was plated on Vero E6-heACE2 cells and incubated for 14 days at 37°C and 5% CO<sub>2</sub>; cell lysate and supernatant from these cultures were blind passaged on fresh monolayers of Vero E6-heACE2 cells and again incubated for 14 days at 37°C and 5% CO<sub>2</sub>.