

Certificate of Analysis for NR-55583

Klebsiella pneumoniae, Strain MRSN 614201

Catalog No. NR-55583

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (K. pneumoniae), strain MRSN 614201 was isolated in 2018 from an environmental sample in Africa as part of a global surveillance program. NR-55583 was deposited as a multidrug-resistant (MDR) strain, sensitive to amikacin, ceftazidime/avibactam, ceftolozane/tazobactam, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline, tobramycin and trimethoprim/sulfamethoxazole and resistant to ampicillin/sulbactam, aztreonam, ceftazidime, ceftriaxone and cefepime. NR-55583 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051603 Manufacturing Date: 06APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology	Report results	Circular, convex, entire, smooth and cream (Figure 1)
Motility (wet mount)	Report results	Non-motile
VITEK [®] 2 (GN card)	K. pneumoniae (≥ 89%)	K. pneumoniae (99%)
Antibiotic Susceptibility Profile ^{1,2}		
Amikacin	Sensitive	Sensitive (≤ 2 μg/mL)
Ampicillin/sulbactam	Resistant	Sensitive (4 μg/mL) ³
Aztreonam	Resistant	Sensitive (2 μg/mL) ⁴
Cefepime	Resistant	Sensitive (≤ 1 µg/mL) ⁵
Ceftazidime	Resistant	Sensitive (4 μg/mL) ⁶
Ceftazidime/avibactam	Sensitive	Sensitive (0.094 µg/mL)
Ceftolozane/tazobactam	Sensitive	Sensitive (0.125 µg/mL)
Ceftriaxone	Resistant	Resistant (16 µg/mL)
Ciprofloxacin	Sensitive	Sensitive (0.016 µg/mL)
Ertapenem	Sensitive	Sensitive (≤ 0.5 μg/mL)
Gentamicin	Sensitive	Sensitive (≤ 1 μg/mL)
Imipenem	Sensitive	Sensitive (0.25 µg/mL)
Levofloxacin	Sensitive	Sensitive (≤ 0.12 μg/mL)
Meropenem	Sensitive	Sensitive (≤ 0.25 μg/mL)
Piperacillin/tazobactam	Sensitive	Sensitive (≤ 4 μg/mL)
Tetracycline	Sensitive	Sensitive (≤ 1 μg/mL)
Tigecycline	Sensitive	Sensitive (≤ 0.5 μg/mL) ⁷
Tobramycin	Sensitive	Sensitive (≤ 1 μg/mL)
Trimethoprim/sulfamethoxazole	Sensitive	Sensitive (≤ 20 μg/mL)
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene	≥ 99% sequence identity to	99.9% sequence identity to
(1470 base pairs)	K. pneumoniae, strain MRSN 614201 (GenBank: JAGYCC010000052.1)	K. pneumoniae, strain MRSN 614201 (GenBank: JAGYCC010000052.1)

BEI Resources
www.beiresources.org

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898



Certificate of Analysis for NR-55583

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar		Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

Figure 1: Colony Morphology



/Sonia Bjorum Brower/ Sonia Bjorum Brower

23 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org Tel: 800-359-7370

Fax: 703-365-2898

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK®2 GN74 and ETEST®.

³K. pneumoniae, strain MRSN 614201 was deposited as resistant to ampicillin/sulbactam but showed a MIC of 4 μg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

4K. pneumoniae, strain MRSN 614201 was deposited as resistant to aztreonam but showed a MIC of 2 µg per mL (interpreted as sensitive) for this

antibiotic during QC testing. Testing was performed in duplicate.

⁵K. pneumoniae, strain MRSN 614201 was deposited as resistant to cefepime but showed a MIC of ≤ 1 μg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁶K. pneumoniae, strain MRSN 614201 was deposited as resistant to ceftazidime but showed a MIC of 4 μg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁷MIC Interpretation Guideline: EUCAST Version 8.0 (2018)