

Certificate of Analysis for NR-55533

Klebsiella pneumoniae, Strain MRSN 25112

Catalog No. NR-55533

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (K. pneumoniae), strain MRSN 25112 was isolated in 2014 from a human urine sample in North America as part of a global surveillance program. NR-55533 was deposited as a multidrug-resistant strain, sensitive to amikacin, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftriaxone, ceftolozane/tazobactam, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, tigecycline, tobramycin and tetracycline, intermediately resistant to piperacillin/tazobactam and resistant to ampicillin/sulbactam and trimethoprim/sulfamethoxazole. NR-55533 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70050362 Manufacturing Date: 25FEB2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS RESULTS				
Phenotypic Analysis					
Cellular morphology	Gram-negative rods	Gram-negative rods			
Colony morphology	Report results	Circular, convex, entire, smooth mucoid and cream			
Motility (wet mount)	Report results	Non-motile			
VITEK [®] 2 (GN card)	K. pneumoniae (≥ 89%)	K. pneumoniae (99%)			
Antibiotic Susceptibility Profile ^{1,2}					
Amikacin	Sensitive	Sensitive (≤ 2 µg/mL)			
Ampicillin/sulbactam	Resistant	Resistant (≥ 32 µg/mL)			
Aztreonam	Sensitive	Sensitive (≤ 1 μg/mL)			
Cefepime	Sensitive	Sensitive (≤ 1 µg/mL)			
Ceftazidime	Sensitive	Sensitive (≤ 1 µg/mL)			
Ceftazidime/avibactam	Sensitive	Sensitive (0.094 µg/mL)			
Ceftolozane/tazobactam	Sensitive	Sensitive (0.064 µg/mL)			
Ceftriaxone	Sensitive	Sensitive (≤ 1 µg/mL)			
Ciprofloxacin	Sensitive	Sensitive (0.016 µg/mL)			
Ertapenem	Sensitive	Sensitive (≤ 0.5 µg/mL)			
Gentamicin	Sensitive	Sensitive (≤ 1 µg/mL)			
Imipenem	Sensitive	Sensitive (0.19 µg/mL)			
Levofloxacin	Sensitive	Sensitive (≤ 0.12 μg/mL)			
Meropenem	Sensitive	Sensitive (≤ 0.25 μg/mL)			
Piperacillin/tazobactam	Intermediate	Sensitive (3 µg/mL) ³			
Tetracycline	Sensitive	Sensitive (≤ 1 µg/mL)			
Tigecycline	Sensitive	Sensitive (1 µg/mL) ⁴			
Tobramycin	Sensitive	Sensitive (≤ 1 μg/mL)			
Trimethoprim/sulfamethoxazole	Resistant	Resistant (≥ 320 µg/mL)			
Genotypic Analysis					
Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to K. pneumoniae, strain MRSN 25112 (GenBank: JAGYDZ010000105.1)				

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898



Certificate of Analysis for NR-55533

TEST	SPECIFICATIONS			RESULTS				
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	_	consistent morphology	with	expected	_	consistent morphology	with	expected
Viability	Growth				Growth			

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

/Heather Couch/

Heather Couch 06 MAY 2022

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Tel: 800-359-7370 Fax: 703-365-2898

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK®2 GN74 and ETEST®.

³K. pneumoniae, strain MRSN 15219 was deposited as intermediately resistant to piperacillin/tazobactam, but showed a MIC of 3 μg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed twice in duplicate.

⁴MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁵Also consistent with other *Klebsiella* species