

Certificate of Analysis for NR-55593

Klebsiella pneumoniae, Strain MRSN 702325

Catalog No. NR-55593

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (K. pneumoniae), strain MRSN 702325 was isolated in 2019 from a human wound sample in Africa as part of a global surveillance program. NR-55593 was deposited as a multidrug-resistant strain, sensitive to amikacin, aztreonam, cefepime, ceftazidime, ceftriaxone, ciprofloxacin, ceftazidime/avibactam, ceftolozane/tazobactam, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, tobramycin, tigecycline and trimethoprim/sulfamethoxazole, intermediately resistant to piperacillin/tazobactam and resistant to ampicillin/sulbactam and tetracycline. NR-55593 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70051624 Manufacturing Date: 14APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology	Report results	Circular, low convex, entire, smooth and cream (Figure 1)
Motility (wet mount)	Report results	Non-motile
VITEK® 2 (GN card)	K. pneumoniae (≥ 89%)	K. pneumoniae (99%)
Antibiotic Susceptibility Profile ^{1,2}		
Amikacin	Sensitive	Sensitive (≤ 2 μg/mL)
Ampicillin/sulbactam	Resistant	Sensitive (6 µg/mL) ³
Aztreonam	Sensitive	Sensitive (≤ 1 μg/mL)
Cefepime	Sensitive	Sensitive (≤ 1 μg/mL)
Ceftazidime	Sensitive	Sensitive (≤ 1 μg/mL)
Ceftazidime/avibactam	Sensitive	Sensitive (0.19 μg/mL)
Ceftolozane/tazobactam	Sensitive	Sensitive (0.125 µg/mL)
Ceftriaxone	Sensitive	Sensitive (≤ 1 µg/mL)
Ciprofloxacin	Sensitive	Sensitive (0.032 µg/mL)
Ertapenem	Sensitive	Sensitive (≤ 0.5 μg/mL)
Gentamicin	Sensitive	Sensitive (≤ 1 µg/mL)
Imipenem	Sensitive	Sensitive (0.19 to 0.25 µg/mL)
Levofloxacin	Sensitive	Sensitive (≤ 0.12 μg/mL)
Meropenem	Sensitive	Sensitive (≤ 0.25 μg/mL)
Piperacillin/tazobactam	Intermediate	Sensitive (1.5 µg/mL) ⁴
Tetracycline	Resistant	Resistant (≥ 16 µg/mL)
Tigecycline	Sensitive	Sensitive (1 µg/mL) ⁵
Tobramycin	Sensitive	Sensitive (≤ 1 µg/mL)
Trimethoprim/sulfamethoxazole	Sensitive	Sensitive (≤ 20 μg/mL)
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene (~ 1470 base pairs)	≥ 99% sequence identity to K. pneumoniae, strain MRSN 702325 (GenBank: JAGYBS010000047.1)	99.5% sequence identity to K. pneumoniae, strain MRSN 702325 (GenBank: JAGYBS010000047.1) ⁶

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898



Certificate of Analysis for NR-55593

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

Figure 1: Colony Morphology



/Sonia Bjorum Brower/ Sonia Bjorum Brower

30 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK®2 GN74 and ETEST®.

³K. pneumoniae, strain MRSN 702325 was deposited as resistant to ampicillin/sulbactam, but showed a MIC of 6 µg per mL (interpreted as sensitive)

for this antibiotic during QC testing. Testing was performed in duplicate.

4κ. pneumoniae, strain MRSN 702325 was deposited as intermediately resistant to piperacillin/tazobactam, but showed a MIC of 1.5 μg per mL (interpreted as sensitive) for this antibiotic during QC testing. Testing was performed in duplicate.

⁵MIC Interpretation Guideline: EUCAST Version 8.0 (2018)

⁶Also consistent with other *Klebsiella* species