

Certificate of Analysis for NR-55601

Klebsiella pneumoniae, Strain MRSN 750999

Catalog No. NR-55601

This reagent is the tangible property of the U.S. Government.

Product Description:

Klebsiella pneumoniae (K. pneumoniae), strain MRSN 750999 was isolated in 2020 from a human urine sample in North America as part of a global surveillance program. NR-55601 was deposited as a susceptible strain, sensitive to amikacin, ampicillin/sulbactam, aztreonam, cefepime, ceftazidime, ceftazidime/avibactam, ceftolozane/tazobactam, ceftriaxone, ciprofloxacin, ertapenem, gentamicin, imipenem, levofloxacin, meropenem, piperacillin/tazobactam, tetracycline, tigecycline, tobramycin, and trimethoprim/sulfamethoxazole. NR-55601 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed on Tryptic Soy agar under propagation conditions unless otherwise noted.

Lot: 70051640 Manufacturing Date: 06APR2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: contact@BEIResources.org. We try to respond to feedback within 24 hours.

| TEST | SPECIFICATIONS | RESULTS |
|--|---|---|
| Phenotypic Analysis | | |
| Cellular morphology | Gram-negative rods | Gram-negative rods |
| Colony morphology | Report results | Circular, convex, entire, smooth and cream |
| Motility (wet mount) | Report results | Non-motile |
| VITEK® 2 (GN card) | K. pneumoniae (≥ 89%) | K. pneumoniae (98%) |
| Antibiotic Susceptibility Profile ^{1,2} | | |
| Amikacin | Sensitive | Sensitive (≤ 2 µg/mL) |
| Ampicillin/sulbactam | Sensitive | Sensitive (4 µg/mL) |
| Aztreonam | Sensitive | Sensitive (≤ 1 µg/mL) |
| Cefepime | Sensitive | Sensitive (≤ 1 μg/mL) |
| Ceftazidime | Sensitive | Sensitive (≤ 1 µg/mL) |
| Ceftazidime/avibactam | Sensitive | Sensitive (0.125 µg/mL) |
| Ceftolozane/tazobactam | Sensitive | Sensitive (0.125 µg/mL) |
| Ceftriaxone | Sensitive | Sensitive (≤ 1 µg/mL) |
| Ciprofloxacin | Sensitive | Sensitive (0.047 µg/mL) |
| Ertapenem | Sensitive | Sensitive (≤ 0.5 μg/mL) |
| Gentamicin | Sensitive | Sensitive (≤ 1 μg/mL) |
| Imipenem | Sensitive | Sensitive (0.25 to 0.38 µg/mL) |
| Levofloxacin | Sensitive | Sensitive (≤ 0.12 μg/mL) |
| Meropenem | Sensitive | Sensitive (≤ 0.25 μg/mL) |
| Piperacillin/tazobactam | Sensitive | Sensitive (≤ 4 μg/mL) |
| Tetracycline | Sensitive | Sensitive (≤ 1 μg/mL) |
| Tigecycline | Sensitive | Sensitive (≤ 0.5 μg/mL) ³ |
| Tobramycin | Sensitive | Sensitive (≤ 1 μg/mL) |
| Trimethoprim/sulfamethoxazole | Sensitive | Sensitive (≤ 20 μg/mL) |
| Genotypic Analysis | | |
| Sequencing of 16S ribosomal RNA gene | ≥ 99% sequence identity to | 99.6% sequence identity to |
| (~ 1470 base pairs) | K. pneumoniae, strain MRSN 750999 (GenBank: JAGYBK010000105.1) | K. pneumoniae, strain MRSN 750999 (GenBank: JAGYBK010000105.1) |

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898




Certificate of Analysis for NR-55601

| TEST | SPECIFICATIONS | RESULTS |
|---|---------------------------------------|---|
| Purity 7 days at 37°C in an aerobic atmosphere with 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood | · · · · · · · · · · · · · · · · · · · | Growth consistent with expected colony morphology |
| Viability | Growth | Growth |

¹Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

Figure 1: Colony Morphology



/Sonia Bjorum Brower/ Sonia Bjorum Brower

20 JUN 2022

Lead Technical Writer or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

Support Provided by NIAID

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898

²Antibiotic susceptibility was tested using a combination of bioMérieux VITEK®2 GN74 and ETEST®.

³MIC Interpretation Guideline: EUCAST Version 8.0 (2018)