

## **Certificate of Analysis for NR-55349**

## Genomic RNA from SARS-Related Coronavirus 2, Isolate hCoV-19/USA/MD-HP01542/2021 (Lineage B.1.351)

Catalog No. NR-55349

## **Product Description:**


Genomic RNA was isolated from a preparation of cell lysate and supernatant from *Homo sapiens* lung adenocarcinoma cells (Calu-3; ATCC® HTB-55™) infected with SARS-CoV-2, isolate hCoV-19/USA/MD-HP01542/2021 (BEI Resources NRC-55282 lot 70043444) using QIAamp® Viral RNA Mini Kit (Qiagen® 52906). The viral genomic RNA is in a background of cellular nucleic acid and carrier RNA.

Lot: 70045171 Manufacturing Date: 01JUN2021

TEST	SPECIFICATIONS	RESULTS
Genotypic Analysis Sequencing of species-specific region (~ 840 nucleotides)	≥ 98% identity with SARS-CoV-2, hCoV-19/USA/MD-HP01542/2021 (GISAID: EPI_ISL_890360)	100% identity with SARS-CoV-2, hCoV-19/USA/MD-HP01542/2021 (GISAID: EPI_ISL_890360)
Functional Activity by RT-PCR Amplification <sup>1</sup> ORF1ab	~ 950 base pair amplicon	~ 1000 base pair amplicon (Figure 1)
Estimated Concentration (post-dilution) by Qubit® Measurement (Viral, Cellular and Carrier) <sup>2</sup>	Report results	0.5 ng per 100 μL (0.0046 μg per mL)
Estimated Amount per Vial <sup>3</sup>	Report results	0.5 ng
Genome Copy Number Using BioRad QX200 Droplet Digital PCR (ddPCR™) System (Post vial; 18 replicates)	Report results	7.8 × 10 <sup>5</sup> genome copies per mL
Virus Inactivation 10% of total yield inoculated on Calu-3 cells and evaluated for cytopathic effect and by RT-PCR after serial passage	No viable virus detected	No viable virus detected

<sup>&</sup>lt;sup>1</sup>Amplified using iTaq™ Universal SYBR Green One-step Kit (Bio-Rad<sup>®</sup> 172-5151) with 5 μL of NR-55349 in a 50 μL reaction

Figure 1: Functional Activity of NR-55349 by RT-PCR Amplification



Lane 1: Invitrogen™ TrackIt™ 1 Kb Plus DNA Ladder Lane 2: PCR product from 1 µL of NR-55349

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

<sup>&</sup>lt;sup>2</sup>Measurement is determined pre-vial prior to dilution due to the limit of detection of the quantification method

<sup>&</sup>lt;sup>3</sup>Use of the QIAamp<sup>®</sup> Viral RNA Mini Kit has been demonstrated to consistently inactivate 100% of SARS-CoV-2 as shown by the absence of cytopathic effect (CPE) and viral RNA expression by RT-PCR after plating the entire extract on virus-susceptible cells for two passages.



## Certificate of Analysis for NR-55349

/Sonia Bjorum Brower/ Sonia Bjorum Brower

02 AUG 2022

Technical Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898