

Certificate of Analysis for NR-55216

Staphylococcus aureus, Strain AJUL2

Catalog No. NR-55216

Product Description:

Staphylococcus aureus (S. aureus), strain AJUL2 is deposited as a gentamycin, kanamycin and tobramycin-resistant strain derived from S. aureus, strain SH1000 through introduction of plasmid pSK5487M containing the gene aacA-aphD (encoding aminoglycoside N-acetyltransferase) and a chloramphenicol resistance gene (cat) for selection. NR-55216 was produced by resuspension of a lyophilized vial of deposited material in Tryptic Soy broth. Broth inoculum was added to Tryptic Soy broth containing 25 μg per mL chloramphenicol and grown for 1 day at 37°C in an aerobic atmosphere. The material from the initial growth was added to Tryptic Soy broth containing 25 μg per mL chloramphenicol, which was grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed using Tryptic Soy broth or agar containing 25 μg per mL chloramphenicol under propagation conditions unless otherwise noted.

Lot: 70052524 Manufacturing Date: 05MAY2022

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-positive cocci	Gram-positive cocci
Colony morphology	Report results	Circular, convex, entire, smooth and cream (Figure 1)
Motility (wet mount)	Report results	Non-motile
Hemolysis	Report results	β-hemolytic
Catalase	Positive	Positive
VITEK® MS (MALDI-TOF)	S. aureus	S. aureus (99.9%)
Antibiotic Susceptibility Profile ¹		,
Etest® antibiotic test strips		
1 day at 35°C in an aerobic atmosphere on Mueller Hinton agar		
Gentamycin	Resistant	Resistant (8 µg per mL) ²
Kanamycin	Resistant	Resistant (64 µg per mL)
Tobramycin	Resistant	Resistant (6 µg per mL)
Genotypic Analysis		
Digital DNA-DNA hybridization (dDDH) ³	≥ 70% for species identification	S. aureus (99.4%) ⁴
Next-Generation Sequencing (NGS) analysis for		, ,
antimicrobial resistance genes ⁵		
Gentamycin	Resistant	Resistant
Kanamycin	Resistant	Resistant
Tobramycin	Resistant	Resistant
Purity (post-freeze) 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar with 5% defibrinated sheep blood		Growth consistent with expected colony morphology
Viability (post-freeze)	Growth	Growth

¹Minimum Inhibitory Concentration (MIC); Clinical & Laboratory Standards Institute (CLSI) M100-S28 (2018) or MIC Interpretation Guideline: European Committee on Antimicrobial Susceptibility Testing (EUCAST) (2023).

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898


²The susceptibility result for this antibiotic is within one doubling dilution of specification, which is considered an equivalent result.

³Relatedness between bacterial strains has traditionally been determined using DDH. For additional information refer to Auch, A.F., et al. "Digital DNA-DNA Hybridization for Microbial Species Delineation by Means of Genome-to-Genome Sequence Comparison." <u>Stand. Genomic Sci.</u> 2 (2010): 117-134. PubMed: 21304684.



Certificate of Analysis for NR-55216

Figure 1: Colony Morphology



/Sonia Bjorum Brower/ Sonia Bjorum Brower

07 FEB 2023

Technical Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org

E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898

⁴The whole genome of *S. aureus*, strain AJUL2 (contig total length approximately 2.66 megabase pairs) was sequenced using the Illumina® MiSeq® system.

⁵In silico analysis of NGS data for antimicrobial resistance genes was performed using the Bacterial and Viral Bioinformatics Resource Center (BV-BRC), ResFinder and Pathogenwatch genome analysis tools.