

***Plasmodium falciparum*, Strain LA476-1Δhrp2**

**Catalog No. MRA-1331**

**Product Description:**

*Plasmodium falciparum* (*P. falciparum*), strain LA476-1Δhrp2 is deletion mutant derived from the progenitor strain LA476-1 (MRA-1330). Strain LA476-1 is a clone of *P. falciparum*, strain LA476, which was originally isolated in 2008 from a patient in Malawi. Strain LA476-1Δhrp2 was generated by the deletion of the *P. falciparum* histidine-rich protein 2 (*hrp2*) gene, located outside of the telomeric region of chromosomes 8, using CRISPR-Cas9 technology. MRA-1331 was produced by cultivation of deposited material in fresh human erythrocytes suspended in RPMI 1640 medium adjusted to contain 10% (v/v) heat-inactivated human serum (Type A), 25 mM HEPES, 2 mM L-glutamine, 2 g/L D-glucose, 27 µg/mL hypoxanthine and 5 µg/mL gentamicin. The culture was propagated in human Type O erythrocytes at 37°C in sealed flasks outgassed with a blood-gas atmosphere (90% N<sub>2</sub>, 5% CO<sub>2</sub>, 5% O<sub>2</sub>) and monitored for parasitemia for 22 days. Every 1 to 4 days, uninfected, leukocyte-filtered, Type O erythrocytes in complete culture medium were added dropwise as needed to maintain 1% to 3% hematocrit. The culture was harvested when the total parasitemia reached ≥ 2% to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

**Lot: 70064018**

**Manufacturing Date: 17OCT2023**

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: [Contact@BEIResources.org](mailto:Contact@BEIResources.org). We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
<b>Identification by Giemsa Stain Microscopy<sup>1</sup></b>	Blood-stage parasites present	Blood-stage parasites present
<b>Genotypic Analysis<sup>1</sup></b> Sequencing of Merozoite Surface Protein 2 (MSP2) gene (~ 670 base pairs)	Consistent with <i>P. falciparum</i>	Consistent with <i>P. falciparum</i> (Figure 1)
<b>Confirmation of Gene Deletion by PCR Amplification<sup>1,2</sup></b> <i>hrp2</i> (MRA-1331) <i>hrp2</i> (MRA-1330; positive control)	No amplicon ~ 300 base pair amplicon	No amplicon ~ 300 base pair amplicon
<b>Level of Parasitemia by Giemsa Stain Microscopy</b> Pre-freeze (22 days post-infection) <sup>3</sup> Ring-stage parasitemia Total parasitemia Post-freeze (2 days post-infection) <sup>1</sup> Ring-stage parasitemia Total parasitemia	Report results ≥ 2%  Report results ≥ 1%	4.4% 5.1%  4.6% 5.3%
<b>Viability (2 days post-infection)<sup>1</sup></b>	Growth in infected red blood cells	Growth in infected red blood cells
<b>Sterility (21-day incubation)<sup>1</sup></b> Harpo's HTYE broth, 37°C and 26°C, aerobic <sup>4</sup> Trypticase soy broth, 37°C and 26°C, aerobic Sabouraud broth, 37°C and 26°C, aerobic DMEM with 10% FBS, 37°C, aerobic Sheep blood agar, 37°C, aerobic Sheep blood agar, 37°C, anaerobic Thioglycollate broth, 37°C, anaerobic	No growth No growth No growth No growth No growth No growth No growth	No growth No growth No growth No growth No growth No growth No growth
<b>Mycoplasma Contamination<sup>1</sup></b> DNA detection by PCR	None detected	None detected

<sup>1</sup>Testing completed on vial, post-freeze material.

<sup>2</sup>Primer sequences and conditions for PCR are available upon request.

<sup>3</sup>Testing completed on bulk material prior to vialing and freezing.

<sup>4</sup>Atlas, Ronald M. *Handbook of Microbiological Media*. 3rd ed. Ed. Lawrence C. Parks. Boca Raton: CRC Press, 2004, p. 798.

**Figure 1: MRA-1331 MSP2 Sequence**

```
CAATGCTTATAATATGAGTATAAGGAGAAGTATGGCAAATGAAGGTTCTAATACTAATAAGTGTAGGTGCAGATGCTCCAAAAGCTGA
TACTATTGCTAGTGGAAGTCAAAGTAGTACAAATAGTGCAAGTACTAGTACTACTAATAATGGAGAATCACAACTACTACTCCTAC
CGCTGCTGATACCCCTACTGCTACAGAAAGTAATTCACCTTCACCACCCATCACTACTACAAAAAGTAATTCACCTTCACCACCCAT
CACTACTACAAAAAGTAATTCACCTTCACCACCCATCACTACTACAAAAAGTAATTCACCTTCACCACCCATCACTACTACAGAAAG
TTCAAGTTCTGGCAATGCACCAAATAAAACAGACGGTAAAGGAGAAGAGAGTGAAAAACAAAATGAATTAAATGAATCAACTGAAGA
AGGACCCAAAGCTCCACAAGAACCTCAAACGGCAGAAAATGAAAATCCTGCTGCACCAGAGAATAAAGGTACAGGACAACATGGACA
TATGCATGGTTCTAGAAATAATCATCCACAAAATACTTCTGATAGTCAAAAAGAATGTACCGATGGTAACAAAGAAAACCTGTGGAGC
AGCAACATCCCTCTTAAATAACTCTAGTAATATTGCTTCAATAAATAAATTTGTTGTTTTAATT
```

/Sonia Bjorum Brower/

Sonia Bjorum Brower

02 APR 2024

Technical Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

