SUPPORTING INFECTIOUS DISEASE RESEARCH

## Genomic RNA from Influenza B Virus, B/Macha/I1080/2021

### Catalog No. NR-59590

#### **Product Description:**

Genomic RNA was isolated from a preparation of cell lysate and supernatant from Madin-Darby canine kidney SIAT-1 (MDCK-SIAT1) cells infected with influenza B virus, B/Macha/I1080/2021 using QIAamp<sup>®</sup> Viral RNA Mini Kit (Qiagen<sup>®</sup> 52906). The viral genomic RNA is in a background of cellular nucleic acid and carrier RNA.

#### Lot: 70063759

## Manufacturing Date: 10JAN2024

| TEST   | SPECIFICATIONS               | RESULTS                                   |
|--|------------------------------|---|
| Genotypic Analysis<br>Sequencing of species-specific region<br>Hemagglutinin gene (~ 880 nucleotides)  | Consistent with source virus | Consistent with source virus <sup>1</sup> |
| Functional Activity by RT-PCR Amplification <sup>2</sup><br>Hemagglutinin gene   | ~ 1000 base pair amplicon    | ~ 1000 base pair amplicon                 |
| Estimated Concentration (post-dilution) by Qubit <sup>®</sup><br>Measurement (Viral, Cellular and Carrier) <sup>3</sup>                                    | Report results               | 4.6 ng per 100 μL (0.0462 μg/mL)          |
| Estimated Amount per Vial <sup>3</sup>   | Report results               | 4.6 ng                                    |
| Genome Copy Number Using BioRad QX200 Droplet<br>Digital PCR (ddPCR™) System<br>(Post vial; 9 replicates)  | Report results               | 4.8 × 10 <sup>8</sup> NDU/mL <sup>4</sup> |
| Virus Inactivation<br>10% of total yield inoculated on MDCK-SIAT1 cells and<br>evaluated for cytopathic effect and HA after serial<br>passage <sup>5</sup> | No viable virus detected     | No viable virus detected                  |

<sup>1</sup>Sequence information for influenza B virus B/Macha/I1080/2021 is not available in the NCBI database; nucleotide sequence obtained for NR-59590 lot 70063759 is identical to the source virus.

<sup>2</sup>Amplified using iTaq™ Universal SYBR Green One-step Kit (Bio-Rad<sup>®</sup> 172-5151) with 5 µL of NR-59590 in a 50 µL reaction

<sup>3</sup>Measurement is determined pre-vial prior to dilution due to the limit of detection of the quantification method

<sup>4</sup>NDU; NAAT-detectable units

<sup>5</sup>Use of the QIAamp<sup>®</sup> Viral RNA Mini Kit has been demonstrated to consistently inactivate 100% of influenza B viruses as shown by the absence of cytopathic effect (CPE) and HA after plating the entire extract on virus-susceptible cells for two passages.

# /Sonia Bjorum Brower/

11 APR 2024

Technical Manager or designee, ATCC Federal Solutions

ATCC<sup>®</sup>, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC<sup>®</sup>'s knowledge.



ATCC<sup>®</sup> is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.