SUPPORTING INFECTIOUS DISEASE RESEARCH

Staphylococcus aureus, Strain HIP15178

Catalog No. NR-46415

For research use only. Not for human use.

Contributor:

Jean Patel, PhD, Deputy Director, Office of Antimicrobial Resistance, Centers for Disease Control and Prevention, Atlanta, Georgia, USA

Manufacturer:

BEI Resources

Product Description:

Bacteria Classification: Staphylococcaceae, Staphylococcus Species: Staphylococcus aureus Strain: HIP15178 (also referred to as VRSA-5) NARSA Catalog Number: VRS5

- <u>Original Source</u>: *Staphylococcus aureus* (*S. aureus*), strain HIP15178 was isolated in 2005 in Michigan, USA from a surgical site infection of a 58-year-old female who was recently treated with an eight-week course of vancomycin.^{1,2}
- Comments: S. aureus, strain HIP15178 is a vancomycinresistant *S. aureus* (VRSA) strain.^{1,2} *S. aureus*, strain HIP15178 was deposited as positive for mec (subtype II) and vanA; negative for vanB, vanC1, vanC2, vanD, vanE, PVL and arginine catabolic mobile element (ACME); pulsed-field type USA100; spa repeats TJMBMDMGMK; Ridom spa type t002.1-5 S. aureus, strain HIP15178 is a USA100 isolate. USA100 isolates have the same MLST profile (ST 5), SCCmec (subtype II) and spa motif (MDMGMK) and are usually resistant to erythromycin and spectinomycin as well as being multiresistant to other commonly used therapeutic agents. USA100 is the most prevalent U.S health care-associated pulsed-field type and is endemic in many U.S. hospitals.³ The complete genome sequence of S. aureus, strain HIP15178 is available (GenBank: AHBO0000000).

S. aureus is a Gram-positive, cluster-forming coccus that normally inhabits human nasal passages, skin and mucus membranes. It is also a human pathogen and causes a variety of pus-forming infections as well as septicemia and endocarditis. S. aureus infections are difficult to treat due to resistance to numerous antibiotics. The development and dissemination of methicillin-resistant S. aureus (MRSA) strains has proven to be particularly difficult to contain and treat.6 Vancomycin has been the preferred antibiotic of choice for the treatment of MRSA infections, however, there have now been MRSA strains isolated that are also resistant to vancomycin.^{6,7} It is believed that this resistance results from either mutations that ultimately lead to a reduction of vancomycin at its site of action or from the acquisition of the vancomycin resistance gene, vanA, from Enterococcus faecalis.⁶⁻⁸ The vanA gene is carried by the Tn1546

transposon that resides on a plasmid in all VRSA strains.⁷ For VRSA strains carrying both *mecA* and *vanA*, β -lactams and glycopeptides seem to have a synergistic effect against these strains, both *in vitro* and in an animal model.^{8,9} Combination therapy, therefore, may be a more effective treatment option for VRSA infections than monotherapy with either antibiotic.^{8,9}

Material Provided:

Each vial contains approximately 0.5 mL of bacterial culture in Brain Heart Infusion broth supplemented with 6 μ g/mL vancomycin and 10% glycerol.

<u>Note</u>: If homogeneity is required for your intended use, please purify prior to initiating work.

Packaging/Storage:

NR-46415 was packaged aseptically in cryovials. The product is provided frozen and should be stored at -60°C or colder immediately upon arrival. For long-term storage, the vapor phase of a liquid nitrogen freezer is recommended. Freeze-thaw cycles should be avoided.

Growth Conditions:

<u>Note</u>: For stability purposes, it is recommended that the strain is subcultured in the presence of vancomycin.

Media:

Brain Heart Infusion broth or Tryptic Soy broth or equivalent

Brain Heart Infusion agar or Tryptic Soy agar with 5% defibrinated sheep blood or equivalent

Incubation:

Temperature: 37°C

Atmosphere: Aerobic

Propagation:

- 1. Keep vial frozen until ready for use, then thaw.
- 2. Transfer the entire thawed aliquot into a single tube of broth.
- 3. Use several drops of the suspension to inoculate an agar slant and/or plate.
- 4. Incubate the tube, slant and/or plate at 37°C for 18 to 24 hours.

Citation:

Acknowledgment for publications should read "The following reagent was provided by the Network on Antimicrobial Resistance in *Staphylococcus aureus* (NARSA) for distribution by BEI Resources, NIAID, NIH: *Staphylococcus aureus*, Strain HIP15178, NR-46415."

Biosafety Level: 2

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. <u>Biosafety in Microbiological and Biomedical Laboratories</u>. 5th ed. Washington, DC: U.S. Government Printing Office, 2009; see www.cdc.gov/biosafety/publications/bmbl5/index.htm.

E-mail: <u>contact@beiresources.org</u> Tel: 800-359-7370 Fax: 703-365-2898 **b**|**e**|**i** resources

SUPPORTING INFECTIOUS DISEASE RESEARCH

Disclaimers:

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at <u>www.beiresources.org</u>.

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC[®] nor the U.S. Government makes any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC[®] nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. ATCC[®] and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, ATCC[®], their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

Use Restrictions:

This material is distributed for internal research, noncommercial purposes only. This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale.

References:

- 1. NARSA, VRS5
- Sievert, D. M., et al. "Vancomycin-Resistant Staphylococcus aureus in the United States, 2002-2006." <u>Clin. Infect. Dis.</u> 46 (2008): 668-674. PubMed: 18257700.
- McDougal, L. K., et al. "Pulsed-Field Gel Electrophoresis Typing of Oxacillin-Resistant *Staphylococcus aureus* Isolates from the United States: Establishing a National Database." <u>J. Clin. Microbiol.</u> 41 (2003): 5113-5120. PubMed: 14605147.
- Saravolatz, L. D., J. Pawlak and L. B. Johnson. "In vitro Susceptibilities and Molecular Analysis of Vancomycin-Intermediate and Vancomycin-Resistant *Staphylococcus aureus* Isolates." <u>Clin. Infect. Dis.</u> 55 (2012): 582-586. PubMed: 22615331.
- Kos, V. N., et al. "Comparative Genomics of Vancomycin-Resistant Staphylococcus aureus Strains and Their Positions within the Clade Most Commonly Associated with Methicillin-Resistant S. aureus Hospital-Acquired Infection in the United States." <u>MBio.</u> 3 (2012): e00112-12.

PubMed: 22617140.

- Howden, B. P., et al. "Reduced Vancomycin Susceptibility in *Staphylococcus aureus*, Including Vancomycin-Intermediate and Heterogeneous Vancomycin-Intermediate Strains: Resistance Mechanisms, Laboratory Detection, and Clinical Implications." <u>Clin. Microbiol. Rev.</u> 23 (2010): 99-139. PubMed: 20065327.
- Courvalin P. "Vancomycin-Resistance in Gram-Positive Cocci." <u>Clin. Infect. Dis.</u> 42 (2006): S25-34. PubMed: 16323116.
- Severin, A., et al. "High Level Oxacillin and Vancomycin Resistance and Altered Cell Wall Composition in *Staphylococcus aureus* Carrying the Staphylococcal *mecA* and the Enterococcal *vanA* Gene Complex." <u>J.</u> <u>Biol. Chem.</u> 30 (2004): 3398-3407. PubMed: 14613936.
- Fox, P. M., et al. "Successful Therapy of Experimental Endocarditis Caused by Vancomycin-Resistant Staphylococcus aureus with a Combination of Vancomycin and Beta-Lactam Antibiotics." <u>Antimicrob.</u> <u>Agents Chemother.</u> 50 (2006): 2951-2956. PubMed: 16940087.

 $\mathsf{ATCC}^{\circledast}$ is a trademark of the American Type Culture Collection.

