



NIH AIDS Reagent Program

20301 Century Boulevard
Building 6, Suite 200
Germantown, MD 20874
USA

Phone: 240 686 4740
Fax: 301 515 4015
aidsreagent.org

DATA SHEET

Reagent:	Een217 T Cell Clone
Catalog Number:	755
Lot Number:	130236
Release Category:	A
Provided:	5 x 10 ⁶ cells/vial. Viability at cryopreservation time was 80%. Cells are cryopreserved in RPMI supplemented with 50% FBS and 10% DMSO.
Propagation Medium:	RPMI 1640 supplemented with 10% heat-inactivated fetal bovine serum, 4 mM L-glutamine, 50 U/ml penicillin, 50 µg/ml streptomycin, 50 U/ml recombinant human IL-2
Freeze Medium:	Fetal calf serum, 90%; DMSO, 10%.
Growth Characteristics:	For continued growth, these cells must be re-stimulated every 7-14 days with PHA and irradiated allogeneic peripheral blood mononuclear cells. Please see cell propagation instructions attached to data sheet.
Sterility:	Negative for aerobic and anaerobic bacteria, mycoplasma, fungi, and yeast.
Description:	CD4+ gp120-specific human cytolytic T-cell clone.
Special Characteristics:	This clone was derived from an HIV-seronegative donor by <i>in vitro</i> stimulation with recombinant gp120 followed by soft agar cloning. It recognizes aa 410-429 of HIV-1 pV22 gp120 in association with certain subtypes of DR4 (Dw10 and Dw15), and is cytolytic. Cells are 99% CD4 ⁺ , 99% CD3 ⁺ .
	<u>PROPAGATION OF HUMAN T-CELL CLONES</u>
Recommended Storage:	Liquid nitrogen.

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

Contributor: Dr. Robert F. Siliciano.

References: Preparation of irradiated PBMCs: Enhanced culture assay for detection and quantitation of latently infected, resting CD4+ T-cells carrying replication-competent virus in HIV-1-infected individuals. Siliciano JD, Siliciano RF. *Methods Mol Biol.* 2005;304:3-15.[Abstract](#)

Siliciano RF, Lawton T, Knal C, Karr R, Berman P, Gregory T, Reinherz E. Analysis of host-virus interactions in AIDS with anti-gp120 T cell clones: Effects of HIV sequence variation and a mechanism for CD4+ cell depletion. *Cell* **54**:561-575, 1988.

Callahan K, Fort M, Obah E, Reinherz E, Siliciano, R. Genetic variability in HIV-1 gp120 affects interactions with HLA molecules and T cell receptor. *J Immunol* **144**:3341-3346, 1990.

Polydefkis M, Koenig S, Flexner C, Obah E, Gebo K, Chakrabarti S, Earl P, Moss B, Siliciano R. Anchor sequence-dependent endogenous processing of human immunodeficiency virus 1 envelope glycoprotein gp160 for CD4+ T cell recognition. *J Exp Med* **171**:875-887, 1990.

NOTE: Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: Een217 T Cell Clone from Dr. Robert F. Siliciano (cat# 755)." Also include the references cited above in any publications.

Last Updated October 01, 2018

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.